

Miltex[®]
Redefining Excellence



MILTEX[®] DISPOSABLES
for ENT, GYN, Dermal and General Surgery



MILTEX® DISPOSABLES

Enjoy the convenience of disposables with the quality you want. Miltex offers an extensive line of single-use instrumentation manufactured to provide the feel and performance of stainless steel instruments.

Contact us for a free trial sample of our Miltex® Disposable items to truly appreciate the quality and performance of our products.

TABLE OF CONTENTS

Biopsy Punches	2-3
Disposable Scalpels	4-5
Surgical Disposables & Accessories	6-7
GYN Disposables	8-10
ENT Disposables	11-13
Sparta Products	14
General Instruments	15
Ear Curettes	16
Mirrors	16

Illustrations and content provide general description only and may be subject to minor change.

Note: Some products in this catalog are not available for sale in the EU.



HOW MILTEX CREATES A SUPERBLY SHARP PUNCH... EVERYTIME ...AND EVERYTIME... A PERFECT BIOPSY!

Miltex produces the **SHARPEST** seamless edge of any biopsy punch when compared to other tested brands!



- **25% more measured cutting sharpness (KgF)* than other tested punches**
- **Significant metal hardness (HV)* begins the process of creating sharpness**
- **A steeper edge angle (by degree)* is a sharpness design unique to Miltex**
- **A precise punch diameter* for a true and accurate biopsy**

* KgF (Kilogram Force). The Miltex® Punch requires less weighted "force" or push from the clinician to achieve the desired biopsy depth as compared to other brands measured KgF.

* HV (Vickers Hardness). The Miltex® Punch (metal) tip "measures" harder than other tested brands employing this metallurgical standard and results in an extremely sharper edge.

* Degree. The Miltex® Punch is designed at a steeper angle than tested competitive brands to ensure a finely honed, sharper edge.

* Diameter. The Miltex® Punch is manufactured to accurately measure to the specified size.

* Data on file


Biopsy Punch, Sterile

Now you can have the convenience of disposable punches along with the quality you expect from Miltex. Miltex offers an extensive product line of sterile disposable biopsy punches that are manufactured from the finest grade surgical steel to ensure superior cutting every time.

- Seamless, stainless steel, razor-sharp cutting edges eliminate torn or pulled biopsies
- Ribbed handle provides comfort and control
- Convenience of disposables, single-use
- Sterile, 5-year shelf life
- Sizes clearly marked on each punch
- Individual color-coded blister packs

Biopsy Punch with Plunger System


- Effortless ejection of lodged specimens
- Available in boxes of 25
- Color-coded by size - 5 sizes available



REF	Size
33-31AA-P/25	1.0mm w/plunger
33-31A-P/25	1.5mm w/plunger
33-31-P/25	2.0mm w/plunger
33-32-P/25	3.0mm w/plunger
33-34-P/25	4.0mm w/plunger

Biopsy Punch, Standard

- Each dispenser box contains 50 punches
- Ten sizes for a wide range of applications
- Assortment box contains 10 each of the following sizes: 2mm, 3mm, 4mm, 5mm, and 6mm



REF	Size
33-31AA	1.0mm
33-31A	1.5mm
33-31	2.0mm
33-31B	2.5mm
33-32	3.0mm
33-33	3.5mm
33-34	4.0mm
33-35	5.0mm
33-36	6.0mm
33-37	8.0mm
33-38	Assorted sizes*

Each box contains 50 sterile color-coded blister packs.

*10 each of 2mm, 3mm, 4mm, 5mm, and 6mm.

Safety Scalpels, Sterile

- Complies with OSHA directive for the Needlestick Safety & Prevention Act
- Stainless steel, razor-sharp cutting edge for superior cutting and clean incisions
- Pre-assembled, fully retractable scalpel blade is user activated
- Sterile, 5-year shelf life
- Individually packaged in sealed peel-apart pouches
- Convenience of disposables, single-use
- Metric ruler on each scalpel to assist in determining incision length
- Sizes clearly marked on each scalpel blade
- Each dispenser box contains 10 safety scalpels
- Six sizes available

REF	Size
4-510	10
4-511	11
4-512	12
4-515	15
4-515C	15C Angled Blade
4-522	22



4-510



4-511



4-515



4-515C




4-512



4-522

Directions For Use

Gently push slider to right and completely forward in one simple motion. Blade locks into place with audible confirmation. To retract, push slider to right and back to return blade to original protected position.



Scalpels, Standard, Sterile

- Stainless steel, razor-sharp cutting edge for superior cutting and clean incisions
- Safety sleeve for blade protects you and the blade
- Sterile, 5-year shelf life
- Individually packaged in sealed peel-apart pouches
- Convenience of disposables, single-use
- Sizes clearly marked on each scalpel blade
- Each dispenser box contains 10 scalpels
- Eight sizes available

REF	Size
4-410	10
4-411	11
4-412	12
4-415	15
4-420	20
4-421	21
4-422	22
4-423	23



4-410



4-411



4-412



4-415



4-420



4-421



4-422



4-423

Miltex BiopBlade™, Sterile

- Ideal for skin/tissue saucerization and shave biopsy
- Used in dermatology, podiatry, family practice and plastic surgery
- Incorporates a super sharp blade with rear finger guard
- Flexible design allows the blade to be positioned at any angle for the intended procedure
- Disposable, one-time use only
- Sterile, 5-year shelf life




REF
33-100 Bx/50

Dermal Curettes, Sterile

- Seamless, stainless steel, semi-sharp cutting edge featuring a unique symmetrically ground beveled edge
- Textured handle provides better comfort and control
- Sterile, 5-year shelf life
- Individually packaged in peel-apart pouches
- Convenience of disposables, single-use
- Sizes clearly marked on each curette
- Each dispenser box contains 50 curettes
- Five sizes available for a wide range of applications




REF	Size
33-52	2.0mm
33-53	3.0mm
33-54	4.0mm
33-55	5.0mm
33-57	7.0mm



Surgical Blades, Sterile

Carbon or Stainless Steel

- Razor-sharp cutting edge for superior cutting and clean incisions
- Sterile, 5-year shelf life
- Individually packaged in foil peel packs
- Convenience of disposables, single-use
- Sizes clearly marked on each blade
- Available in carbon and stainless steel
- Each dispenser box contains 100 blades
- Eleven sizes available



CS*	SS**	Size
REF	REF	
4-110	4-310	10
4-111	4-311	11
4-112	4-312	12
4-112B ^A	4-312B ^A	12B ^A
4-115	4-315	15
4-115C ^{AA}	4-315C ^{AA}	15C ^{AA}

*Carbon Steel ** Stainless Steel

^AB is a double edged blade

^{AA}C is an angled blade

CS*	SS**	Size
REF	REF	
4-120	4-320	20
4-121	4-321	21
4-122	4-322	22
4-123	4-323	23
4-125	4-325	25

*Carbon Steel ** Stainless Steel

Scalpel Handles

For your convenience, Miltex, Inc. offers a wide array of scalpel blade handles and the Miltex® Blade Safe Remover to aid in the safe removal of blades and to help protect against accidental cuts and punctures.




4-100

REF	Description	Fits Blades 10-15C	Fits Blades 20-25
4-7	No. 3, Stainless Steel 5" graduated in mm and cm	x	
4-7R	No. 3, Nickel 5" graduated in mm and cm	x	
4-8	No. 4, Stainless Steel 5.375" graduated in mm and cm		x
4-8R	No. 4, Nickel 5.375" graduated in mm and cm		x
4-9	No. 7, Stainless Steel 6.5"	x	
4-10*	No. 3L, Stainless Steel 8.5"	x	
4-11*	No. 4L, Stainless Steel 8.5"		x
4-12*	No. 3LA, Stainless Steel, angled 8.5"	x	
4-16	No. 9, Stainless Steel 4.5"	x	
4-17	SIEGEL, Stainless Steel, round handle 5.75"	x	
4-18	No. 5, plastic handle 5", green	x	
4-20	No. 6, plastic handle 5.25", red		x
4-30**	No. 3, fitting blades 5.5", grey	x	
4-31**	No. 4, fitting blades 5.5", blue		x
4-400***	Miniature Blade Handle w/chuck 3-3/4", round knurled (3K type)		
4-401***	Miniature Blade Handle w/self-locking chuck, 3-3/4", round knurled (3K type)		
4-405***	Self-locking chisel blade handle, Stainless Steel 5.31"		
4-100	Blade Safe surgical blade remover for any style handle		

* For deep surgery

** EZ-OFF reusable plastic knife handles, autoclavable, permitting quick and safe removal of blade, 10/box

*** Fits Mini Blade (blades not available from Miltex, Inc.)



4-7 4-8 4-9 4-10 4-11 4-12 4-16 4-17 4-18 4-20 4-30 4-31 4-400 4-401 4-405
4-7R 4-8R

Miltex® Disposable Electrodes


FEATURES

- Range of most popular patterns and sizes
- Loop and Square working ends manufactured with Tungsten Wire to maintain shape throughout procedure
- 5.1" (13 cm) shaft
- 3/32" (2.4 mm) stainless steel connector shaft fits most standard electrosurgical pencils and generators
- Each Loop and Electrode individually packaged and sterilized in peel pouch
- Disposable, one-time use only
- Sterile, 5-year shelf life
- Packaged 5 per box

REF	Description*
31-5000	Ball (3mm Diameter)
31-5002	Ball (5mm Diameter)
31-5022	Square (10mm width x 8mm depth)
31-5054	Loop (10mm width x 8mm depth)
31-5056	Loop (10mm width x 10mm depth)
31-5070	Loop (15 mm width x 5mm depth)
31-5072	Loop (15mm width x 8mm depth)
31-5074	Loop (15mm width x 10mm depth)
31-5090	Loop (20mm width x 8mm depth)
31-5092	Loop (20mm width x 10mm depth)
31-5094	Loop (20mm width x 12mm depth)
31-5096	Loop (20mm width x 15mm depth)
31-5098	Loop (20mm width x 20mm depth)

*Sold in box of 5

Images signify general design only.




Contact Miltex Customer Service for more information and a Brochure of our Insulated Stainless Steel Instruments for Electrosurgical Procedures

Disposable Os Locator, Canal Dilator & Fundus Sound

Three convenient disposable and latex free patterns. The Os Locator is used to locate the external Os, the Canal Dilator is designed to dilate the cervical canal and the Fundus Sound assists with determining uterine depth. Each pattern is packaged individually for specific usage.

- Packaged 25 per box
- Disposable, one-time use only
- Sterile, 3-year shelf life
- Latex-free

REF	Description
30-5001	Os Locator (25/Box), sterile
30-5002	Canal Dilator (25/Box), sterile
30-5003	Fundus Sound (25/Box), sterile




Uterine Sound

Miltex® Disposable Uterine Sounds are flexible and sterile and are measured in metric units. The soft design is a helpful alternative to traditional patterns manufactured with stainless steel, silver or chrome plating.

- Packaged 25 per box
- Disposable, one-time use only
- Sterile, 3-year shelf life
- Latex-free

REF	Description
30-6000	Uterine Sound (25/Box), sterile



Endometrial Sampling Set

The Miltex® Endometrial Sampling Set is cost effective and accurate for screening the endometrium. It is ideal for efficient sampling for both pre- and post-menopausal patients and affords reduced physical discomfort often associated with endometrial tissue collection.

Slight Curve on Curette

Optimizes access to the uterine contour and structure for obtaining a specimen from the endometrial cavity.

3 mm O.D Curette

Allows insertion and permits gentle curette placement for sampling while reducing the risk of uterine perforation.

FEATURES

- Set includes one each of Curette and Luer Lock Syringe
- 10 ml standard Twist-and-Lock Syringe
 - Creates 21" Hg continuous negative pressure for optimized suction and sample capture
 - Prevents backflow of sample into syringe
 - Luer Lock fitting has winged flanges for ease of grasp and attachment
- Curette working length = 8-1/8" (20.64 cm)
- Disposable and latex-free
- Packaged 10 per box
- Sterile, 3-year shelf life

Semi-Rigid Curette with Tapered Tip

For controlled insertion. Ideal for patients where stenosis or atrophic endometrium is a concern.

Randall-Type Sharp Cutting Edge

Efficient sample collection.

REF	Description
30-3015	Endometrial Sampling Set (includes Disposable Curette with a Twist-and-Lock Syringe), sterile, single-use, box of 10



Endometrial Biopsy Cannula

Endometrial sampling is the most widely used screening technique of the female patient population for a number of indications. The Miltex® SoftFlex™ and FirmFlex™ Endometrial Biopsy Cannula are ideal for accurate as well as economical in-office specimen collection of endometrial tissue.

TWO STYLES

SoftFlex™

- For standard/routine sampling
- When a more flexible cannula is desired

FirmFlex™

- For sampling when stenosis or atrophic endometrium is a concern
- For post menopausal women where stenosis of the cervix is a concern
- When a more rigid cannula is needed for controlled insertion

FEATURES

- One step sampling – no extra syringe or suction required
- Cannula and inner sleeve design generate sufficient vacuum for maximum sample yield without the need for external suction device
- 23.5 cm length allows for insertion depth to accommodate a range of patients
- Soft and Firm flexibility choices optimize formation to the natural uterine contour and structure for obtaining specimens from the endometrial cavity
- Disposable and latex-free
- Packaged 25 per box
- Sterile, 3-year shelf life


REF	Description
30-3010	SoftFlex™ Endometrial Biopsy Cannula, sterile, single-use, box of 25
30-3012	FirmFlex™ Endometrial Biopsy Cannula, sterile, single-use, box of 25

3 mm O.D.

Appropriate size for ease of insertion without dilation in most cases.

3 cm - 10 cm Graduated Markings & Tapered Tip

Aids proper insertion and placement of the cannula during specimen collection while reducing risk of uterine perforation.



ENT Packing Line

The Miltex® ENT Packing line is a full offering of Polyvinyl Alcohol dressings for nasal, sinus and ear surgeries. The soft pre-compressed micro-pore size material allows for quick insertion and atraumatic removal for patient comfort. The Miltex® ENT Packing line is 100% lint and fiber free, easily trimmed for a custom fit and most products are available with a drawstring for patient safety.

FEATURES

- Hydrophilic and fast wicking
- Absorbs up to 25 times its weight in fluid
- Supplied compressed for easier insertion
- Disposable, one-time use only
- Sterile, 4-year shelf life

EAR PACKING

- Excellent wicking capabilities for administering otic medication in the treatment of Swimmer's Ear and other post-otic surgery.
- Fits into swollen ear canal and continually expands during healing process; helps provide contact and medication delivery to the infection site.

REF	Description
19-2000	Short Otitis Externa Wick (9 mm x 15 mm), (1/pkt - 50 pkts/box)
19-2002	Long Otitis Externa Wick (9 mm x 24 mm), (1/pkt - 10 pkts/box)
19-2004	Ear Packing with string (12 mm x 15 mm), (1/pkt - 10 pkts/box)



19-2000



19-2002



19-2004

EPISTAXIS

- Highly absorbent and effective control of epistaxis
- Inserts in seconds vs. traditional gauze packing method requiring 4 to 6 feet of gauze.
- Removal is easy and atraumatic.

REF	Description
20-3020	Small Epistaxis Dressing w/string, (5.5 cm x 1.5 cm x 2.5 cm), (1/pkt - 20 pkts/box)
20-3022	Large Epistaxis Dressing w/string, (10 cm x 1.5 cm x 2.5 cm), (1/pkt - 10 pkts/box)



20-3020



20-3022

NASAL

- Nasal products can be used for absorption and drainage for Rhinoplasty, Septoplasty and other nasal post-operative procedures.
- Non-stick skinned surface on dressing helps with atraumatic removal and minimizes re-bleeding.
- Dressings with airway (cannulated) tubes help ensure patient comfort.

REF	Description
20-3000	Small Standard Nasal Dressing w/string (4.5 cm x 1.5 cm x 2 cm), (1/pkt - 10 pkts/box)
20-3002	Small Standard Nasal Dressing with airway tube and string (4.5 cm x 1.5 cm x 2 cm), (1/pkt - 10 pkts/box)
20-3004	Large Standard Nasal Dressing w/string (8 cm x 1.5 cm x 2 cm), (1/pkt - 10 pkts/box)
20-3006	Large Standard Nasal Dressing with airway tube and string (8 cm x 1.5 cm x 2 cm), (1/pkt - 10 pkts/box)
20-3010	Thin Anatomical Nasal Dressing w/string (8 cm x 0.75 cm x 3 cm), (1/pkt - 10 pkts/box)
20-3012	Standard Anatomical Nasal Dressing with string (8 cm x 1.5 cm x 3 cm), (1/pkt - 10 pkts/box)
20-3014	Standard Anatomical Nasal Dressing with airway tube and string (8 cm x 1.5 cm x 3 cm), (1/pkt - 10 pkts/box)



20-3000



20-3002



20-3004



20-3006



20-3010



20-3012



20-3014

SINUS

- Designed to prevent lateralization of the middle turbinate after Functional Endoscopic Sinus Surgery.
- Compressed configuration for easy small cavity placement.
- Absorbs post-operative drainage and helps clear clots and debris upon removal, leaving a clean cavity.

REF	Description
20-3030	Sinus Dressing B-Type w/string (2.5 cm x 0.9 cm x 2 cm), (2/pkt - 10 pkts/box)
20-3032	Sinus Dressing K-Type w/string (3.5 cm x 0.9 cm x 1.2 cm), (2/pkt - 10 pkts/box)
20-3034	Sinus Dressing T-Type w/string (4.5 cm x 0.9 cm x 2 cm), (2/pkt - 10 pkts/box)
20-3036	Slim Sinus Dressing B-Type w/string (2.5 cm x 0.6 cm x 2 cm), (2/pkt - 10 pkts/box)
20-3038	Slim Sinus Dressing K-Type w/string (3.5 cm x 0.6 cm x 1.2 cm), (2/pkt - 10 pkts/box)
20-3040	Slim Sinus Dressing T-Type w/string (4.5 cm x 0.6 cm x 2 cm), (2/pkt - 10 pkts/box)
20-3042	Sinus Dressing Kit - Assortment of Sinus Dressing, Nasal Dressing/Airway Tube. (1/pkt - 10 pkts/box)



20-3030



20-3032



20-3034



20-3036



20-3038



20-3040



20-3042

NASOSTAT™ Treatment of Epistaxis

- Nasostat™ products may be used in conjunction with anterior packing
- Will control bleeding while minimizing patient discomfort associated with nasal packs
- Is a safe and easily tolerated alternative to other treatments for epistaxis, e.g. cautery
- Sterile, 1-year shelf life

REF	Description
41900	Nasostat Regular, Sterile, 5/Bx
41901	Nasostat Slimline, Sterile, 5/Bx
41902	Nasostat Large, Sterile, 5/Bx
41903	Nasostat Small, Sterile, 5/Bx



SURGI-CLAMP, PLASTIC

- Non-sterile operator clamp
- Multi-functional
- Steam sterilization safe
- Designed for reuse
- No rust or stain problems

REF	Description	
374201	Surgi-Clamp, Plastic	100/bx



SURGI-PREP, DEPILATORY

- Fast-acting removal of body hair; dissolves unwanted hair below skin surface
- Lasts days longer than shaving
- Pain free formula; removes even coarse, stubborn hair
- No nicks, cuts, or prickly stubble compared to shave preparation
- Works in 10 minutes

REF	Description	
600900	Surgi-Prep Depilatory	12/bx



SURGI-BRUSH, STERILE

- Sterile brush used for cleaning surgical instruments (in the OR)
- Double-ended nylon bristle brush
- One end has 3 rows of bristles for general decontamination
- One end contains a single row of bristles for fine, delicate cleaning
- Disposable, one-time use only
- Sterile, 1-year shelf life

REF	Description	
405043	Surgi-Brush, Sterile	20/bx



SURGICAL INSTRUMENTS, FLOOR GRADE QUALITY, STERILE

SINGLE-USE INSTRUMENTS


- Convenience and cost savings provide an exceptional value
- Eliminate the time and expense of re-processing soiled instruments
- Includes all of the most common instrument patterns for routine surgical needs
- Sterile, 2-year shelf life



REF	Description	
ST5-120	MAYO, 5-1/2", ST, B/B	50/Bx
ST5-124	MAYO, 6-3/4", ST, B/B	50/Bx
ST5-16	OR Scissors, 5-1/2", ST, S/B	50/Bx
ST5-180	METZENBAUM, 5-1/2", Curved	50/Bx
ST5-182	METZENBAUM, 7", Curved	50/Bx
ST5-304	IRIS Scissors, 4-1/2", Straight	50/Bx
ST5-306	IRIS Scissors, 4-1/2", Curved	50/Bx
ST5-504SS	LISTER Bandage Scissors, 5-1/2"	50/Bx
ST9-102	LITTAUER Stitch Scissors, 4-1/2"	50/Bx
SCISSORS		
ST6-118	ADSON Forceps, 4-3/4"	50/Bx
ST6-120	ADSON Forceps, 4-3/4", 1X2	50/Bx
ST6-44	TISSUE Forceps, 5-1/2", 1X2	50/Bx
ST6-8	DRESSING Forceps, 5-1/2"	50/Bx
ST7-138	ROCHESTER-PEAN, 6-1/4", Curved	50/Bx
ST7-2	HALSTEAD Forceps, 5", Straight	50/Bx
ST7-4	HALSTEAD Forceps, 5", Curved	50/Bx
ST7-36	KELLY Forceps, 5-1/2", Straight	50/Bx
ST7-38	KELLY Forceps, 5-1/2", Curved	50/Bx
ST7-42	CRILE Forceps, 5-1/2", Straight	50/Bx
ST7-44	CRILE Forceps, 5-1/2", Curved	50/Bx
FORCEPS & HEMOSTATS		
ST8-44	MAYO-HEGAR Needle Holder, 6"	50/Bx
ST8-50	CRILE-WOOD Needle Holder, 6"	50/Bx
ST8-6	WEBSTER Needle Holder, 5"	50/Bx
NEEDLE HOLDERS		

EAR CURETTES, DISPOSABLE

- Each curette is individually wrapped, non-sterile
- Available in boxes of 50
- Color-coded and available in 4 tip configurations



REF	Description
19-320	Ring Tip, Pink
19-321	Loop Tip, White
19-322	Spoon Tip, Gray
19-323	Large Spoon Tip, Peach

DISPOSABLE MIRRORS

- Disposable plane surface mirrors
- Individually wrapped, non-sterile
- Used for air abrasion procedures
- Available in boxes of 72



017-26810



017-26811

REF	Description
017-26810	Plane Surface Mirror, Size 4, Economical, 72/box
017-26811	Plane Surface Mirror, Size 4, Economical Deluxe, 72/box

DENLITE® ILLUMINATED MIRROR

The DenLite® Illuminated Mirror System is designed to give Dental and Medical Professionals a shadow-free view of the mouth during routine oral and dental exams. It provides a bright white light which eliminates the need for constant adjustment of cumbersome overhead lights during everyday procedures. The DenLite® system is the only cordless, illuminated instrument that can be used as a mirror while also being strong enough to be also used as a retractor during an examination.

- Cordless illuminated mirror system
- Provides shadow-free view of mouth
- Strong enough to be used as a retractor
- Mirrors sold separately in boxes of 50, non-sterile



DP5000

DP5030

REF	Description
DP5000	DenLite LED Kit
DP5150	DenLite LED Light Cartridge
DP5200	DenLite Handle Set
DP5400	DenLite Charger Set
DP5030	DenLite Disposable Mirrors, single use (50 mirrors per box - used with DP5000)
DP5250	DenLite Handpiece Collet Assembly
DP5070	DenLite Illuminator Tip
DP1030-3	DenLite Disposable Mirror, 110V, 3 boxes of 72 mirrors (used with DP1000 AND DP1500)
DP1020	DenLite Lamp for DP1000
DP1520	DenLite Lamp for DP1500



www.miltex.com
toll-free phone 866-854-8300
toll-free fax 866-854-8400
phone 717-840-9335
fax 717-840-9347

Please call Customer Service today for additional information, an in-office demonstration, to place an order, or to receive other Miltex catalogs.

Miltex, Miltex Redefining Excellence and DenLite are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. SoftFlex, FirmFlex and Miltex BiopBlade are trademarks of Miltex, Inc. The Integra wave logo is a trademark of Integra LifeSciences Corporation. Any other trademarks identified in this brochure are the property of their respective owners. Copyright© 2009 Integra LifeSciences Corporation

