

MeisterHand *By Miltex*



OR-GRADE SURGICAL INSTRUMENTS



## I N T R O D U C T I O N

MeisterHand® Instruments are expertly crafted from surgical grade stainless steel German forgings and meet the very highest Operating Room quality standards and specifications.

This new instrument line is advantageously priced while offering greater quality and value compared to instruments from any other source.

A wide range of General Instruments are included in the MeisterHand® line as well as an extensive selection of Specialty Instrumentation ideal for a variety of surgical procedures and settings.





## WARRANTY


Miltex hereby warrants all MeisterHand® instruments purchased directly from Miltex or an Authorized Dealer in accordance with the following terms and conditions:

All MeisterHand® instruments are guaranteed to be free from defects in material and workmanship when used for their intended surgical purpose. Any MeisterHand® instrument that proves defective in workmanship or material will either be repaired or replaced, at the discretion of Miltex, without charge.

The above is a limited warranty and is the only warranty made by Miltex. Miltex makes no other warranty, express or implied, including any warranty of merchantability or fitness for a particular purpose.

All instruments shown in this catalog are made of stainless steel, unless otherwise indicated.

Illustrations and content provide general description only and may be subject to change.



## TABLE OF CONTENTS

Scalpel Handles . . . . .	6	Ear Instruments . . . . .	104
Scissors . . . . .	7	Nasal Instruments . . . . .	113
Tungsten Carbide		Plastic Surgery	
Scissors . . . . .	26	Instruments . . . . .	121
SuperCut Scissors . . . . .	29	Dermal Instruments . . . . .	152
Dressing & Tissue		Mouth & Throat	
Forceps . . . . .	32	Instruments . . . . .	168
Tungsten Carbide		Tracheal Instruments . . . . .	176
Dressing Forceps . . . . .	44	Cardiovascular & Thoracic	
Splinter Forceps . . . . .	45	Instruments . . . . .	178
Hemostatic Forceps . . . . .	46	Neurosurgical	
Towel Clamps . . . . .	55	Instruments . . . . .	189
Sponge Forceps . . . . .	56	Orthopedic Instruments . . . . .	206
Needle Holders . . . . .	57	Podiatry Instruments . . . . .	219
Tungsten Carbide		Rectal Instruments . . . . .	220
Needle Holders . . . . .	61	Urological Instruments . . . . .	222
Suction Tubes . . . . .	65	Obstetrics &	
Retractors . . . . .	67	Gynecology . . . . .	224
Hooks . . . . .	81	Miscellaneous	
Micro Surgery		Instruments . . . . .	239
Instruments . . . . .	86	Alphabetical Index . . . . .	241
Ophthalmic Instruments . . . . .	89	Numerical Index . . . . .	251



MH4-7



MH4-8



MH4-9



MH4-10

REF	Description
MH4-7*	Knife Handle no. 3, fitting surgical blades nos. 10 thru 15c, extra fine
MH4-8*	Knife Handle no. 4, fitting surgical blades nos. 20 thru 25, extra fine
MH4-9	Knife Handle no. 7, fitting surgical blades nos. 10 thru 15c, extra fine
MH4-10	Knife Handle no. 3L, for deep surgery, fitting surgical blades nos. 10 thru 15c, extra fine.

\*Handles have mm and cm graduations.



MH5-2 to MH5-8



MH5-12 to MH5-18



MH5-22 to MH5-26

### OPERATING SCISSORS

REF	Description
MH5-2	Operating Scissors, straight, 4-1/2" (11.4 cm), sharp-sharp points
MH5-4	Operating Scissors, straight, 5" (12.7 cm), sharp-sharp points
MH5-6	Operating Scissors, straight, 5-1/2" (14 cm), sharp-sharp points
MH5-8	Operating Scissors, straight, 6-1/2" (16.5 cm), sharp-sharp points
MH5-12	Operating Scissors, straight, 4-1/2" (11.4 cm), sharp-blunt points
MH5-14	Operating Scissors, straight, 5" (12.7 cm), sharp-blunt points
MH5-16	Operating Scissors, straight, 5-1/2" (14 cm), sharp-blunt points
MH5-18	Operating Scissors, straight, 6-1/2" (16.5 cm), sharp-blunt points
MH5-22	Operating Scissors, straight, 4-1/2" (11.4 cm), blunt-blunt points
MH5-24	Operating Scissors, straight, 5" (12.7 cm), blunt-blunt points
MH5-26	Operating Scissors, straight, 5-1/2" (14 cm), blunt-blunt points





MH5-36




MH5-46




MH5-56

### OPERATING SCISSORS

REF	Description
MH5-36	Operating Scissors, curved, 5-1/2" (14 cm), sharp-sharp points
MH5-46	Operating Scissors, curved, 5-1/2" (14 cm), sharp-blunt points
MH5-56	Operating Scissors, curved, 5-1/2" (14 cm), blunt-blunt points




MH5-120  
MH5-124  
MH5-128



MH5-122  
MH5-126  
MH5-130

REF	Description
MH5-120	MAYO Dissecting Scissors, 5-1/2" (14 cm), straight, standard beveled blades
MH5-122	MAYO Dissecting Scissors, 5-1/2" (14 cm), curved, standard beveled blades
MH5-124	MAYO Dissecting Scissors, 6-3/4" (17.1 cm), straight, standard beveled blades
MH5-126	MAYO Dissecting Scissors, 6-3/4" (17.1 cm), curved, standard beveled blades
MH5-128	MAYO Dissecting Scissors, 9" (22.9 cm), straight, standard beveled blades
MH5-130	MAYO Dissecting Scissors, 9" (22.9 cm), curved, standard beveled blades



MH5-179  
MH5-180

MH5-181

MH5-182

MH5-288

REF	Description
MH5-179	METZENBAUM (LAHEY) Scissors, Standard Pattern, 5-1/2" (14 cm), straight, blunt points
MH5-180	METZENBAUM (LAHEY) Scissors, Standard Pattern, 5-1/2" (14 cm), curved, blunt points
MH5-181	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), straight, blunt points
MH5-182	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), curved, blunt points
MH5-288	METZENBAUM Scissors, Delicate Pattern, 7" (17.8 cm), curved, blunt points



MH5-270  
MH5-272

MH5-274  
MH5-276

MH5-278

MH5-280

## PLASTIC SURGERY AND DELICATE SCISSORS

REF	Description
MH5-270	Delicate Scissors, 4-3/4" (12.1 cm), straight, sharp-sharp points
MH5-272	Delicate Scissors, 4-3/4" (12.1 cm), straight, sharp-blunt points
MH5-274	Delicate Scissors, 4-3/4" (12.1 cm), straight, blunt-blunt points
MH5-276	Delicate Scissors, 4-3/4" (12.1 cm), curved, sharp-sharp points
MH5-278	Delicate Scissors, 4-3/4" (12.1 cm), curved, sharp-blunt points
MH5-280	Delicate Scissors, 4-3/4" (12.1 cm), curved, blunt-blunt points



**MH5-300**  
**MH5-304**



**MH5-302**  
**MH5-306**



**MH18-1396**



**MH18-1398**

REF	Description
MH5-300	Iris Scissors, 4" (10.2 cm), straight
MH5-302	Iris Scissors, 4" (10.2 cm), curved
MH5-304	Iris Scissors, 4-1/2" (11.4 cm), straight
MH5-306	Iris Scissors, 4-1/2" (11.4 cm), curved
MH18-1396	Iris Scissors, 3-1/2" (8.9 cm), straight, with 20 mm blades, delicate
MH18-1398	Iris Scissors, 3-1/2" (8.9 cm), curved, with 20 mm blades, delicate



MH18-1415

MH18-1416

MH18-1417

MH18-1418

REF	Description
MH18-1415	Iris Scissors, 4" (10.2 cm), straight, standard pattern, ribbon-type
MH18-1416	Iris Scissors, 4" (10.2 cm), curved, standard pattern, ribbon-type
MH18-1417	Iris Scissors, 3-3/4" (9.5 cm), straight, miniature pattern with 14 mm blades, ribbon-type
MH18-1418	Iris Scissors, 3-3/4" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon-type



MH18-1420



MH18-1422



MH18-1426



MH18-1428

REF	Description
MH18-1420	KNAPP Iris Scissors, 4" (10.2 cm), straight, sharp-blunt points
MH18-1422	KNAPP Iris Scissors, 4" (10.2 cm), straight, blunt-blunt points
MH18-1426	KNAPP Iris Scissors, 4" (10.2 cm), curved, sharp-blunt points
MH18-1428	KNAPP Iris Scissors, 4" (10.2 cm), curved, blunt-blunt points



MH5-312  
MH5-314

MH5-264

REF	Description
MH5-312	Strabismus Scissors, 4" (10.2 cm), straight
MH5-314	Strabismus Scissors, 4" (10.2 cm), curved
MH5-264	DEAN Dissecting Scissors, 6-3/4" (17.1 cm), blades angled on flat, one serrated blade





MH18-1460



MH18-1462 to MH18-1466



MH18-1472 to MH18-1476

### STEVENS TENOTOMY SCISSORS

REF	Description
MH18-1460	4-1/8" (10.5 cm), straight, short blades, sharp points
MH18-1462	4-1/8" (10.5 cm), straight, short blades, blunt points
MH18-1464	4-1/8" (10.5 cm), curved, short blades, sharp points
MH18-1466	4-1/8" (10.5 cm), curved, short blades, blunt points
MH18-1472	4-1/2" (11.4 cm), straight, long blades, blunt points
MH18-1476	4-1/2" (11.4 cm), curved, long blades, blunt points



MH18-1572



MH18-1576



MH18-1578



MH18-1620



MH18-1622

REF	Description
MH18-1572	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), angled blades, 11 mm, blunt points
MH18-1576	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), curved, sharp points
MH18-1578	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), curved, blunt points
MH18-1620	VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), straight, extra delicate sharp points
MH18-1622	VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), curved, extra delicate sharp points



MH20-160



MH21-600



MH21-602



MH21-604

REF	Description
MH20-160	KNIGHT Nasal Scissors, 6-3/4" (17.1 cm), angled on side, standard pattern
MH21-600	FOMON Saber-Back Scissors, 5" (12.7 cm), slightly curved, semi-sharp outer edges used for dissection
MH21-602	FOMON Lower Lateral Scissors, 5" (12.7 cm), strong curve, dissecting blades
MH21-604	FOMON Dorsal Scissors, 5-1/2" (14 cm), light pattern, angular shanks



MH21-610



MH21-630

MH21-632

REF	Description
MH21-610	COTTLE Dorsal Scissors, 6-1/2" (16.5 cm), heavy pattern with rounded blades, angular shanks
MH21-630	JOSEPH Scissors, 5-3/4" (14.6 cm), straight, sharp points
MH21-632	JOSEPH Scissors, 5-3/4" (14.6 cm), curved, sharp points



MH27-1000



MH30-2190  
MH30-2195

REF	Description
MH27-1000	MARTIN Cartilage Scissors, 8" (20.3 cm), curved, serrated non-slip cutting edges
MH30-2190	BRAUN Episiotomy Scissors, 5-1/2" (14 cm), angled on side, guarded lower blade
MH30-2195	BRAUN Episiotomy Scissors, 8-1/2" (21.6 cm), angled to side, guarded lower blade



MH9-98

MH9-100

MH9-102  
MH9-104

MH9-110

REF	Description
MH9-98	NEWS Suture Scissors, 5-1/2" (14 cm), angled on flat
MH9-100	SPENCER Stitch Scissors, 3-1/2" (8.9 cm), delicate
MH9-102	LITTAUER JR. Stitch Scissors, 4-1/2" (11.4 cm), light pattern
MH9-104	LITTAUER Stitch Scissors, 5-1/2" (14 cm), standard pattern
MH9-110	O'BRIEN Stitch Scissors, 3-3/4" (9.5 cm), angled, sharp points



MH9-124



MH18-1486

REF	Description
MH9-124	Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade
MH18-1486	WESTCOTT Stitch Scissors, 4-1/2" (11.4 cm), curved, sharp points



MH18-1652



MH18-1653



MH21-536

REF	Description
MH18-1652	Eye Suture (GRADLE) Scissors, 3-3/4" (9.5 cm), slightly curved, sharp points
MH18-1653	Eye Suture Scissors, 4" (10.2 cm), ribbon type, slightly curved, sharp points
MH21-536	LITTLER Suture Carrying Scissors, 4-5/8" (11.8 cm), with suture hole in blades, curved





MH5-512 to  
MH5-516



MH5-550



MH5-560

REF	Description
MH5-512	LISTER Bandage Scissors, 4-1/2" (11.4 cm)
MH5-514	LISTER Bandage Scissors, 5-1/2" (14 cm)
MH5-516	LISTER Bandage Scissors, 7-1/4" (18.4 cm)
MH5-550	LISTER Bandage Scissors, 8" (20.3 cm), with one large finger ring
MH5-560	KNOWLES Bandage Scissors, 5-1/2" (14 cm), straight



MH5-568



MH27-3070

REF	Description
MH5-568	ESMARCH Heavy Duty Bandage and Cast Shears, 8" (20.3 cm)
MH27-3070	BRUNS Plaster Shears, 9-1/4" (23.5 cm), serrated blade



MH5-120TC  
MH5-124TC  
MH5-128TC



MH5-122TC  
MH5-126TC  
MH5-130TC

REF	Description
MH5-120TC	MAYO Dissecting Scissors, 5-1/2" (14 cm), straight, standard beveled blades, Tungsten Carbide blades
MH5-122TC	MAYO Dissecting Scissors, 5-1/2" (14 cm), curved, standard beveled blades, Tungsten Carbide blades
MH5-124TC	MAYO Dissecting Scissors, 6-3/4" (17.1 cm), straight, standard beveled blades, Tungsten Carbide blades
MH5-126TC	MAYO Dissecting Scissors, 6-3/4" (17.1 cm), curved, standard beveled blades, Tungsten Carbide blades
MH5-128TC	MAYO Dissecting Scissors, 9" (22.9 cm), straight, standard beveled blades, Tungsten Carbide blades
MH5-130TC	MAYO Dissecting Scissors, 9" (22.9 cm), curved, standard beveled blades, Tungsten Carbide blades

\*For easy identification, MeisterHand Tungsten Carbide scissors have two gold-plated finger rings.



MH5-179TC

MH5-180TC  
MH5-182TC

MH5-182A-TC

REF	Description
MH5-179TC	METZENBAUM Scissors, Standard Pattern, 5-1/2" (14 cm), straight, blunt points, Tungsten Carbide blades
MH5-180TC	METZENBAUM Scissors, Standard Pattern, 5-1/2" (14 cm), curved, blunt points, Tungsten Carbide blades
MH5-182TC	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), curved, blunt points, Tungsten Carbide blades
MH5-182A-TC	METZENBAUM Scissors, Delicate Pattern, 7" (17.8 cm), curved, blunt points, Tungsten Carbide blades

\*For easy identification, MeisterHand Tungsten Carbide scissors have two gold-plated finger rings.



MH5-304TC  
MH5-306TC



MH5-312TC  
MH5-314TC



MH9-124TC

REF	Description
MH5-304TC	Iris Scissors, 4-1/2" (11.4 cm), straight, Tungsten Carbide blades
MH5-306TC	Iris Scissors, 4-1/2" (11.4 cm), curved, Tungsten Carbide blades
MH5-312TC	Strabismus Scissors, 4-1/2" (11.4 cm), straight, Tungsten Carbide blades
MH5-314TC	Strabismus Scissors, 4-1/2" (11.4 cm), curved, Tungsten Carbide blades
MH9-124TC	Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade, Tungsten Carbide blades

\*For easy identification, MeisterHand Tungsten Carbide scissors have two gold-plated finger rings.



MH5-SC-120



MH5-SC-122



MH5-SC-124



MH5-SC-126

REF	Description
MH5-SC-120	SuperCut MAYO Scissors, 5-1/2" (14 cm), straight
MH5-SC-122	SuperCut MAYO Scissors, 5-1/2" (14 cm), curved
MH5-SC-124	SuperCut MAYO Scissors, 6-3/4" (17.1 cm), straight
MH5-SC-126	SuperCut MAYO Scissors, 6-3/4" (17.1 cm), curved

\*For easy identification, MeisterHand SuperCut scissors have two black finger rings.



MH5-SC-128



MH5-SC-130



MH5-SC-180



MH5-SC-182

REF	Description
MH5-SC-128	SuperCut MAYO Scissors, 9" (22.9 cm), straight
MH5-SC-130	SuperCut MAYO Scissors, 9" (22.9 cm), curved
MH5-SC-180	SuperCut METZENBAUM Scissors, 5-1/2" (14 cm), curved
MH5-SC-182	SuperCut METZENBAUM Scissors, 7" (17.8 cm), curved

\*For easy identification, MeisterHand SuperCut scissors have two black finger rings.



MH5-SC-304



MH5-SC-306


REF	Description
MH5-SC-304	SuperCut Iris Scissors, 4-1/2" (11.4 cm), straight
MH5-SC-306	SuperCut Iris Scissors, 4-1/2" (11.4 cm), curved

\*For easy identification, MeisterHand SuperCut scissors have two black finger rings.




**MH6-4 to MH6-14**
**MH6-40 to MH6-48**
**MH6-104**
**MH6-106**

REF	Description
MH6-4	Dressing Forceps, 4-1/2" (11.4 cm), serrated tips, standard pattern serrated handles
MH6-6	Dressing Forceps, 5" (12.7 cm), serrated tips, standard pattern serrated handles
MH6-8	Dressing Forceps, 5-1/2" (14 cm), serrated tips, standard pattern serrated handles
MH6-10	Dressing Forceps, 6" (15.2 cm), serrated tips, standard pattern serrated handles
MH6-14	Dressing Forceps, 8" (20.3 cm), serrated tips, standard pattern serrated handles
MH6-40	Tissue Forceps, 4-1/2" (11.4 cm), 1 X 2 teeth, standard pattern serrated handles
MH6-42	Tissue Forceps, 5" (12.7 cm), 1 X 2 teeth, standard pattern serrated handles
MH6-44	Tissue Forceps, 5-1/2" (14 cm), 1 X 2 teeth, standard pattern serrated handles
MH6-46	Tissue Forceps, 6" (15.2 cm), 1 X 2 teeth, standard pattern serrated handles
MH6-48	Tissue Forceps, 8" (20.3 cm), 1 X 2 teeth, standard pattern serrated handles
MH6-104	SEMKEN Dressing Forceps, 5" (12.7 cm), serrated
MH6-106	SEMKEN Tissue Forceps, 5" (12.7 cm), 1 X 2 teeth



MH6-118

MH6-120  
MH6-122

MH6-121

MH6-123

MH6-124

REF	Description
MH6-118	ADSON Dressing Forceps, 4-3/4" (12.1 cm), delicate, serrated
MH6-120	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, delicate
MH6-121	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, cross serrated tips, delicate
MH6-122	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 2 X 3 teeth, delicate
MH6-123	ADSON Tissue & Suture Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, with tying platform
MH6-124	BROWN-ADSON Tissue Forceps, 4-3/4" (12.1 cm), 7 X 7 side grasping teeth



MH17-2500  
MH17-2510



MH6-142  
MH6-144  
MH6-145



MH6-148



MH6-150

REF	Description
MH17-2500	Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, tips 0.5 mm wide
MH17-2510	Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated tips, 0.5 mm wide
MH6-142	RUSSIAN Tissue Forceps, 6" (15.2 cm)
MH6-144	RUSSIAN Tissue Forceps, 8" (20.3 cm)
MH6-145	RUSSIAN Tissue Forceps, 10" (25.4 cm)
MH6-148	BONNEY Tissue Forceps, 7" (17.8 cm), 1 X 2 teeth, serrated tips
MH6-150	BONNEY Tissue Forceps, 7" (17.8 cm), 2 X 3 teeth, serrated tips



MH6-154

MH6-160

MH6-180

MH6-182

MH6-214

REF	Description
MH6-154	POTTS-SMITH Dressing Forceps, 7" (17.8 cm), serrated
MH6-160	POTTS-SMITH Tissue Forceps, 7" (17.8 cm), 1 X 2 teeth, serrated tips
MH6-180	GERALD Dressing Forceps, 7" (17.8 cm), straight, serrated
MH6-182	GERALD Tissue Forceps, 7" (17.8 cm), 1 X 2 teeth, straight
MH6-214	SINGLEY Tissue Forceps, 9" (22.9 cm), fenestrated serrated jaws



MH6-190



MH6-192



MH6-196



MH6-197



MH6-198

REF	Description
MH6-190	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), bayonet, serrated, scraper end
MH6-192	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), 1 X 2 teeth, bayonet, scraper end
MH6-196	GRUENWALD (ADSON) Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated
MH6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps, 6-1/2" (16.5 cm), serrated
MH6-198	ADSON Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated, delicate tips



MH18-780



MH18-782



MH18-784



MH18-786



MH18-788



MH18-790

REF	Description
MH18-780	Eye Dressing Forceps, 4" (10.2 cm), straight, standard pattern, 0.8 mm wide serrated tips
MH18-782	Eye Dressing Forceps, 4" (10.2 cm), half curved, standard pattern, 0.8 mm wide serrated tips
MH18-784	Eye Dressing Forceps, 4" (10.2 cm), full curved, standard pattern, 0.8 mm wide serrated tips
MH18-786	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, straight, standard pattern, 0.8 mm wide tips
MH18-788	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, half curved, standard pattern, 0.8 mm wide tips
MH18-790	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, full curved, standard pattern, 0.8 mm wide tips



MH18-835



MH18-836



MH18-866



MH18-867



MH18-868

REF	Description
MH18-835	BONN Micro Suturing Forceps, 2-3/4" (7 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-836	BONN Micro Suturing Forceps, 3-3/4" (9.5 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-866	BISHOP-HARMON Dressing Forceps, 3-3/8" (8.6 cm), fine cross serrated tips 0.5 mm wide
MH18-867	BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth
MH18-868	BISHOP-HARMON Tissue Forceps, 3-3/8" (8.6 cm), standard 0.7 mm 1 X 2 teeth



MH18-948



MH18-949



MH18-954



MH18-955

REF	Description
MH18-948	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), straight, smooth tying platform 0.3 mm X 5 mm, no teeth
MH18-949	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), angled, smooth tying platform 0.3 mm X 5 mm, no teeth
MH18-954	CASTROVIEJO Suturing Forceps, 4" (10.2 cm), with 11 mm wide handles, 1 X 2 teeth 0.6 mm, with tying platform
MH18-955	CASTROVIEJO Suturing Forceps, 4" (10.2 cm), with 11 mm wide handles, 1 X 2 teeth 0.9 mm, with tying platform





MH19-350



MH19-370



MH26-956  
MH26-958

REF	Description
MH19-350	WILDE Dressing Forceps, 5" (12.7 cm), angular with serrated tips
MH19-370	LUCAE Dressing Forceps, 5-1/2" (14 cm), bayonet shape with serrated tips
MH26-956	FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 1 X 2 teeth
MH26-958	FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 2 X 3 teeth



MH16-4



MH16-6 to MH16-12

REF	Description
MH16-4	Baby ALLIS Tissue Forceps, 5-1/2" (14 cm), 4 X 5 teeth, delicate jaws 4 mm wide
MH16-6	ALLIS Tissue Forceps, 6" (15.2 cm), 4 X 5 teeth
MH16-8	ALLIS Tissue Forceps, 6" (15.2 cm), 5 X 6 teeth
MH16-10	ALLIS Tissue Forceps, 7-1/4" (18.4 cm), 5 X 6 teeth
MH16-12	ALLIS Tissue Forceps, 10" (25.4 cm), 5 X 6 teeth



MH16-30



MH16-31




MH16-38

REF	Description
MH16-30	ALLIS-ADAIR Tissue Forceps, 6-1/2" (16.5 cm), 9 X 10 teeth
MH16-31	BOYS-ALLIS Tissue Forceps, 6" (15.2 cm), 5 X 6 teeth, heavy pattern
MH16-38	THOMS-ALLIS Tissue Forceps, 8" (20.3 cm), 6 X 7 teeth



MH16-42 to MH16-48



MH16-51

REF	Description
MH16-42	Baby BABCOCK Intestinal Forceps, 5-1/2" (14 cm), delicate jaws 6 mm wide
MH16-44	BABCOCK Intestinal Forceps, 6-1/4" (15.9 cm), jaws 9 mm wide
MH16-46	BABCOCK Intestinal Forceps, 8-1/4" (21 cm), jaws 10 mm wide
MH16-48	BABCOCK Intestinal Forceps, 9-1/2" (24.1 cm), jaws 15 mm wide
MH16-51	ADAIR Tissue Holding Tenaculum Forceps, 6-1/2" (16.5 cm)



MH6-118TC



MH6-154TC



MH6-165TC



REF	Description
MH6-118TC	ADSON Dressing Forceps, 4-3/4" (12.1 cm), cross serrated tips, Tungsten Carbide jaws
MH6-154TC	POTTS-SMITH Dressing Forceps, 7" (17.8 cm), cross serrated tips, Tungsten Carbide jaws
MH6-165TC	EVANS Pelvic Dressing Forceps, 10" (25.4 cm), cross serrated tips, Tungsten Carbide jaws



MH6-300  
MH6-304



MH6-316  
MH6-318



MH6-400

REF	Description
MH6-300	Plain Splinter Forceps, 3-1/2" (8.9 cm)
MH6-304	Plain Splinter Forceps, 4-1/2" (11.4 cm)
MH6-316	FEILCHENFELD Splinter Forceps, 3" (7.6 cm)
MH6-318	FEILCHENFELD Splinter Forceps, 4-1/2" (11.4 cm)
MH6-400	Splinter Liberator, 5" (12.7 cm), trocar point in reversible handle




MH7-2  
MH7-4

MH7-8  
MH7-10

MH7-24  
MH7-26

MH7-36  
MH7-38

REF	Description
MH7-2	HALSTED Mosquito Forceps, 5" (12.7 cm), straight
MH7-4	HALSTED Mosquito Forceps, 5" (12.7 cm), curved
MH7-8	HALSTED Mosquito Forceps, 5" (12.7 cm), straight, extra delicate
MH7-10	HALSTED Mosquito Forceps, 5" (12.7 cm), curved, extra delicate
MH7-24	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), straight
MH7-26	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), curved
MH7-36	KELLY Forceps, 5-1/2" (14 cm), straight
MH7-38	KELLY Forceps, 5-1/2" (14 cm), curved



MH7-42

MH7-44

MH7-46

MH7-48

MH7-50

MH7-52

MH7-66

MH7-68

REF	Description
MH7-42	CRILE Forceps, 5-1/2" (14 cm), straight
MH7-44	CRILE Forceps, 5-1/2" (14 cm), curved
MH7-46	CRILE Forceps, 6-1/4" (15.9 cm), straight
MH7-48	CRILE Forceps, 6-1/4" (15.9 cm), curved
MH7-50	Baby CRILE Forceps, 5-1/2" (14 cm), extra delicate, straight
MH7-52	Baby CRILE Forceps, 5-1/2" (14 cm), extra delicate, curved
MH7-66	KOCHER Forceps, 5-1/2" (14 cm), 1 X 2 teeth, straight
MH7-68	KOCHER Forceps, 5-1/2" (14 cm), 1 X 2 teeth, curved






MH7-92  
MH7-94

MH7-102

MH7-106

MH7-118 to  
MH7-144

REF	Description
MH7-92	RANKIN Forceps, 6-1/4" (15.9 cm), straight
MH7-94	RANKIN Forceps, 6-1/4" (15.9 cm), curved
MH7-102	DANDY Forceps, 5-1/2" (14 cm), curved sideways
MH7-106	PRATT T-Shape Forceps, 6" (15.2 cm)
MH7-118	ROCHESTER-PEAN Forceps, 5-1/2" (14 cm), straight
MH7-120	ROCHESTER-PEAN Forceps, 6-1/4" (15.9 cm), straight
MH7-122	ROCHESTER-PEAN Forceps, 7-1/4" (18.4 cm), straight
MH7-124	ROCHESTER-PEAN Forceps, 8" (20.3 cm), straight
MH7-126	ROCHESTER-PEAN Forceps, 9" (22.9 cm), straight
MH7-136	ROCHESTER-PEAN Forceps, 5-1/2" (14 cm), curved
MH7-138	ROCHESTER-PEAN Forceps, 6-1/4" (15.9 cm), curved
MH7-140	ROCHESTER-PEAN Forceps, 7-1/4" (18.4 cm), curved
MH7-142	ROCHESTER-PEAN Forceps, 8" (20.3 cm), curved
MH7-144	ROCHESTER-PEAN Forceps, 9" (22.9 cm), curved



MH7-150 to  
MH7-164



MH7-170 to  
MH7-182

REF	Description
MH7-150	ROCHESTER-OCHSNER Forceps, 6-1/4" (15.9 cm), straight
MH7-152	ROCHESTER-OCHSNER Forceps, 7-1/4" (18.4 cm), straight
MH7-154	ROCHESTER-OCHSNER Forceps, 8" (20.3 cm), straight
MH7-160	ROCHESTER-OCHSNER Forceps, 6-1/4" (15.9 cm), curved
MH7-162	ROCHESTER-OCHSNER Forceps, 7-1/4" (18.4 cm), curved
MH7-164	ROCHESTER-OCHSNER Forceps, 8" (20.3 cm), curved
MH7-170	ROCHESTER-CARMALT Forceps, 6-1/4" (15.9 cm), straight
MH7-172	ROCHESTER-CARMALT Forceps, 6-1/4" (15.9 cm), curved
MH7-180	ROCHESTER-CARMALT Forceps, 8" (20.3 cm), straight
MH7-182	ROCHESTER-CARMALT Forceps, 8" (20.3 cm), curved



MH7-202

MH7-206

MH7-220

MH7-221

MH7-240

MH7-242


REF	Description
MH7-202	MIXTER Forceps, 7-1/4" (18.4 cm), fully curved
MH7-206	MIXTER Forceps, 9" (22.9 cm), fully curved
MH7-220	PENNINGTON Hemostatic Forceps, 6" (15.2 cm), triangular jaws
MH7-221	PENNINGTON Hemostatic Forceps, 8" (20.3 cm), triangular jaws
MH7-240	ADSON Hemostatic Forceps, 7-1/4" (18.4 cm), straight
MH7-242	ADSON Hemostatic Forceps, 7-1/4" (18.4 cm), curved



MH7-255 to  
MH7-259



MH14-48



MH16-160  
MH16-162

REF	Description
MH7-255	GEMINI-MIXTER Hemostatic Forceps, 5-1/2" (14 cm), delicate, full curved jaws
MH7-256	GEMINI-MIXTER Hemostatic Forceps, 7" (17.8 cm), delicate, full curved jaws
MH7-257	GEMINI-MIXTER Hemostatic Forceps, 8" (20.3 cm), delicate, full curved jaws
MH7-258	GEMINI-MIXTER Hemostatic Forceps, 9" (22.9 cm), delicate, full curved jaws
MH7-259	GEMINI-MIXTER Hemostatic Forceps, 11" (27.9 cm), delicate, full curved jaws
MH14-48	LAHEY Gall Duct Forceps, 7-1/2" (19.1 cm)
MH16-160	DOYEN Intestinal Forceps, 9" (22.9 cm), straight, flexible blades with longitudinal serrations
MH16-162	DOYEN Intestinal Forceps, 9" (22.9 cm), curved, flexible blades with longitudinal serrations



MH17-2600  
MH17-2602

MH22-662

MH22-668

MH25-800

REF	Description
MH17-2600	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), straight, extremely delicate
MH17-2602	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), curved, extremely delicate
MH22-662	SCHNIDT Tonsil Hemostatic Forceps, 7-1/2" (19.1 cm), slight curve, closed rings
MH22-668	SAWTELL Tonsil Hemostatic Forceps, 7-1/2" (19.1 cm), full curve, closed rings
MH25-800	SAROT Intra-Thoracic Artery Forceps, 10" (25.4 cm), with 2-1/2" (6.4 cm) long and slightly curved jaws



MH25-814



MH25-836

MH25-837

REF	Description
MH25-814	MIXTER Thoracic Forceps, 11" (27.9 cm), right angle jaws, longitudinal serrations
MH25-836	KANTROWITZ Thoracic Forceps, 7-1/2" (19.1 cm), delicate right angle jaws
MH25-837	KANTROWITZ Thoracic Forceps, 9-1/2" (24.1 cm), delicate right angle jaws




MH30-1025



MH30-1700

MH30-1710




MH30-1716

MH30-1718

REF	Description
MH30-1025	JACOBS Uterine Vulsellum Forceps, 8-1/2" (21.6 cm), straight, serrated jaws, 2 X 2 teeth
MH30-1700	HEANEY Hysterectomy Forceps, 8-1/4" (21 cm), heavy pattern, single tooth, curved
MH30-1710	HEANEY Hysterectomy Forceps, 8-1/4" (21 cm), heavy pattern, double tooth, curved
MH30-1716	HEANEY-BALLENTINE Hysterectomy Forceps, 8-1/4" (21 cm), longitudinal serrations, single tooth, curved
MH30-1718	BALLENTINE Hysterectomy Forceps, 8-1/2" (21.6 cm), longitudinal serrations, single tooth, straight



MH7-504  
MH7-506



MH7-514  
MH7-515



MH7-516

REF	Description
MH7-504	BACKHAUS Towel Clamp, 4" (10.2 cm)
MH7-506	BACKHAUS Towel Clamp, 5-1/4" (13.3 cm)
MH7-514	LORNA Non-Perforating Towel Clamp, 4" (10.2 cm)
MH7-515	LORNA Non-Perforating Towel Clamp, 5-1/4" (13.3 cm)
MH7-516	PEERS Non-Perforating Towel Clamp, 5-3/4" (14.6 cm)





MH7-600  
MH7-602  
MH7-604



MH7-610  
MH7-614



MH7-622

REF	Description
MH7-600	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), straight, serrated jaws
MH7-602	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), straight, smooth jaws
MH7-604	FOERSTER Sponge Forceps, 7" (17.8 cm), straight, serrated jaws
MH7-610	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), curved, serrated jaws
MH7-614	FOERSTER Sponge Forceps, 7" (17.8 cm), curved, serrated jaws
MH7-622	BOZEMAN Sponge Forceps, 10-1/2" (26.7 cm), double curved



MH8-2

MH8-4

MH8-6

MH8-8

REF	Description
MH8-2	COLLIER Needle Holder, 5" (12.7 cm), fenestrated jaws
MH8-4	DERF Needle Holder, 4-3/4" (12.1 cm)
MH8-6	WEBSTER Needle Holder, 5" (12.7 cm), delicate, smooth jaws, standard pattern
MH8-8	HALSEY Needle Holder, 5" (12.7 cm), smooth jaws



MH8-14

MH8-15

MH8-20

MH8-22

MH8-16

MH8-17

REF	Description
MH8-14	OLSEN-HEGAR Needle Holder with Suture Scissors, 4-3/4" (12.1 cm), serrated jaws, extra delicate pattern
MH8-15	OLSEN-HEGAR Needle Holder with Suture Scissors, 5-1/2" (14 cm), serrated jaws
MH8-16	OLSEN-HEGAR Needle Holder with Suture Scissors, 6-1/2" (16.5 cm), serrated jaws
MH8-17	OLSEN-HEGAR Needle Holder with Suture Scissors, 7-1/4" (18.4 cm), serrated jaws
MH8-20	CONVERSE Needle Holder, 4-1/4" (10.8 cm), longitudinal serrations
MH8-22	BROWN Needle Holder, 5-1/4" (13.3 cm), convex jaws, fine serrations



MH8-40



MH8-42 to MH8-48



MH8-50 to MH8-52

REF	Description
MH8-40	BAUMGARTNER Needle Holder, 5-1/2" (14 cm)
MH8-42	MAYO-HEGAR Needle Holder, 5" (12.7 cm)
MH8-44	MAYO-HEGAR Needle Holder, 6-1/4" (15.9 cm)
MH8-46	MAYO-HEGAR Needle Holder, 7-1/4" (18.4 cm)
MH8-48	MAYO-HEGAR Needle Holder, 8" (20.3 cm)
MH8-50	CRILE-WOOD Needle Holder, 6" (15.2 cm)
MH8-51	CRILE-WOOD Needle Holder, 7" (17.8 cm)
MH8-52	CRILE-WOOD Needle Holder, 8" (20.3 cm)



MH8-100



MH8-130



MH18-1828

MH18-1832

REF	Description
MH8-100	HEANEY Needle Holder, 8-1/2" (21.6 cm), curved jaws
MH8-130	MAYO-HEGAR (MASSON) Needle Holder, 10-1/2" (26.7 cm)
MH18-1828	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), straight, with lock
MH18-1832	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), curved, with lock



MH8-4TC

MH8-6TC

MH8-8A-TC

REF	Description
MH8-4TC	DERF Needle Holder, 4-3/4" (12.1 cm), serrated jaws, Tungsten Carbide jaws
MH8-6TC	WEBSTER Needle Holder, 5-1/4" (13.3 cm), with smooth jaws, standard pattern, delicate, Tungsten Carbide jaws
MH8-8A-TC	HALSEY Needle Holder, 5" (12.7 cm) serrated jaws, Tungsten Carbide jaws

\*For easy identification, MeisterHand® Tungsten Carbide Needle Holders have two gold-plated finger rings.



MH8-15TC  
MH8-17TC




MH8-16TC



MH8-40TC

REF	Description
MH8-15TC	OLSEN-HEGAR Needle Holder with Suture Scissors, 5-1/2" (14 cm), serrated, Tungsten Carbide jaws
MH8-16TC	OLSEN-HEGAR Needle Holder with Suture Scissors, 6-1/2" (16.5 cm), serrated, Tungsten Carbide jaws
MH8-17TC	OLSEN-HEGAR Needle Holder with Suture Scissors, 7-1/2" (19.1 cm), serrated, Tungsten Carbide jaws
MH8-40TC	BAUMGARTNER Needle Holder, 5-1/2" (14 cm), serrated, Tungsten Carbide jaws

\*For easy identification, MeisterHand® Tungsten Carbide Needle Holders have two gold-plated finger rings.



MH8-42TC to  
MH8-48TC



MH8-50TC to  
MH8-52TC



MH8-100TC


REF	Description
MH8-42TC	MAYO-HEGAR Needle Holder, 5" (12.7 cm), serrated, Tungsten Carbide jaws
MH8-44TC	MAYO-HEGAR Needle Holder, 6-1/4" (15.9 cm), serrated, Tungsten Carbide jaws
MH8-46TC	MAYO-HEGAR Needle Holder, 7-1/4" (18.4 cm), serrated, Tungsten Carbide jaws
MH8-48TC	MAYO-HEGAR Needle Holder, 8" (20.3 cm), serrated, Tungsten Carbide jaws
MH8-50TC	CRILE-WOOD Needle Holder, 6" (15.2 cm), serrated, Tungsten Carbide jaws
MH8-51TC	CRILE-WOOD Needle Holder, 7" (17.8 cm), serrated, Tungsten Carbide jaws
MH8-52TC	CRILE-WOOD Needle Holder, 8" (20.3 cm), serrated, Tungsten Carbide jaws
MH8-100TC	HEANEY Needle Holder, 8-1/2" (21.6 cm), curved, serrated, Tungsten Carbide jaws

\*For easy identification, MeisterHand® Tungsten Carbide Needle Holders have two gold-plated finger rings.





MH8-101TC to  
MH8-105TC



MH8-130TC



MH18-1828TC

REF	Description
MH8-101TC	RYDER Needle Holder, 5" (12.7 cm), serrated, Tungsten Carbide jaws
MH8-103TC	RYDER Needle Holder, 6" (15.2 cm), serrated, Tungsten Carbide jaws
MH8-104TC	RYDER Needle Holder, 7" (17.8 cm), serrated, Tungsten Carbide jaws
MH8-105TC	RYDER Needle Holder, 8" (20.3 cm), serrated, Tungsten Carbide jaws
MH8-130TC	MAYO-HEGAR (MASSON) Needle Holder, 10-1/2" (26.7 cm), serrated, Tungsten Carbide jaws
MH18-1828TC	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), smooth jaws, straight, with lock, Tungsten Carbide jaws

\*For easy identification, MeisterHand® Tungsten Carbide Needle Holders have two gold-plated finger rings.



MH2-104SS  
MH2-109



MH2-110SS



MH10-310



MH10-312

REF	Description
MH2-104SS	YANKAUER Suction Tube, 11-3/4" (29.8 cm), with removable tip and tubing connection
MH2-106SS	Tip only for YANKAUER Suction Tube (MH2-104SS)
MH2-109	YANKAUER Pediatric Suction Tube, 8-1/2" (21.6 cm), with removable tip, delicate pattern
MH2-110SS	ANDREWS-PYNCHON Suction Tube 10" (25.4 cm), delicate pattern
MH10-310	POOLE Suction Tube, 23 French, (7.6 mm), curved
MH10-312	POOLE Suction Tube, 30 French, (9.9 mm), straight



MH19-570 to  
MH19-577



MH19-580 to  
MH19-584

REF	Description
MH19-570	FRAZIER Suction Tube, 6 French, (2 mm), with finger cut-off, angled
MH19-572	FRAZIER Suction Tube, 7 French, (2.3 mm), with finger cut-off, angled
MH19-573	FRAZIER Suction Tube, 8 French, (2.6 mm), with finger cut-off, angled
MH19-574	FRAZIER Suction Tube, 9 French, (3 mm), with finger cut-off, angled
MH19-575	FRAZIER Suction Tube, 10 French, (3.3 mm), with finger cut-off, angled
MH19-576	FRAZIER Suction Tube, 11 French, (3.6 mm), with finger cut-off, angled
MH19-577	FRAZIER Suction Tube, 12 French, (4 mm), with finger cut-off, angled
MH19-580	BARON Suction Tube, 3 French (1 mm), with finger cut-off, working length 7.5 cm
MH19-582	BARON Suction Tube, 5 French (1.7 mm), with finger cut-off, working length 7.5 cm
MH19-584	BARON Suction Tube, 7 French (2.3 mm), with finger cut-off, working length 7.5 cm



MH11-2

MH11-12  
MH11-13

MH11-14



MH11-16


REF	Description
MH11-2	JANSEN Retractor, 4-1/4" (10.8 cm), spread 2" (5.1 cm), 3 X 3 blunt prongs 10 mm deep, standard model
MH11-12	ALM Retractor, 2-3/4" (7 cm), spread 2-1/4" (5.7 cm), prongs 1/8" (.32 cm) deep
MH11-13	ALM Retractor, 3-7/8" (9.8 cm), spread 3" (7.6 cm), prongs 1/2" (1.3 cm) deep
MH11-14	HEISS Self-Retaining Cross Action Retractor, 4" (10.2 cm), sharp prongs
MH11-16	HEISS Self-Retaining Cross Action Retractor, 4" (10.2 cm), blunt prongs



MH11-22



MH11-73

MH11-74  
MH11-76

MH11-78


REF	Description
MH11-22	GREEN Retractor (Thyroid), 8-1/2" (21.6 cm)
MH11-73	RAGNELL Retractor, 6" (15.2 cm), double end 3 X 8 mm and 5 X 15 mm blades, delicate
MH11-74	SENN Retractor, 6-3/8" (16.2 cm), sharp, double end
MH11-76	SENN Retractor, 6-3/8" (16.2 cm), blunt, double end
MH11-78	CRILE Retractor, 4-1/2" (11.4 cm), double end 3 X 11 mm and 6 X 20 mm blades



MH11-102



MH11-110



MH11-114



MH11-120

REF	Description
MH11-102	CUSHING Vein Retractor, 9" (22.9 cm), blade 11 X 13 mm
MH11-110	FARABEUZ Retractors, 6" (15.2 cm), double end, set of 2
MH11-114	PARKER Retractors, 5" (12.7 cm), double end, set of 2
MH11-120	GOELET Retractor, 7-1/2" (19.1 cm), double end



MH11-122



MH11-132 to  
MH11-146



MH11-162 to  
MH11-166

REF	Description
MH11-122	U.S. ARMY Retractors, 8-1/2" (21.6 cm), double end, set of 2
MH11-132	Ribbon Retractor, 1" (2.5 cm) X 13" (33 cm), malleable
MH11-136	Ribbon Retractor, 1-1/2" (3.8 cm) X 13" (33 cm), malleable
MH11-140	Ribbon Retractor, 2" (5.1 cm) X 13" (33 cm), malleable
MH11-144	Ribbon Retractor, Pediatric, 5/8" (1.6 cm) X 8" (20.3 cm), malleable
MH11-146	Ribbon Retractor, Pediatric, 1" (2.5 cm) X 8" (20.3 cm), malleable
MH11-162	LANGENBECK Retractor, 8-1/2" (21.6 cm), blade 3/8" (1 cm) X 1-1/2" (3.8 cm)
MH11-164	LANGENBECK Retractor, 8-1/2" (21.6 cm), blade 1/2" (1.3 cm) X 1-1/2" (3.8 cm)
MH11-166	LANGENBECK Retractor, 8-1/2" (21.6 cm), blade 5/8" (1.6 cm) X 1-1/2" (3.8 cm)

MH11-180 to  
MH11-184

MH11-185



MH11-202

MH11-219  
MH11-220MH11-206 to  
MH11-210

REF	Description
MH11-180	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1" (2.5 cm) X 1-3/4" (4.4 cm)
MH11-182	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1-1/8" (2.8 cm) X 2-7/8" (7.3 cm)
MH11-184	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1-1/8" (2.8 cm) X 3-1/4" (8.3 cm)
MH11-185	HEANEY-SIMON Retractor, 11-1/2" (29.2 cm), blade 1" (2.5 cm) X 4-1/2" (11.4 cm)
MH11-202	DEAVER Retractor, 1" (2.5 cm) X 12" (30.5 cm)
MH11-206	DEAVER Retractor, 1-1/2" (3.8 cm) X 12" (30.5 cm)
MH11-208	DEAVER Retractor, 2" (5.1 cm) X 12" (30.5 cm)
MH11-210	DEAVER Retractor, 3" (7.6 cm) X 12" (30.5 cm)
MH11-219	DEAVER Retractor, Pediatric, 3/4" (1.9 cm) X 7" (17.8 cm)
MH11-220	DEAVER Retractor, Pediatric, 7/8" (2.2 cm) X 8" (20.3 cm)





MH11-240 to  
MH11-244



MH11-246 to  
MH11-249

REF	Description
MH11-240	MEYERDING Retractor, 9-1/2" (24.1 cm), large blade, 3-1/2" (8.9 cm) X 2" (5.1 cm)
MH11-242	MEYERDING Retractor, 9-1/2" (24.1 cm), medium blade, 3" (7.6 cm) X 1" (2.5 cm)
MH11-244	MEYERDING Retractor, 9-1/2" (24.1 cm), small blade, 2" (5.1 cm) X 5/8" (1.6 cm)
MH11-246	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1" (2.5 cm) X 3/4" (1.9 cm)
MH11-247	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1-1/4" (3.2 cm) X 1" (2.5 cm)
MH11-248	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1-1/2" (3.8 cm) X 1-1/2" (3.8 cm)
MH11-249	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 2" (5.1 cm) X 3/4" (1.9 cm)



MH11-260 to  
MH11-260-B




MH11-266 to  
MH11-269



MH11-314  
MH11-316

REF	Description
MH11-260	RICHARDSON-EASTMAN Retractors, double end, set of 2
MH11-260-A	RICHARDSON-EASTMAN Retractor, double end, 11" (27.9 cm), large size
MH11-260-B	RICHARDSON-EASTMAN Retractor, double end, 10" (25.4 cm), small size
MH11-266	KELLY Retractor, 9-3/4" (24.8 cm), loop handle, 2" (5.1 cm) X 1-1/2" (3.8 cm)
MH11-268	KELLY Retractor, 9-3/4" (24.8 cm), loop handle, 3" (7.6 cm) X 2-1/2" (6.4 cm)
MH11-269	KELLY Retractor, 9-3/4" (24.8 cm), loop handle, 3-1/2" (8.9 cm) X 2-3/4" (7 cm)
MH11-314	HARRINGTON Retractor, 12-1/2" (31.8 cm), blade 2-1/2" (6.4 cm) wide
MH11-316	HARRINGTON Retractor, 12-1/2" (31.8 cm), blade 1-1/2" (3.8 cm) wide



MH11-530 to  
MH11-544



MH11-580  
MH11-582



MH11-590

REF	Description
MH11-530	VOLKMAN Retractor, 8-1/2" (21.6 cm), 2 prongs, sharp, ring handle
MH11-532	VOLKMAN Retractor, 8-1/2" (21.6 cm), 3 prongs, sharp, ring handle
MH11-534	VOLKMAN Retractor, 8-1/2" (21.6 cm), 4 prongs, sharp, ring handle
MH11-536	VOLKMAN Retractor, 8-1/2" (21.6 cm), 6 prongs, sharp, ring handle
MH11-538	VOLKMAN Retractor, 8-1/2" (21.6 cm), 2 prongs, blunt, ring handle
MH11-540	VOLKMAN Retractor, 8-1/2" (21.6 cm), 3 prongs, blunt, ring handle
MH11-542	VOLKMAN Retractor, 8-1/2" (21.6 cm), 4 prongs, blunt, ring handle
MH11-544	VOLKMAN Retractor, 8-1/2" (21.6 cm), 6 prongs, blunt, ring handle
MH11-580	BENNETT Tibia Retractor, 9-1/2" (24.1 cm), small blade 1-3/4" (4.4 cm) wide
MH11-582	BENNETT Tibia Retractor, 9-1/2" (24.1 cm), large blade 2-1/2" (6.4 cm) wide
MH11-590	HIBBS Retractor, 9-1/4" (23.5 cm), blade 1" (2.5 cm) X 3" (7.6 cm) deep



MH11-364 to  
MH11-366

MH11-600  
MH11-602

MH11-608  
MH11-610

MH11-614  
MH11-616

REF	Description
MH11-364	GELPI Retractor, 3-1/2" (8.9 cm), with grip lock, mini
MH11-365	GELPI Retractor, Pediatric, 5-1/2" (14 cm), with grip lock
MH11-366	GELPI Retractor, 7-1/2" (19.1 cm), with grip lock
MH11-600	WEITLANER Retractor, 4" (10.2 cm), sharp, 2 X 3 teeth
MH11-602	WEITLANER Retractor, 4" (10.2 cm), blunt, 2 X 3 teeth
MH11-608	WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth
MH11-610	WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth
MH11-614	WEITLANER Retractor, 6-1/2" (16.5 cm), sharp, 3 X 4 teeth
MH11-616	WEITLANER Retractor, 6-1/2" (16.5 cm), blunt, 3 X 4 teeth



**MH11-630-BL**  
**MH11-630-SH**




**MH11-684**



**MH11-700**  
**MH11-702**


REF	Description
MH11-630-BL	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth, with hinged blades
MH11-630-SH	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth, with hinged blades
MH11-684	BECKMAN-ADSON Retractor (Thyroid), 12-1/2" (31.8 cm), hinged, sharp, 4 X 4 teeth, 3/4" (1.9 cm) X 1" (2.5 cm)
MH11-700	ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth 3/4" (1.9 cm) X 3/4" (1.9 cm), straight arms
MH11-702	ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth 3/4" (1.9 cm) X 3/4" (1.9 cm), angled arms



MH18-100  
MH18-102



MH18-110 to  
MH18-116



MH21-145  
MH21-146

MH21-148

REF	Description
MH18-100	KNAPP (ROLLET) Retractor, 5-1/4" (13.3 cm), with 4 sharp prongs
MH18-102	KNAPP Retractor, 5-1/4" (13.3 cm), with 4 blunt prongs
MH18-110	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 1, 11 mm
MH18-112	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 2, 13 mm
MH18-114	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 3, 15 mm
MH18-116	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 4, 17 mm
MH21-145	COTTLE Retractor, 5-1/2" (14 cm), right ball tip
MH21-146	COTTLE Retractor, 5-1/2" (14 cm), left ball tip
MH21-148	FOMON Retractor, 6-1/2" (16.5 cm), two prongs with ball tips, 11 mm wide



MH21-165



MH21-188



MH22-854



MH23-1049

REF	Description
MH21-165	COTTLE Knife Guide & Retractor, 8" (20.3 cm), double ended, ball tip prongs 10 mm wide, and 5 mm wide slotted hook
MH21-188	AUFRICHT Retractor-Speculum, 6-1/2" (16.5 cm), solid blade, 45 mm long
MH22-854	HURD Dissector and Pillar Retractor, 9" (22.9 cm), standard pattern, 10.2 mm X 13.9 mm blade
MH23-1049	JACKSON Trachea Retractor, 7" (17.8 cm), double ended, 2 blunt prongs 5/16" (2.4 cm) wide, and solid blade 3/8" (1 cm) wide



MH27-213  
MH27-214



MH27-215 to  
MH27-219

REF	Description
MH27-213	MINI-HOHMANN Retractor, 6-1/4" (15.9 cm), blade 6 mm wide
MH27-214	MINI-HOHMANN Retractor, 6-1/4" (15.9 cm), blade 8 mm wide
MH27-215	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 10 mm wide
MH27-216	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 18 mm wide
MH27-217	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 25 mm wide
MH27-218	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 43 mm wide
MH27-219	HOHMANN Retractor, 9-3/4" (24.8 cm), blade 70 mm wide





MH28-68



MH28-69



MH30-395  
MH30-406  
MH30-415

REF	Description
MH28-68	HILL-FERGUSON Retractor, 8-1/2", small blade 2.2 X 6.4 cm
MH28-69	HILL-FERGUSON Retractor, 8-3/4", large blade 3 X 7.6 cm
MH30-395	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 1-3/4" (4.4 cm) X 2-1/2" (6.4 cm)
MH30-406	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 2-1/8" (5.4 cm) X 3" (7.6 cm)
MH30-415	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 2-1/4" (5.7 cm) X 3-1/2" (8.9 cm)



MH18-412  
MH18-414

MH18-416

MH18-440

MH18-446

REF	Description
MH18-412	TYRELL Hook, 5" (12.7 cm), sharp
MH18-414	TYRELL Hook, 5" (12.7 cm), blunt
MH18-416	GUTHRIE Fixation Hook, 5" (12.7 cm), sharp double prongs 1.5 mm wide
MH18-440	STEVENS Tenotomy Hook, 5" (12.7 cm), right angle 5 mm deep
MH18-446	JAMESON Strabismus Hook, 5" (12.7 cm), 10 mm deep



MH18-450  
MH18-452  
MH18-454




MH19-220  
MH19-222



MH21-88  
MH21-89

REF	Description
MH18-450	VON GRAEFE Strabismus Hook, 5-1/4" (13.3 cm), small size 8 mm deep
MH18-452	VON GRAEFE Strabismus Hook, 5-1/2" (14 cm), medium size 10 mm deep
MH18-454	VON GRAEFE Strabismus Hook, 5-3/4" (14.6 cm), large size 11 mm deep
MH19-220	DAY Ear Hook, 6-1/2" (16.5 cm), small size
MH19-222	DAY Ear Hook, 6-1/2" (16.5 cm), large size
MH21-88	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), large
MH21-89	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), small



MH21-90  
MH21-92  
MH21-94

MH21-144


MH21-145  
MH21-146

MH21-152

REF	Description
MH21-90	COTTLE Skin Hook, 5-1/2" (14 cm), small, shallow curve
MH21-92	COTTLE Skin Hook, 5-1/2" (14 cm), small, deep curve
MH21-94	COTTLE Skin Hook, 5-1/2" (14 cm), large, deep curve
MH21-144	COTTLE Double Hook, 5-1/2" (14 cm), two sharp prongs, 12 mm wide
MH21-145	COTTLE Double Hook, 5-1/2" (14 cm), right ball tip
MH21-146	COTTLE Double Hook, 5-1/2" (14 cm), left ball tip
MH21-152	Skin Hook, 4-3/4" (12.1 cm), sharp prong, 2 mm diameter



MH21-153



MH21-154 to  
MH21-160



MH23-1050



MH23-1051

REF	Description
MH21-153	JOSEPH Single Hook, 6-1/4" (15.9 cm), one sharp prong
MH21-154	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 2 mm wide
MH21-156	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 5 mm wide
MH21-158	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 7 mm wide
MH21-160	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 10 mm wide
MH23-1050	JACKSON Tenaculum Hook and Retractor, 5-3/4" (14.6 cm), sharp point
MH23-1051	NEWS Hook, 6" (15.2 cm), single sharp prong



MH26-1084



MH26-1090

REF	Description
MH26-1084	FRAZIER Dura Hook, 5" (12.7 cm), sharp point
MH26-1090	ADSON Dura and Dissecting Hook, 8" (20.3 cm), sharp blade 5 mm at right angle



MH17-2500  
MH17-2510



MH17-2600  
MH17-2602

REF	Description
MH17-2500	Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, tips 0.5 mm wide
MH17-2510	Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated tips, 0.5 mm wide
MH17-2600	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), straight, extremely delicate
MH17-2602	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), curved, extremely delicate



MH18-835



MH18-836



MH18-948



MH18-949

REF	Description
MH18-835	BONN Micro Suturing Forceps, 2-3/4" (7 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-836	BONN Micro Suturing Forceps, 3-3/4" (9.5 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-948	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), straight, smooth tying platform 0.3 mm X 5 mm, no teeth
MH18-949	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), angled, smooth tying platform 0.3 mm X 5 mm, no teeth





MH18-867



MH19-2081



MH19-2110

REF	Description
MH18-867	BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth
MH19-2081	Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, serrated jaws, 0.6 mm wide
MH19-2110	Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, oval cup jaws, 0.8 mm wide



MH7-2

MH7-4

MH7-8

MH7-10

MH7-24

MH7-26

MH17-2600

MH17-2602

REF	Description
MH7-2	HALSTED Mosquito Forceps, 5" (12.7 cm), straight
MH7-4	HALSTED Mosquito Forceps, 5" (12.7 cm), curved
MH7-8	HALSTED Mosquito Forceps, 5" (12.7 cm), straight, extra delicate
MH7-10	HALSTED Mosquito Forceps, 5" (12.7 cm), curved, extra delicate
MH7-24	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), straight
MH7-26	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), curved
MH17-2600	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), straight, extremely delicate
MH17-2602	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), curved, extremely delicate




MH18-1968



MH18-35



MH18-36



MH18-100  
MH18-102



MH18-110 to  
MH18-116

REF	Description
MH18-1968	FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
MH18-35	BARRAQUER Eye Speculum, 1-3/4" (4.4 cm), large blades 15 mm wide
MH18-36	BARRAQUER Eye Speculum, 1-3/4" (4.4 cm), small blades 10 mm wide
MH18-100	KNAPP (ROLLET) Retractor, 5-1/4" (13.3 cm), with 4 sharp prongs
MH18-102	KNAPP Retractor, 5-1/4" (13.3 cm), with 4 blunt prongs
MH18-110	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 1, 11 mm
MH18-112	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 2, 13 mm
MH18-114	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 3, 15 mm
MH18-116	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 4, 17 mm



MH18-412  
MH18-414



MH18-416



MH18-440



MH18-446



MH21-152

REF	Description
MH18-412	TYRELL Hook, 5" (12.7 cm), sharp
MH18-414	TYRELL Hook, 5" (12.7 cm), blunt
MH18-416	GUTHRIE Fixation Hook, sharp double prongs 1.5 mm wide
MH18-440	STEVENS Tenotomy Hook, 5" (12.7 cm), right angle 5 mm deep
MH18-446	JAMESON Strabismus Hook, 5" (12.7 cm), 10 mm deep
MH21-152	Skin Hook, 4-3/4" (12.1 cm), sharp prong, 2 mm diameter



MH18-450 to  
MH18-454

MH18-740 to  
MH18-750

MH18-780

MH18-782

MH18-784

REF	Description
MH18-450	VON GRAEFE Strabismus Hook, 5-1/4" (13.3 cm), small size 8 mm deep
MH18-452	VON GRAEFE Strabismus Hook, 5-1/2" (14 cm), medium size 10 mm deep
MH18-454	VON GRAEFE Strabismus Hook, 5-3/4" (14.6 cm), large size 11 mm deep
MH18-740	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 0000-000
MH18-742	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 00-0
MH18-744	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 1-2
MH18-746	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 3-4
MH18-748	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 5-6
MH18-750	BOWMAN Lacrimal Probe, stainless, 5" (12.7 cm), size 7-8
MH18-780	Eye Dressing Forceps, 4" (10.2 cm), straight, standard pattern, 0.8 mm wide serrated tips
MH18-782	Eye Dressing Forceps, 4" (10.2 cm), half curved, standard pattern, 0.8 mm wide serrated tips
MH18-784	Eye Dressing Forceps, 4" (10.2 cm), full curved, standard pattern, 0.8 mm wide serrated tips

↔



MH18-786

↔




MH18-788

↔



MH18-790

REF	Description
MH18-786	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, straight, standard pattern, 0.8 mm wide tips
MH18-788	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, half curved, standard pattern, 0.8 mm wide tips
MH18-790	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, full curved, standard pattern, 0.8 mm wide tips



MH18-835

MH18-836

MH18-866

MH18-867

MH18-868

REF	Description
MH18-835	BONN Micro Suturing Forceps, 2-3/4" (7 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-836	BONN Micro Suturing Forceps, 3-3/4" (9.5 cm), 1 X 2 teeth, 0.12 mm, with tying platform
MH18-866	BISHOP-HARMON Dressing Forceps, 3-3/8" (8.6 cm), fine cross serrated tips 0.5 mm wide
MH18-867	BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth
MH18-868	BISHOP-HARMON Tissue Forceps, 3-3/8" (8.6 cm), standard 0.7 mm 1 X 2 teeth



MH18-948



MH18-949



MH18-954



MH18-955

REF	Description
MH18-948	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), straight, smooth tying platform 0.3 mm X 5 mm, no teeth
MH18-949	McPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), angled, smooth tying platform 0.3 mm X 5 mm, no teeth
MH18-954	CASTROVIEJO Suturing Forceps, 4" (10.2 cm), with 11 mm wide handles, 1 X 2 teeth 0.6 mm, with tying platform
MH18-955	CASTROVIEJO Suturing Forceps, 4" (10.2 cm), with 11 mm wide handles, 1 X 2 teeth 0.9 mm, with tying platform





**MH5-300**  
**MH5-304**

**MH5-302**  
**MH5-306**

**MH5-312**  
**MH5-314**

**MH18-1396**  
**MH18-1398**

REF	Description
MH5-300	Iris Scissors, 4" (10.2 cm), straight
MH5-302	Iris Scissors, 4" (10.2 cm), curved
MH5-304	Iris Scissors, 4-1/2" (11.4 cm), straight
MH5-306	Iris Scissors, 4-1/2" (11.4 cm), curved
MH5-312	Strabismus Scissors, 4" (10.2 cm), straight
MH5-314	Strabismus Scissors, 4" (10.2 cm), curved
MH18-1396	Iris Scissors, 3-1/2" (8.9 cm), straight, with 20 mm blades, delicate
MH18-1398	Iris Scissors, 3-1/2" (8.9 cm), curved, with 20 mm blades, delicate



MH18-1415

MH18-1416

MH18-1417

MH18-1418

REF	Description
MH18-1415	Iris Scissors, 4" (10.2 cm), straight, standard pattern, ribbon-type
MH18-1416	Iris Scissors, 4" (10.2 cm), curved, standard pattern, ribbon-type
MH18-1417	Iris Scissors, 3-3/4" (9.5 cm), straight, miniature pattern with 14 mm blades, ribbon-type
MH18-1418	Iris Scissors, 3-3/4" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon-type



MH18-1420



MH18-1422



MH18-1426



MH18-1428

REF	Description
MH18-1420	KNAPP Iris Scissors, 4" (10.2 cm), straight, sharp-blunt points
MH18-1422	KNAPP Iris Scissors, 4" (10.2 cm), straight, blunt-blunt points
MH18-1426	KNAPP Iris Scissors, 4" (10.2 cm), curved, sharp-blunt points
MH18-1428	KNAPP Iris Scissors, 4" (10.2 cm), curved, blunt-blunt points



MH18-1460



MH18-1462 to MH18-1466



MH18-1472 to MH18-1476

### STEVENS TENOTOMY SCISSORS

REF	Description
MH18-1460	4-1/8" (10.5 cm), straight, short blades, sharp points
MH18-1462	4-1/8" (10.5 cm), straight, short blades, blunt points
MH18-1464	4-1/8" (10.5 cm), curved, short blades, sharp points
MH18-1466	4-1/8" (10.5 cm), curved, short blades, blunt points
MH18-1472	4-1/2" (11.4 cm), straight, long blades, blunt points
MH18-1476	4-1/2" (11.4 cm), curved, long blades, blunt points



MH18-1572



MH18-1576



MH18-1578



MH18-1620



MH18-1622

REF	Description
MH18-1572	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), angled blades 11 mm, blunt points
MH18-1576	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), curved, sharp points
MH18-1578	CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), curved, blunt points
MH18-1620	VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), straight, extra delicate sharp points
MH18-1622	VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), curved, extra delicate sharp points



MH9-100



MH9-110

REF	Description
MH9-100	SPENCER Stitch Scissors, 3-1/2" (8.9 cm), delicate
MH9-110	O'BRIEN Stitch Scissors, 3-3/4" (9.5 cm), angled, sharp points



MH18-1486



MH18-1652



MH18-1653

REF	Description
MH18-1486	WESTCOTT Stitch Scissors, 4-1/2" (11.4 cm), curved, sharp points
MH18-1652	Eye Suture (GRADLE) Scissors, 3-3/4" (9.5 cm), slightly curved, sharp points
MH18-1653	Eye Suture Scissors, 4" (10.2 cm), ribbon type, slightly curved, sharp points



MH18-1828  
MH18-1832



MH18-1828TC

REF	Description
MH18-1828	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), straight, with lock
MH18-1828TC	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), smooth jaws, straight, with lock, Tungsten Carbide jaws
MH18-1832	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), curved, with lock





MH11-2



MH11-12  
MH11-13



MH11-14



MH11-16

REF	Description
MH11-2	JANSEN Retractor, 4-1/4" (10.8 cm), spread 2" (5.1 cm), 3 X 3 blunt prongs 10 mm deep, standard model
MH11-12	ALM Retractor, 2-3/4" (7 cm), spread 2-1/4" (5.7 cm), prongs 1/8" (.32 cm) deep
MH11-13	ALM Retractor, 3-7/8" (9.8 cm), spread 3" (7.6 cm), prongs 1/2" (1.3 cm) deep
MH11-14	HEISS Self-Retaining Cross Action Retractor, 4" (10.2 cm), sharp prongs
MH11-16	HEISS Self-Retaining Cross Action Retractor, 4" (10.2 cm), blunt prongs



MH11-600  
MH11-602



MH11-608  
MH11-610



MH11-614  
MH11-616



MH11-630-BL  
MH11-630-SH

REF	Description
MH11-600	WEITLANER Retractor, 4" (10.2 cm), sharp, 2 X 3 teeth
MH11-602	WEITLANER Retractor, 4" (10.2 cm), blunt, 2 X 3 teeth
MH11-608	WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth
MH11-610	WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth
MH11-614	WEITLANER Retractor, 6-1/2" (16.5 cm), sharp, 3 X 4 teeth
MH11-616	WEITLANER Retractor, 6-1/2" (16.5 cm), blunt, 3 X 4 teeth
MH11-630-BL	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth, with hinged blades
MH11-630-SH	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth, with hinged blades



MH19-220  
MH19-222


MH19-280 to  
MH19-298

MH19-314

MH19-316

MH19-318

REF	Description
MH19-220	DAY Ear Hook, 6-1/2" (16.5 cm), small size
MH19-222	DAY Ear Hook, 6-1/2" (16.5 cm), large size
MH19-280	BUCK Ear Curette, 6-1/2" (16.5 cm), straight, blunt, size 00
MH19-282	BUCK Ear Curette, 6-1/2" (16.5 cm), straight, blunt, size 0
MH19-284	BUCK Ear Curette, 6-1/2" (16.5 cm), straight, blunt, size 1
MH19-286	BUCK Ear Curette, 6-1/2" (16.5 cm), straight, blunt, size 2
MH19-288	BUCK Ear Curette, 6-1/2" (16.5 cm), straight, blunt, size 3
MH19-290	BUCK Ear Curette, 6-1/2" (16.5 cm), angled, blunt, size 00
MH19-292	BUCK Ear Curette, 6-1/2" (16.5 cm), angled, blunt, size 0
MH19-294	BUCK Ear Curette, 6-1/2" (16.5 cm), angled, blunt, size 1
MH19-296	BUCK Ear Curette, 6-1/2" (16.5 cm), angled, blunt, size 2
MH19-298	BUCK Ear Curette, 6-1/2" (16.5 cm), angled, blunt, size 3
MH19-314	BILLEAU Flexible Ear Loop, 6-1/2" (16.5 cm), small size no. 1
MH19-316	BILLEAU Flexible Ear Loop, 6-1/2" (16.5 cm), medium size no. 2
MH19-318	BILLEAU Flexible Ear Loop, 6-1/2" (16.5 cm), large size no. 3



MH19-700 to  
MH19-720



MH19-793



MH19-796

**SPRATT (BRUN) MASTOID CURETTES, OVAL CUPS**

REF	Description
MH19-700	6-1/2" (16.5 cm), size 5/0
MH19-702	6-1/2" (16.5 cm), size 4/0
MH19-704	6-1/2" (16.5 cm), size 3/0
MH19-706	6-1/2" (16.5 cm), size 2/0
MH19-708	6-1/2" (16.5 cm), size 0
MH19-710	6-1/2" (16.5 cm), size 1

REF	Description
MH19-712	6-1/2" (16.5 cm), size 2
MH19-714	6-1/2" (16.5 cm), size 3
MH19-716	6-1/2" (16.5 cm), size 4
MH19-718	6-1/2" (16.5 cm), size 5
MH19-720	6-1/2" (16.5 cm), size 6

REF	Description
MH19-793	Mallet, 7-1/2" (19.1 cm), lightweight 7.5 oz. (213 g), stainless head 1" (2.5 cm) diameter with 2 replaceable nylon caps
MH19-796	MEAD Mallet, 7-1/2" (19.1 cm), head 8 oz. (227 g), with 2 replaceable nylon caps



MH6-196



MH6-197



MH6-198



MH19-350



MH19-370

REF	Description
MH6-196	GRUENWALD (ADSON) Forceps, 8-1/4" (21 cm), serrated
MH6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps, 6-1/2" (16.5 cm), serrated
MH6-198	ADSON Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated, delicate tips
MH19-350	WILDE Ear Forceps, 5" (12.7 cm), angular with serrated tips
MH19-370	LUCAE Ear Forceps, 5-1/2" (14 cm), bayonet shape with serrated tips



MH19-410



MH19-411



MH19-434

REF	Description
MH19-410	HARTMAN Ear Forceps, 5" (12.7 cm), delicate serrated jaws 2.5 X 6 mm
MH19-411	HARTMAN Alligator Ear Forceps, 5-1/2" (14 cm), serrated jaws
MH19-434	HARTMAN-NOYES Alligator Ear Forceps, 3-1/4" (8.3 cm) shaft, serrated jaws



MH19-2081



MH19-2110

REF	Description
MH19-2081	Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, serrated jaws, 0.6 mm wide
MH19-2110	Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, oval cup jaws, 0.8 mm wide




MH19-570 to  
MH19-577



MH19-580 to  
MH19-584

REF	Description
MH19-570	FRAZIER Suction Tube, 6 French, (2 mm), with finger cut-off, angled
MH19-572	FRAZIER Suction Tube, 7 French, (2.3 mm), with finger cut-off, angled
MH19-573	FRAZIER Suction Tube, 8 French, (2.6 mm), with finger cut-off, angled
MH19-574	FRAZIER Suction Tube, 9 French, (3 mm), with finger cut-off, angled
MH19-575	FRAZIER Suction Tube, 10 French, (3.3 mm), with finger cut-off, angled
MH19-576	FRAZIER Suction Tube, 11 French, (3.6 mm), with finger cut-off, angled
MH19-577	FRAZIER Suction Tube, 12 French, (4 mm), with finger cut-off, angled
MH19-580	BARON Suction Tube, 3 French (1 mm), with finger cut-off, working length 7.5 cm
MH19-582	BARON Suction Tube, 5 French (1.7 mm), with finger cut-off, working length 7.5 cm
MH19-584	BARON Suction Tube, 7 French (2.3 mm), with finger cut-off, working length 7.5 cm




**MH19-844**
**MH19-848**
**MH19-850**
**MH19-1230**
**MH19-849**
**MH19-852**
**MH19-1232**

REF	Description
MH19-844	BEYER Rongeur, 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide
MH19-848	RUSKIN Rongeur, 6" (15.2 cm), double action, straight 3 mm bite jaws
MH19-849	RUSKIN Rongeur, 6" (15.2 cm), double action, curved 3 mm bite jaws
MH19-850	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, straight 5 mm bite jaws
MH19-852	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, curved 5 mm bite jaws
MH19-1230	LEMPERT Rongeur, 6-1/4" (15.9 cm), straight, delicate, 2.5 mm bite
MH19-1232	LEMPERT Rongeur, 6-1/4" (15.9 cm), curved, delicate, 2.5 mm bite



MH6-190

MH6-192

MH6-196

MH6-197

MH6-198

REF	Description
MH6-190	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), bayonet, serrated, scraper end
MH6-192	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), bayonet, 1 X 2 teeth, scraper end
MH6-196	GRUENWALD (ADSON) Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated
MH6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps, 6-1/2" (16.5 cm), serrated
MH6-198	ADSON Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated, delicate tips



MH20-140  
MH20-142



MH20-160

REF	Description
MH20-140	NOYES Alligator Forceps, 5-1/2" (14 cm) shaft, serrated jaws
MH20-142	NOYES Alligator Forceps, 6-1/2" (16.5 cm) shaft, serrated jaws
MH20-160	KNIGHT Nasal Scissors, 6-3/4" (17.1 cm), angled on side, standard pattern



MH20-504  
MH20-506



MH20-522

REF	Description
MH20-504	WATSON-WILLIAMS Nasal Polyp Forceps, 7-3/4" (19.7 cm), small jaws 6 mm wide, shaft 4-1/4" (10.8 cm)
MH20-506	WATSON-WILLIAMS Nasal Polyp Forceps, 7-3/4" (19.7 cm), large jaws 9 mm wide, shaft 4-1/4" (10.8 cm)
MH20-522	KNIGHT Septum Forceps, 6-3/4" (17.1 cm), cup jaws, 15 X 5 mm, shaft 3-1/4" (8.3 cm)



MH19-574



MH20-2



MH20-3



MH20-4

REF	Description
MH19-574	FRAZIER Suction Tube, 9 French, (3 mm), with finger cut-off, angled
MH20-2	VIENNA Nasal Speculum, 5-3/8" (13.7 cm), large size, standard pattern
MH20-3	VIENNA Nasal Speculum, 5-3/8" (13.7 cm), medium size, standard pattern
MH20-4	VIENNA Nasal Speculum, 5-3/8" (13.7 cm), small size, standard pattern



MH20-12



MH20-24



MH20-28 to  
MH20-32



MH20-180

REF	Description
MH20-12	TIECK-HALLE Infant Nasal Speculum, 5-1/2" (14 cm)
MH20-24	BECKMAN Nasal Speculum, 6" (15.2 cm), angled on side
MH20-28	KILLIAN Septum Speculum, 6" (15.2 cm), 2" (5.1 cm) long blades, with set screw
MH20-29	KILLIAN Septum Speculum, 6" (15.2 cm), 2-1/2" (6.4 cm) long blades, with set screw
MH20-30	KILLIAN Septum Speculum, 6" (15.2 cm), 3" (7.6 cm) long blades, with set screw
MH20-32	KILLIAN Septum Speculum, 6" (15.2 cm), 3-1/2" (8.9 cm) long blades, with set screw
MH20-180	KRAUSE Nasal Snare, 10" (25.4 cm), standard pattern with slightly angled cannula


**MH20-200**

**MH20-246**

**MH20-258**

**MH18-1968**

**MH20-278 to  
MH20-284**

REF	Description
MH20-200	FREER Septum Knife, 6-1/4" (15.9 cm), "D" model blade
MH20-246	PENNINGTON Septum Elevator, 8-1/4" (21 cm), sharp, double ended with blades to left and right side
MH20-258	HAJEK-BALLENGER Septum Elevator, 7-1/4" (18.4 cm), double end straight, semi-sharp and blunt blades
MH18-1968	FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
MH20-278	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 2 mm wide
MH20-280	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 4 mm wide
MH20-282	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 6 mm wide
MH20-284	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 8 mm wide



MH20-552  
MH20-553

MH20-554

MH20-556

MH20-570

REF	Description
MH20-552	RUBIN Septal Morselizer, 8" (20.3 cm), double action, deeply serrated jaws 17 X 5 mm, supplied with one slip-on guard
MH20-553	Slip-On Guard for RUBIN Morselizer only, interchangeable
MH20-554	JANSEN-MIDDLETON Septum Forceps, 7-1/2" (19.1 cm), spoon-shaped cups 15 X 5 mm
MH20-556	JANSEN-MIDDLETON Septum Forceps, 7-1/2" (19.1 cm), through-cutting jaws 19 X 4 mm
MH20-570	TAKAHASHI Ethmoid Forceps, 4-1/2" (11.4 cm) shaft, 2.5 X 10 mm bite





MH20-800 to  
MH20-808




MH20-852 to  
MH20-856



MH20-858-1 to  
MH20-858-3

REF	Description
MH20-800	BARNHILL Adenoid Curette, 8-1/2" (21.6 cm), size 0, 11 mm blade
MH20-802	BARNHILL Adenoid Curette, 8-1/2" (21.6 cm), size 1, 13 mm blade
MH20-804	BARNHILL Adenoid Curette, 8-1/2" (21.6 cm), size 2, 15 mm blade
MH20-806	BARNHILL Adenoid Curette, 8-1/2" (21.6 cm), size 3, 17 mm blade
MH20-808	BARNHILL Adenoid Curette, 8-1/2" (21.6 cm), size 4, 19 mm blade
MH20-852	LAFORCE Adenotome, 9-1/2" (24.1 cm), small size 1, cuts 11 mm width
MH20-854	LAFORCE Adenotome, 9-1/2" (24.1 cm), medium size 2, cuts 13 mm width
MH20-856	LAFORCE Adenotome, 9-1/2" (24.1 cm), large size 3, cuts 16 mm width
MH20-858-1	Blades for LAFORCE Adenotome, size 1, 13 mm
MH20-858-2	Blades for LAFORCE Adenotome, size 2, 16 mm
MH20-858-3	Blades for LAFORCE Adenotome, size 3, 19 mm



MH5-179  
MH5-180

MH5-181

MH5-182

MH5-288

REF	Description
MH5-179	METZENBAUM (LAHEY) Scissors, Standard Pattern, 5-1/2" (14 cm), straight, blunt points
MH5-180	METZENBAUM (LAHEY) Scissors, Standard Pattern, 5-1/2" (14 cm), curved, blunt points
MH5-181	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), straight, blunt points
MH5-182	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), curved, blunt points
MH5-288	METZENBAUM Scissors, Delicate Pattern, 7" (17.8 cm), curved, blunt points



MH5-270  
MH5-272

MH5-274  
MH5-276

MH5-278

MH5-280

### PLASTIC SURGERY AND DELICATE SCISSORS

REF	Description
MH5-270	Delicate Scissors, 4-3/4" (12.1 cm), straight, sharp-sharp points
MH5-272	Delicate Scissors, 4-3/4" (12.1 cm), straight, sharp-blunt points
MH5-274	Delicate Scissors, 4-3/4" (12.1 cm), straight, blunt-blunt points
MH5-276	Delicate Scissors, 4-3/4" (12.1 cm), curved, sharp-sharp points
MH5-278	Delicate Scissors, 4-3/4" (12.1 cm), curved, sharp-blunt points
MH5-280	Delicate Scissors, 4-3/4" (12.1 cm), curved, blunt-blunt points



MH5-300  
MH5-304

MH5-302  
MH5-306

MH18-1396

MH18-1398

REF	Description
MH5-300	Iris Scissors, 4" (10.2 cm), straight
MH5-302	Iris Scissors, 4" (10.2 cm), curved
MH5-304	Iris Scissors, 4-1/2" (11.4 cm), straight
MH5-306	Iris Scissors, 4-1/2" (11.4 cm), curved
MH18-1396	Iris Scissors, 3-1/2" (8.9 cm), straight, with 20 mm blades, delicate
MH18-1398	Iris Scissors, 3-1/2" (8.9 cm), curved, with 20 mm blades, delicate



MH18-1415

MH18-1416

MH18-1417

MH18-1418

REF	Description
MH18-1415	Iris Scissors, 4" (10.2 cm), straight, standard pattern, ribbon-type
MH18-1416	Iris Scissors, 4" (10.2 cm), curved, standard pattern, ribbon-type
MH18-1417	Iris Scissors, 3-3/4" (9.5 cm), straight, miniature pattern with 14 mm blades, ribbon-type
MH18-1418	Iris Scissors, 3-3/4" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon-type



MH18-1420

MH18-1422

MH18-1426

MH5-312

MH18-1428

MH5-314

REF	Description
MH18-1420	KNAPP Iris Scissors, 4" (10.2 cm), straight, sharp-blunt points
MH18-1422	KNAPP Iris Scissors, 4" (10.2 cm), straight, blunt-blunt points
MH18-1426	KNAPP Iris Scissors, 4" (10.2 cm), curved, sharp-blunt points
MH18-1428	KNAPP Iris Scissors, 4" (10.2 cm), curved, blunt-blunt points
MH5-312	Strabismus Scissors, 4" (10.2 cm), straight
MH5-314	Strabismus Scissors, 4" (10.2 cm), curved



MH18-1460



MH18-1462 to MH18-1466



MH18-1472 to MH18-1476

### STEVENS TENOTOMY SCISSORS

REF	Description
MH18-1460	4-1/8" (10.5 cm), straight, short blades, sharp points
MH18-1462	4-1/8" (10.5 cm), straight, short blades, blunt points
MH18-1464	4-1/8" (10.5 cm), curved, short blades, sharp points
MH18-1466	4-1/8" (10.5 cm), curved, short blades, blunt points
MH18-1472	4-1/2" (11.4 cm), straight, long blades, blunt points
MH18-1476	4-1/2" (11.4 cm), curved, long blades, blunt points



MH9-100

MH9-102  
MH9-104

MH9-110

REF	Description
MH9-100	SPENCER Stitch Scissors, 3-1/2" (8.9 cm), delicate
MH9-102	LITTAUER JR. Stitch Scissors, 4-1/2" (11.4 cm), light pattern
MH9-104	LITTAUER Stitch Scissors, 5-1/2" (14 cm), standard pattern
MH9-110	O'BRIEN Stitch Scissors, 3-3/4" (9.5 cm), angled, sharp points





MH9-124



MH9-124TC



MH18-1486

REF	Description
MH9-124	Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade
MH9-124TC	Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade, Tungsten Carbide blades
MH18-1486	WESTCOTT Stitch Scissors, 4-1/2" (11.4 cm), curved, sharp points



MH18-1652



MH18-1653



MH21-536

REF	Description
MH18-1652	Eye Suture (GRADLE) Scissors, 3-3/4" (9.5 cm), slightly curved, sharp points
MH18-1653	Eye Suture Scissors, 4" (10.2 cm), ribbon type, slightly curved, sharp points
MH21-536	LITTLER Suture Carrying Scissors, 4-5/8" (11.8 cm), with suture hole in blades, curved



MH21-600




MH21-602



MH21-604

REF	Description
MH21-600	FOMON Scissors, Saber-Back, 5" (12.7 cm), slightly curved, semi-sharp outer edges, used for dissection
MH21-602	FOMON Lower Lateral Scissors, 5" (12.7 cm), strong curve, dissecting blades
MH21-604	FOMON Dorsal Scissors, 5-1/2" (14 cm), light pattern, angular shanks



MH21-610



MH21-630

MH21-632

REF	Description
MH21-610	COTTLE Dorsal Scissors, 6-1/2" (16.5 cm), heavy pattern with rounded blades, angular shanks
MH21-630	JOSEPH Scissors, 5-3/4" (14.6 cm), straight, sharp points
MH21-632	JOSEPH Scissors, 5-3/4" (14.6 cm), curved, sharp points




MH6-104



MH6-106

REF	Description
MH6-104	SEMKEN Dressing Forceps, 5" (12.7 cm), serrated
MH6-106	SEMKEN Tissue Forceps, 5" (12.7 cm), 1 X 2 teeth



MH6-118

MH6-120  
MH6-122

MH6-121

MH6-123

MH6-124

MH17-2500  
MH17-2510

REF	Description
MH6-118	ADSON Dressing Forceps, 4-3/4" (12.1 cm), delicate, serrated
MH6-120	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, delicate
MH6-121	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, cross serrated tips, delicate
MH6-122	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 2 X 3 teeth, delicate
MH6-123	ADSON Tissue & Suture Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, with tying platform
MH6-124	BROWN-ADSON Tissue Forceps, 4-3/4" (12.1 cm), 7 X 7 side-grasping teeth
MH17-2500	Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, tips 0.5 mm wide
MH17-2510	Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated tips, 0.5 mm wide



MH18-780




MH18-782



MH18-784

REF	Description
MH18-780	Eye Dressing Forceps, 4" (10.2 cm), straight, standard pattern, 0.8 mm wide serrated tips
MH18-782	Eye Dressing Forceps, 4" (10.2 cm), half curved, standard pattern, 0.8 mm wide serrated tips
MH18-784	Eye Dressing Forceps, 4" (10.2 cm), full curved, standard pattern, 0.8 mm wide serrated tips



MH18-786

MH18-788

MH18-790

MH18-866

MH18-867

MH18-868

MH21-440

REF	Description
MH18-786	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, straight, standard pattern, 0.8 mm wide tips
MH18-788	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, half curved, standard pattern, 0.8 mm wide tips
MH18-790	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, full curved, standard pattern, 0.8 mm wide tips
MH18-866	BISHOP-HARMON Dressing Forceps, 3-3/8" (8.6 cm), fine cross serrated tips 0.5 mm wide
MH18-867	BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth
MH18-868	BISHOP-HARMON Tissue Forceps, 3-3/8" (8.6 cm), standard 0.7 mm 1 X 2 teeth
MH21-440	COTTLE Columella Forceps, 4-1/2" (11.4 cm), jaws 10 mm wide






MH7-2  
MH7-4




MH7-8  
MH7-10



MH7-24  
MH7-26

REF	Description
MH7-2	HALSTED Mosquito Forceps, 5" (12.7 cm), straight
MH7-4	HALSTED Mosquito Forceps, 5" (12.7 cm), curved
MH7-8	HALSTED Mosquito Forceps, 5" (12.7 cm), straight, extra delicate
MH7-10	HALSTED Mosquito Forceps, 5" (12.7 cm), curved, extra delicate
MH7-24	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), straight
MH7-26	HARTMAN Mosquito Forceps, 3-1/2" (8.9 cm), curved




MH7-50  
MH7-52




MH17-2600  
MH17-2602

REF	Description
MH7-50	Baby CRILE Forceps, 5-1/2" (14 cm), extra delicate, straight
MH7-52	Baby CRILE Forceps, 5-1/2" (14 cm), extra delicate, curved
MH17-2600	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), straight, extremely delicate
MH17-2602	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), curved, extremely delicate



REF	Description
MH8-2	COLLIER Needle Holder, 5" (12.7 cm), fenestrated jaws
MH8-4	DERF Needle Holder, 4-3/4" (12.1 cm)
MH8-6	WEBSTER Needle Holder, 5" (12.7 cm), delicate, smooth jaws, standard pattern
MH8-8	HALSEY Needle Holder, 5" (12.7 cm), smooth jaws



MH8-14

MH8-15  
MH8-16

MH8-20

MH8-22

REF	Description
MH8-14	OLSEN-HEGAR Needle Holder with Suture Scissors, 4-3/4" (12.1 cm), serrated jaws, extra delicate pattern
MH8-15	OLSEN-HEGAR Needle Holder with Suture Scissors, 5-1/2" (14 cm), serrated jaws
MH8-16	OLSEN-HEGAR Needle Holder with Suture Scissors, 6-1/2" (16.5 cm), serrated jaws
MH8-20	CONVERSE Needle Holder, 4-1/4" (10.8 cm), longitudinal serrations
MH8-22	BROWN Needle Holder, 5-1/4" (13.3 cm), convex jaws, fine serrations


**MH8-40**
**MH8-42 to  
MH8-46**
**MH8-50  
MH8-51**
**MH18-1828  
MH18-1828TC  
MH18-1832**

REF	Description
MH8-40	BAUMGARTNER Needle Holder, 5-1/2" (14 cm)
MH8-42	MAYO-HEGAR Needle Holder, 5" (12.7 cm)
MH8-44	MAYO-HEGAR Needle Holder, 6-1/4" (15.9 cm)
MH8-46	MAYO-HEGAR Needle Holder, 7-1/4" (18.4 cm)
MH8-50	CRILE-WOOD Needle Holder, 6" (15.2 cm)
MH8-51	CRILE-WOOD Needle Holder, 7" (17.8 cm)
MH18-1828	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), straight, with lock
MH18-1828TC	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), straight, with lock, Tungsten Carbide
MH18-1832	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), curved, with lock



MH11-73

MH11-74


MH11-76

MH18-100

MH18-102

MH18-110

REF	Description
MH11-73	RAGNELL Retractor, 6" (15.2 cm), double end 3 X 8 mm and 5 X 15 mm blades, delicate
MH11-74	SENN Retractor, 6-3/8" (16.2 cm), sharp, double end
MH11-76	SENN Retractor, 6-3/8" (16.2 cm), blunt, double end
MH18-100	KNAPP (ROLLET) Retractor, 5-1/4" (13.3 cm), with 4 sharp prongs
MH18-102	KNAPP Retractor, 5-1/4" (13.3 cm), with 4 blunt prongs
MH18-110	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 1, 11 mm



MH21-144



MH21-145

MH21-146

MH21-145  
MH21-146



MH21-148



MH23-1049



MH26-1084

REF	Description
MH21-144	COTTLE Double Hook, 5-1/2" (14 cm), two sharp prongs, 12 mm wide
MH21-145	COTTLE Retractor, 5-1/2" (14 cm), right ball tip
MH21-146	COTTLE Retractor, 5-1/2" (14 cm), left ball tip
MH21-148	FOMON Retractor, 6-1/2" (14 cm), two prongs with ball tips, 11 mm wide
MH23-1049	JACKSON Trachea Retractor, 7" (17.8 cm), double ended, 2 blunt prongs 5/16" (2.4 cm) wide, and solid blade 3/8" (1 cm) wide
MH26-1084	FRAZIER Dura Hook, 5" (12.7 cm), sharp point



MH18-412

MH18-414

MH18-416

MH21-88

MH21-89

MH21-90 to

MH21-94

REF	Description
MH18-412	TYRELL Hook, 5" (12.7 cm), sharp
MH18-414	TYRELL Hook, 5" (12.7 cm), blunt
MH18-416	GUTHRIE Fixation Hook, 5" (12.7 cm), sharp double prongs 1.5 mm wide
MH21-88	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), large
MH21-89	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), small
MH21-90	COTTLE Skin Hook, 5-1/2" (14 cm), small, shallow curve
MH21-92	COTTLE Skin Hook, 5-1/2" (14 cm), small, deep curve
MH21-94	COTTLE Skin Hook, 5-1/2" (14 cm), large, deep curve




**MH21-152**

**MH21-153**

**MH21-154 to  
MH21-160**

**MH21-165**

REF	Description
MH21-152	Skin Hook, 4-3/4" (12.1 cm), sharp prong, 2 mm diameter
MH21-153	JOSEPH Hook, 6-1/4" (15.9 cm), one sharp prong
MH21-154	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 2 mm wide
MH21-156	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 5 mm wide
MH21-158	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 7 mm wide
MH21-160	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 10 mm wide
MH21-165	COTTLE Knife Guide & Retractor, 8" (20.3 cm), double ended, ball tip prongs 10 mm wide, and 5 mm wide slotted hook



MH19-793



MH19-796

REF	Description
MH19-793	Mallet, 7-1/2" (19.1 cm), lightweight 7.5 oz. (213 g), stainless steel head 1" (2.5 cm) diameter with 2 replaceable nylon caps
MH19-796	MEAD Mallet, 7-1/2" (19.1 cm), head 8 oz. (227 g), with 2 replaceable nylon caps



MH20-278 to  
MH20-284



MH21-20



MH21-24



MH21-32

REF	Description
MH20-278	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 2 mm wide
MH20-280	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 4 mm wide
MH20-282	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 6 mm wide
MH20-284	BALLENGER V-Shaped Chisel, 6-1/4" (15.9 cm), 8 mm wide
MH21-20	JOSEPH Button End Knife, 6" (15.2 cm), straight blade 4 X 30 mm
MH21-24	CONVERSE Bistouri, 6" (15.2 cm), with button end, curved blade 4 X 35 mm
MH21-32	COTTLE Nasal Knife, 5-1/2" (14 cm), straight blade, 4 mm wide



MH18-1968



MH21-60



MH21-64



MH21-72

REF	Description
MH18-1968	FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
MH21-60	COTTLE Septum Elevator, 9" (22.9 cm), double end, graduated in cm
MH21-64	JOSEPH Periosteal Elevator, 7" (17.8 cm), strong curved blade 6.5 mm wide
MH21-72	COTTLE Skin Elevator, 7-1/2" (19.1 cm), slightly curved blade 9 mm wide


**MH21-188**

**MH21-205-2 to  
MH21-205-12**

**MH21-206-4 to  
MH21-206-12**

**MH21-232 to  
MH21-234**

REF	Description
MH21-188	AUFRICHT Retractor-Speculum, 6-1/2" (16.5 cm), solid blade 45 mm long
MH21-205-2	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 2 mm wide
MH21-205-4	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 4 mm wide
MH21-205-6	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 6 mm wide
MH21-205-8	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 8 mm wide
MH21-205-10	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 10 mm wide
MH21-205-12	SHEEHAN Osteotome, 6-1/4" (15.9 cm), 12 mm wide
MH21-206-4	COTTLE Osteotome, 7" (17.8 cm), 4 mm wide, rounded corners
MH21-206-7	COTTLE Osteotome, 7" (17.8 cm), 7 mm wide, rounded corners
MH21-206-9	COTTLE Osteotome, 7" (17.8 cm), 9 mm wide, rounded corners
MH21-206-12	COTTLE Osteotome, 7" (17.8 cm), 12 mm wide, rounded corners
MH21-232	NEIVERT Osteotome, 8" (20.3 cm), straight, with single guard
MH21-233	NEIVERT-ANDERSON Osteotome, 8" (20.3 cm), curved right, with single guard
MH21-234	NEIVERT-ANDERSON Osteotome, 8" (20.3 cm), curved left, with single guard



MH21-280



MH21-300



MH21-330

REF	Description
MH21-280	JOSEPH Nasal Saw, 7-1/2" (19.1 cm), straight
MH21-300	COTTLE Cartilage Crusher, 16 mm wide anvil, with removing forceps
MH21-330	JOSEPH Rasp, 6-1/4" (15.9 cm), fine cross serrations



MH21-338



MH21-342



MH21-346  
MH21-348



MH21-368 to  
MH21-370

REF	Description
MH21-338	COTTLE Rasp, 8" (20.3 cm), backward cutting
MH21-342	MALTZ Nasal Rasp, 7" (17.8 cm), backward cutting
MH21-346	AUFRICHT Glabella Rasp, 8-1/4" (21 cm), curved blade, forward cutting
MH21-348	AUFRICHT Glabella Rasp, 8-1/4" (21 cm), curved blade, backward cutting
MH21-368	FOMON Rasp, 8-1/4" (21 cm), double ended, 4-sided, fine & coarse teeth, both in flat and convex
MH21-369	FOMON Rasp, 8-1/4" (21 cm), double ended, coarse teeth, flat and convex
MH21-370	FOMON Rasp, 8-1/4" (21 cm), double ended, fine teeth, flat and convex



MH21-626  
MH21-627



MH21-660

REF	Description
MH21-626	KLEINERT-KUTZ (RUSKIN) Bone Splitting Forceps, 6" (15.2 cm), straight
MH21-627	KLEINERT-KUTZ (RUSKIN) Bone Splitting Forceps, 6" (15.2 cm), curved
MH21-660	ASH Septum Straightening Forceps, 9" (22.9 cm), angled blades flat inside





**MH5-179**  
**MH5-180**



**MH5-300**  
**MH5-304**



**MH5-302**  
**MH5-306**



**MH5-312**  
**MH5-314**

REF	Description
MH5-179	METZENBAUM (LAHEY) Scissors, Standard Pattern, 5-1/2" (14 cm), straight, blunt points
MH5-180	METZENBAUM (LAHEY) Scissors, Standard Pattern 5-1/2" (14 cm), curved, blunt points
MH5-300	Iris Scissors, 4" (10.2 cm), straight
MH5-302	Iris Scissors, 4" (10.2 cm), curved
MH5-304	Iris Scissors, 4-1/2" (11.4 cm), straight
MH5-306	Iris Scissors, 4-1/2" (11.4 cm), curved
MH5-312	Strabismus Scissors, 4" (10.2 cm), straight
MH5-314	Strabismus Scissors, 4" (10.2 cm), curved



MH18-1396

MH18-1398

MH18-1415

MH18-1416

REF	Description
MH18-1396	Iris Scissors, 3-1/2" (8.9 cm), straight, with 20 mm blades, delicate
MH18-1398	Iris Scissors, 3-1/2" (8.9 cm), curved, with 20 mm blades, delicate
MH18-1415	Iris Scissors, 4" (10.2 cm), straight, standard pattern, ribbon-type
MH18-1416	Iris Scissors, 4" (10.2 cm), curved, standard pattern, ribbon-type



MH18-1417



MH18-1418



MH18-1420

MH18-1422



MH18-1426

MH18-1428

REF	Description
MH18-1417	Iris Scissors, 3-3/4" (9.5 cm), straight, miniature pattern with 14 mm blades, ribbon-type
MH18-1418	Iris Scissors, 3-3/4" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon-type
MH18-1420	KNAPP Iris Scissors, 4" (10.2 cm), straight, sharp-blunt points
MH18-1422	KNAPP Iris Scissors, 4" (10.2 cm), straight, blunt-blunt points
MH18-1426	KNAPP Iris Scissors, 4" (10.2 cm), curved, sharp-blunt points
MH18-1428	KNAPP Iris Scissors, 4" (10.2 cm), curved, blunt-blunt points



MH18-1460



MH18-1462 to MH18-1466



MH18-1472 to MH18-1476

### STEVENS TENOTOMY SCISSORS

REF	Description
MH18-1460	4-1/8" (10.5 cm), straight, short blades, sharp points
MH18-1462	4-1/8" (10.5 cm), straight, short blades, blunt points
MH18-1464	4-1/8" (10.5 cm), curved, short blades, sharp points
MH18-1466	4-1/8" (10.5 cm), curved, short blades, blunt points
MH18-1472	4-1/2" (11.4 cm), straight, long blades, blunt points
MH18-1476	4-1/2" (11.4 cm), curved, long blades, blunt points



MH21-630  
MH21-632



MH18-1486

REF	Description
MH21-630	JOSEPH Scissors, 5-3/4" (14.6 cm), straight, sharp points
MH21-632	JOSEPH Scissors, 5-3/4" (14.6 cm), curved, sharp points
MH18-1486	WESTCOTT Stitch Scissors, 4-1/2" (11.4 cm), curved, sharp points



MH18-1652

MH18-1653

MH9-100

MH9-102

MH9-110

MH9-104

REF	Description
MH18-1652	Eye Suture (GRADLE) Scissors, 3-3/4" (9.5 cm), slightly curved, sharp points
MH18-1653	Eye Suture Scissors, 4" (10.2 cm), ribbon type, slightly curved, sharp points
MH9-100	SPENCER Stitch Scissors, 3-1/2" (8.9 cm), delicate
MH9-102	LITTAUER JR. Stitch Scissors, 4-1/2" (11.4 cm), light pattern
MH9-104	LITTAUER Stitch Scissors, 5-1/2" (14 cm), standard pattern
MH9-110	O'BRIEN Stitch Scissors, 3-3/4" (9.5 cm), angled, sharp points



MH6-104



MH6-106



MH18-780




MH18-782



MH18-784

REF	Description
MH6-104	SEMKEN Dressing Forceps, 5" (12.7 cm), serrated
MH6-106	SEMKEN Tissue Forceps, 5" (12.7 cm), 1 X 2 teeth
MH18-780	Eye Dressing Forceps, 4" (10.2 cm), straight, standard pattern, 0.8 mm wide serrated tips
MH18-782	Eye Dressing Forceps, 4" (10.2 cm), half curved, standard pattern, 0.8 mm wide serrated tips
MH18-784	Eye Dressing Forceps, 4" (10.2 cm), full curved, standard pattern, 0.8 mm wide serrated tips



MH18-786

MH18-788

MH18-790


MH18-866

MH18-867

MH18-868

REF	Description
MH18-786	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, straight, standard pattern, 0.8 mm wide tips
MH18-788	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, half curved, standard pattern, 0.8 mm wide tips
MH18-790	Iris Tissue Forceps, 4" (10.2 cm), 1 X 2 teeth, full curved, standard pattern, 0.8 mm wide tips
MH18-866	BISHOP-HARMON Dressing Forceps, 3-3/8" (8.6 cm), fine cross serrated tips 0.5 mm wide
MH18-867	BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth
MH18-868	BISHOP-HARMON Tissue Forceps, 3-3/8" (8.6 cm), standard 0.7 mm 1 X 2 teeth




**MH6-118**
**MH6-120**  
**MH6-122**
**MH6-121**
**MH6-123**
**MH6-124**
**MH17-2500**  
**MH17-2510**

REF	Description
MH6-118	ADSON Dressing Forceps, 4-3/4" (12.1 cm), delicate, serrated
MH6-120	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, delicate
MH6-121	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, cross serrated tips, delicate
MH6-122	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 2 X 3 teeth, delicate
MH6-123	ADSON Dressing, Tissue & Suture Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, with tying platform
MH6-124	BROWN-ADSON Tissue Forceps, 4-3/4" (12.1 cm), 7 X 7 side grasping teeth
MH17-2500	Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, tips 0.5 mm wide
MH17-2510	Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated tips, 0.5 mm wide



MH8-4



MH8-6



MH8-8

REF	Description
MH8-4	DERF Needle Holder, 4-3/4" (12.1 cm)
MH8-6	WEBSTER Needle Holder, 5" (12.7 cm), delicate, smooth jaws, standard pattern
MH8-8	HALSEY Needle Holder, 5" (12.7 cm), smooth jaws



MH8-14



MH8-40



MH18-1828

MH18-1832

REF	Description
MH8-14	OLSEN-HEGAR Needle Holder with Suture Scissors, 4-3/4" (12.1 cm), serrated jaws, extra delicate pattern
MH8-40	BAUMGARTNER Needle Holder, 5-1/2" (14 cm)
MH18-1828	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), straight, with lock
MH18-1832	CASTROVIEJO Needle Holder, 5-1/2" (14 cm), curved, with lock



MH11-12  
MH11-13




MH11-14




MH11-16

REF	Description
MH11-12	ALM Retractor, 2-3/4" (7 cm), spread 2-1/4" (5.7 cm), prongs 1/8" (.32 cm) deep
MH11-13	ALM Retractor, 3-7/8" (9.8 cm), spread 3" (7.6 cm), prongs 1/2" (1.3 cm) deep
MH11-14	HEISS Self-Retaining, Cross Action Retractor, 4" (10.2 cm), sharp prongs
MH11-16	HEISS Self-Retaining, Cross Action Retractor, 4" (10.2 cm), blunt prongs



MH11-74

MH11-76



MH18-100

MH18-102

REF	Description
MH11-74	SENN Retractor, 6-3/8" (16.2 cm), sharp, double end
MH11-76	SENN Retractor, 6-3/8" (16.2 cm), blunt, double end
MH18-100	KNAPP (ROLLET) Retractor, 5-1/4" (13.3 cm), with 4 sharp prongs
MH18-102	KNAPP Retractor, 5-1/4" (13.3 cm), with 4 blunt prongs



MH18-110 to  
MH18-116



MH18-412  
MH18-414



MH18-416



MH26-1084

REF	Description
MH18-110	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 1, 11 mm
MH18-112	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 2, 13 mm
MH18-114	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 3, 15 mm
MH18-116	DESMARRES Lid Retractor, 5-1/2" (14 cm), size 4, 17 mm
MH18-412	TYRELL Hook, 5" (12.7 cm), sharp
MH18-414	TYRELL Hook, 5" (12.7 cm), blunt
MH18-416	GUTHRIE Fixation Hook, 5" (12.7 cm), sharp double prongs 1.5 mm wide
MH26-1084	FRAZIER Dura Hook, 5" (12.7 cm), sharp point



MH21-88  
MH21-89



MH21-90 to  
MH21-94



MH21-144



MH21-145  
MH21-146


REF	Description
MH21-88	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), large
MH21-89	GILLIES (CONVERSE) Skin Hook, 7" (17.8 cm), small
MH21-90	COTTLE Skin Hook, 5-1/2" (14 cm), small, shallow curve
MH21-92	COTTLE Skin Hook, 5-1/2" (14 cm), small, deep curve
MH21-94	COTTLE Skin Hook, 5-1/2" (14 cm), large, deep curve
MH21-144	COTTLE Double Hook, 5-1/2" (14 cm), two sharp prongs, 12 mm wide
MH21-145	COTTLE Double Hook, 5-1/2" (14 cm), right ball tip
MH21-146	COTTLE Double Hook, 5-1/2" (14 cm), left ball tip



MH21-152



MH21-153

MH21-154 to  
MH21-160MH33-140  
MH33-142MH33-200  
MH33-201 and  
MH33-202 (not shown)

REF	Description
MH21-152	Skin Hook, 4-3/4" (12.1 cm), sharp prong, 2 mm diameter
MH21-153	JOSEPH Hook, 6-1/4" (15.9 cm), one sharp prong
MH21-154	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 2 mm wide
MH21-156	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 5 mm wide
MH21-158	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 7 mm wide
MH21-160	JOSEPH Double Hook, 6-1/4" (15.9 cm), two sharp prongs, 10 mm wide
MH33-140	Finger Ring Cutter, 6-1/2" (16.5 cm), with circular saw blade
MH33-142	Replacement Blade for Finger Ring Cutter, chrome
MH33-200	SCHAMBERG Comedone Extractor, 3-3/4" (9.5 cm), with square loop ends
MH33-201	SCHAMBERG Comedone Extractor, 3-3/4" (9.5 cm), standard pattern with crimped small loop
MH33-202	SCHAMBERG Comedone Extractor, 3-3/4" (9.5 cm), standard pattern with extra fine loops, small loop crimped





MH2-104SS  
MH2-109



MH2-110SS

REF	Description
MH2-104SS	YANKAUER Suction Tube, 11-3/4" (29.8 cm), with removable tip and tubing connection
MH2-106SS	Tip only for YANKAUER Suction Tube (MH2-104SS)
MH2-109	YANKAUER Pediatric Suction Tube, 8-1/2" (21.6 cm), with removable tip, delicate pattern
MH2-110SS	ANDREWS-PYNCHON Suction Tube, 10" (25.4 cm), delicate pattern



MH2-120 to MH2-136

MH2-180

REF	Description
MH2-120	DAVIS Mouth Gag, 6-1/4" (15.9 cm), double bite, left, without blades
MH2-122	DAVIS Mouth Gag, 6-1/4" (15.9 cm), double bite, right, without blades
MH2-180	McIVOR Mouth Gag, 5-3/4" (14.6 cm), complete with 3 tongue blades, one each Nos. MH2-132, MH2-134, MH2-136
MH2-132	Blades for DAVIS or McIVOR Mouth Gags, size No. 2, 1" X 2-3/4" (2.5 X 7 cm) tongue blade, without Ether tube
MH2-134	Blades for DAVIS or McIVOR Mouth Gags, size No. 3, 1" X 3-1/4" (2.5 X 8.3 cm), tongue blade, without Ether tube
MH2-136	Blades for DAVIS or McIVOR Mouth Gags, size No. 4, 1" X 3-3/4" (2.5 X 9.5 cm), tongue blade, without Ether tube



MH2-140 to  
MH2-147

REF	Description
* MH2-140	Ring Tongue Retractor Blades, size No. 2, 7/8" X 2-1/4" (2.2 X 5.7 cm), left
* MH2-141	Ring Tongue Retractor Blades, size No. 3, 1-1/8" X 3" (2.9 X 7.6 cm), left
* MH2-142	Ring Tongue Retractor Blades, size No. 4, 1-1/4" X 3-1/2" (3.2 X 8.9 cm), left
* MH2-143	Ring Tongue Retractor Blades, size No. 5, 1-1/2" X 3-3/4" (3.8 X 9.5 cm), left
* MH2-144	Ring Tongue Retractor Blades, size No. 2, 7/8" X 2-1/4" (2.2 X 5.7 cm), right
* MH2-145	Ring Tongue Retractor Blades, size No. 3, 1-1/8" X 3" (2.9 X 7.6 cm), right
* MH2-146	Ring Tongue Retractor Blades, size No. 4, 1-1/4" X 3-1/2" (3.2 X 8.9 cm), right
* MH2-147	Ring Tongue Retractor Blades, size No. 5, 1-1/2" X 3-3/4" (3.8 X 9.5 cm), right

\*For use with Davis or McIvor mouth gags. Grooved blade holds endotracheal tube firmly in place.



MH2-149  
MH2-150



MH2-160 to  
MH2-164



MH18-1968

REF	Description
MH2-149	DOYEN-JANSEN Mouth Gag, 4-1/2" (11.4 cm), with grip lock, Child
MH2-150	DOYEN-JANSEN Mouth Gag, 5-1/2" (14 cm), with grip lock, Adult
MH2-160	JENNINGS Mouth Gag, adult size, 6" (15.2 cm), wide
MH2-161	JENNINGS Mouth Gag, adult size, 5" (12.7 cm), wide
MH2-162	JENNINGS Mouth Gag, child size, 4-3/4" (12.1 cm), wide
MH2-164	JENNINGS Mouth Gag, infant size, 3-3/4" (9.5 cm), wide
MH18-1968	FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide



MH19-844



MH22-2

MH22-4

REF	Description
MH19-844	BEYER Rongeur, 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide
MH22-2	WIEDER Tongue Depressor, small blade 28 mm wide, serrated
MH22-4	WIEDER Tongue Depressor, large blade 36 mm wide, serrated



MH5-182



MH5-264



MH7-604

MH7-610

MH7-614

REF	Description
MH5-182	METZENBAUM Scissors, Standard Pattern, 7" (17.8 cm), curved, blunt points
MH5-264	DEAN Dissecting Scissors, 6-3/4" (17.1 cm), blades angled on flat, one serrated blade
MH7-604	FOERSTER Sponge Forceps, 7" (17.8 cm), straight, serrated jaws
MH7-610	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), curved, serrated jaws
MH7-614	FOERSTER Sponge Forceps, 7" (17.8 cm), curved, serrated jaws



MH22-550  
MH22-552



MH22-662



MH22-668

REF	Description
MH22-550	WHITE Tonsil Seizing Forceps, 7-1/2" (19.1 cm), medium curve, one open ring
MH22-552	WHITE Tonsil Seizing Forceps, 9-1/2" (24.1 cm), medium curve, one open ring
MH22-662	SCHNIDT Tonsil Hemostatic Forceps, 7-1/2" (19.1 cm), slight curve, closed rings
MH22-668	SAWTELL Tonsil Hemostatic Forceps, 7-1/2" (19.1 cm), full curve, closed rings



MH22-808



MH22-854



MH22-1000

REF	Description
MH22-808	FISHER Tonsil Knife and Dissector, 8-1/2" (21.6 cm), semi-sharp tip, with teeth
MH22-854	HURD Dissector and Pillar Retractor, 9" (22.9 cm), standard pattern, 10.2 mm X 13.9 mm blade
MH22-1000	EVE Tonsil Snare, 11" (27.9 cm), standard model, without ratchet




**MH6-214**
**MH11-73**
**MH11-74  
MH11-76**
**MH12-80**
**MH23-1049**

REF	Description
MH6-214	SINGLEY Tissue Forceps, 9" (22.9 cm), fenestrated serrated jaws
MH11-73	RAGNELL Retractor, 6" (15.2 cm), double end 3 X 8 mm and 5 X 15 mm blades, delicate
MH11-74	SENN Retractor, 6-3/8" (16.2 cm), sharp, double end
MH11-76	SENN Retractor, 6-3/8" (16.2 cm), blunt, double end
MH12-80	LAHEY Traction Forceps, 6-1/4" (15.9 cm), 3 X 3 teeth
MH23-1049	JACKSON Trachea Retractor, 7" (17.8 cm), double ended, 2 blunt prongs 5/16" (2.4 cm) wide, and solid blade 3/8" (1 cm) wide




MH23-1050



MH3-1051



MH23-1420



MH23-1430

REF	Description
MH23-1050	JACKSON Tenaculum Hook Retractor, 5-3/4" (14.6 cm), sharp point
MH23-1051	NEWS Hook, 6" (15.2 cm), single sharp prong
MH23-1420	TROUSSEAU Trachea Dilator, 5-1/4" (13.3 cm), standard size, jaws 3 cm long X 5 mm wide at tip
MH23-1430	LaBORDE Trachea Dilator, 5-1/2" (14 cm), three bladed



MH24-800 to  
MH24-808



MH24-1180  
MH24-1190  
MH24-1200

REF	Description
MH24-800	DeBAKEY Bulldog Clamp, 5" (12.7 cm), ring handle, straight
MH24-802	DeBAKEY Bulldog Clamp, 4-3/4" (12.1 cm), ring handle, curved
MH24-804	DeBAKEY Bulldog Clamp, 4" (10.2 cm), ring handle, 90 degree angle
MH24-806	DeBAKEY Bulldog Clamp, 4-7/8" (12.4 cm), ring handle, 45 degree angle
MH24-808	DeBAKEY Bulldog Clamp, 4-3/4" (12.1 cm), ring handle, S-curved
MH24-1180	DeBAKEY Peripheral Vascular Clamp, 7" (17.8 cm), angled jaws 6.5 cm
MH24-1190	DeBAKEY Peripheral Vascular Clamp, 7" (17.8 cm), S-shaped
MH24-1200	DeBAKEY Peripheral Vascular Clamp, 6" (15.2 cm), curved jaws



MH24-3200  
MH24-3210  
MH24-3220



MH24-3230  
MH24-3240  
MH24-3250

REF	Description
* MH24-3200	COOLEY Pediatric Vascular Clamp, 5-1/2" (14 cm), angled jaws 25 mm, straight handles/shank
* MH24-3210	COOLEY Pediatric Vascular Clamp, 5-1/2" (14 cm), straight jaws 25 mm, angled handles/shank
* MH24-3220	COOLEY Pediatric Vascular Clamp, 5-1/2" (14 cm), straight jaws 25 mm, straight handles/shank
* MH24-3230	COOLEY Pediatric Vascular Clamp, 4-3/4" (12.1 cm), right angle jaws, 17 mm, straight handles/shank
* MH24-3240	COOLEY Pediatric Vascular Clamp, 4-3/4" (12.1 cm), right angle jaws, 25 mm, straight handles/shank
* MH24-3250	COOLEY Pediatric Vascular Clamp, 4-3/4" (12.1 cm), right angle jaws, 35 mm, straight handles/shank

\*5 mm calibrations on outer sides of jaws




MH24-3260



MH24-3270

REF	Description
MH24-3260	COOLEY-SATINSKY Clamp, 5-1/2" (14 cm)
MH24-3270	COOLEY Vascular Clamp, 5-1/2" (14 cm), spoon shape



REF	Description
MH5-128	MAYO Dissecting Scissors, 9" (22.9 cm), straight, standard beveled blades
MH5-130	MAYO Dissecting Scissors, 9" (22.9 cm), curved, standard beveled blades
MH6-142	RUSSIAN Tissue Forceps, 6" (15.2 cm)
MH6-144	RUSSIAN Tissue Forceps, 8" (20.3 cm)
MH6-145	RUSSIAN Tissue Forceps, 10" (25.4 cm)
MH6-154	POTTS-SMITH Dressing Forceps, 7" (17.8 cm), serrated
MH6-160	POTTS-SMITH Tissue Forceps, 7" (17.8 cm), 1 X 2 teeth, serrated tips



MH7-206



MH7-255 to  
MH7-259

REF	Description
MH7-206	MIXTER Forceps, 9" (22.9 cm), fully curved
MH7-255	GEMINI-MIXTER Hemostatic Forceps, 5-1/2" (14 cm), delicate, full curved jaws
MH7-256	GEMINI-MIXTER Hemostatic Forceps, 7" (17.8 cm), delicate, full curved jaws
MH7-257	GEMINI-MIXTER Hemostatic Forceps, 8" (20.3 cm), delicate, full curved jaws
MH7-258	GEMINI-MIXTER Hemostatic Forceps, 9" (22.9 cm), delicate, full curved jaws
MH7-259	GEMINI-MIXTER Hemostatic Forceps, 11" (27.9 cm), delicate, full curved jaws



MH25-800

MH25-814

MH25-836  
MH25-837


MH25-1160

REF	Description
MH25-800	SAROT Intra-Thoracic Artery Forceps, 10" (25.4 cm), with 2-1/2" (6.4 cm) long jaws
MH25-814	MIXTER Thoracic Forceps, 11" (27.9 cm), right angle jaws, longitudinal serrations
MH25-836	KANTROWITZ Thoracic Forceps, 7-1/2" (19.1 cm), delicate right angle jaws
MH25-837	KANTROWITZ Thoracic Forceps, 9-1/2" (24.1 cm), delicate right angle jaws
MH25-1160	CRAFOORD Coarctation Clamp, 9-1/2" (24.1 cm), slightly curved jaws, 65 mm with longitudinal serrations





MH8-51TC  
MH8-52TC



MH8-101TC to  
MH8-105TC



MH11-246 to  
MH11-249

REF	Description
MH8-51TC	CRILE-WOOD Needle Holder, 7" (17.8 cm), Tungsten Carbide jaws
MH8-52TC	CRILE-WOOD Needle Holder, 8" (20.3 cm), Tungsten Carbide jaws
MH8-101TC	RYDER Needle Holder, 5" (12.7 cm), serrated jaws, Tungsten Carbide jaws
MH8-103TC	RYDER Needle Holder, 6" (15.2 cm), serrated jaws, Tungsten Carbide jaws
MH8-104TC	RYDER Needle Holder, 7" (17.8 cm), serrated jaws, Tungsten Carbide jaws
MH8-105TC	RYDER Needle Holder, 8" (20.3 cm), serrated jaws, Tungsten Carbide jaws
MH11-246	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1" (2.5 cm) X 3/4" (1.9 cm)
MH11-247	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1-1/4" (3.2 cm) X 1" (2.5 cm)
MH11-248	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 1-1/2" (3.8 cm) X 1-1/2" (3.8 cm)
MH11-249	RICHARDSON Retractor, 10" (25.4 cm), loop handle, 2" (5.1 cm) X 3/4" (1.9 cm)

\*For easy identification, MeisterHand® Tungsten Carbide Needle Holders have two gold-plated finger rings.



MH11-314  
MH11-316



MH11-590



MH16-12

REF	Description
MH11-314	HARRINGTON Retractor, 12-1/2" (31.8 cm), blade 2-1/2" (6.4 cm) wide
MH11-316	HARRINGTON Retractor, 12-1/2" (31.8 cm), blade 1-1/2" (3.8 cm) wide
MH11-590	HIBBS Retractor, 9-1/4" (23.5 cm), blade 1" (2.5 cm) X 3" (7.6 cm) deep
MH16-12	ALLIS Tissue Forceps, 10" (25.4 cm), 5 X 6 teeth



MH16-46  
MH16-48



MH19-850  
MH19-852



MH25-384

REF	Description
MH16-46	BABCOCK Intestinal Forceps, 8-1/4" (21 cm), jaws 10 mm wide
MH16-48	BABCOCK Intestinal Forceps, 9-1/2" (24.1 cm), jaws 15 mm wide
MH19-850	RUSKIN Rongeur, 7-1/4" (18.4 cm) long, double action, straight 5 mm bite jaws
MH19-852	RUSKIN Rongeur, 7-1/4" (18.4 cm) long, double action, curved 5 mm bite jaws
MH25-384	RUSKIN Bone Cutting Forceps, 7-1/2" (19.1 cm), straight standard blades



MH25-390

MH25-392

MH25-396

MH25-454

MH25-456

MH25-458

REF	Description
MH25-390	STILLE-LISTON Bone Cutting Forceps, 10-3/4" (27.3 cm), straight
MH25-392	STILLE-LISTON Bone Cutting Forceps, 10-3/4" (27.3 cm), curved
MH25-396	STILLE-HORSLEY Bone Cutting Forceps, 10-1/2" (26.7 cm), angled blades
MH25-454	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), straight jaws 9 X 15 mm
MH25-456	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), curved jaws 9 X 15 mm
MH25-458	STILLE-LUER Duckbill Rongeur, 9" (22.9 cm), angled on side, jaws 6 X 12 mm



MH25-480



MH25-510

MH25-511

REF	Description
MH25-480	SAUERBRUCH Rongeur, 12" (30.5 cm), straight, square jaws 17 mm wide
MH25-510	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 3 mm bite
MH25-511	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 8 mm bite




MH6-104




MH6-106

REF	Description
MH6-104	SEMKEN Dressing Forceps, 5" (12.7 cm), serrated
MH6-106	SEMKEN Tissue Forceps, 5" (12.7 cm), 1 X 2 teeth


**MH6-118**
**MH6-120**  
**MH6-122**
**MH6-121**
**MH6-123**
**MH6-124**
**MH17-2500**  
**MH17-2510**

REF	Description
MH6-118	ADSON Dressing Forceps, 4-3/4" (12.1 cm), delicate, serrated
MH6-120	ADSON Tissue Forceps 4-3/4" (12.1 cm), 1 X 2 teeth, delicate
MH6-121	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, cross serrated tips, delicate
MH6-122	ADSON Tissue Forceps, 4-3/4" (12.1 cm), 2 X 3 teeth, delicate
MH6-123	ADSON Tissue & Suture Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, with tying platform
MH6-124	BROWN-ADSON Tissue Forceps, 4-3/4" (12.1 cm), 7 X 7 side grasping teeth
MH17-2500	Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), 1 X 2 teeth, tips 0.5 mm wide
MH17-2510	Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated tips, 0.5 mm wide



MH26-956  
MH26-958



MH6-180



MH6-182



MH6-190



MH6-192

REF	Description
MH26-956	FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 1 X 2 teeth
MH26-958	FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 2 X 3 teeth
MH6-180	GERALD Dressing Forceps, 7" (17.8 cm), straight, serrated
MH6-182	GERALD Tissue Forceps, 7" (17.8 cm), 1 X 2 teeth, straight
MH6-190	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), bayonet, serrated, scraper end
MH6-192	CUSHING Dressing Forceps, 7-1/4" (18.4 cm), 1 X 2 teeth, bayonet, scraper end





MH6-196



MH6-197



MH6-198

REF	Description
MH6-196	GRUENWALD (ADSON) Forceps, 8-1/4" (21 cm), serrated
MH6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps, 6-1/2" (16.5 cm), serrated
MH6-198	ADSON Bayonet Dressing Forceps, 8-1/4" (21 cm), serrated, delicate tips



MH7-8  
MH7-10



MH7-102

REF	Description
MH7-8	HALSTED Mosquito Forceps, 5" (12.7 cm), straight, extra delicate
MH7-10	HALSTED Mosquito Forceps, 5" (12.7 cm), curved, extra delicate
MH7-102	DANDY Forceps, 5-1/2" (14 cm), curved sideways



MH7-240  
MH7-242



MH17-2600  
MH17-2602

REF	Description
MH7-240	ADSON Hemostatic Forceps, 7-1/4" (18.4 cm), straight
MH7-242	ADSON Hemostatic Forceps, 7-1/4" (18.4 cm), curved
MH17-2600	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), straight, extremely delicate
MH17-2602	JACOBSON Micro Mosquito Forceps, 5" (12.7 cm), curved, extremely delicate



MH11-2



MH11-102

MH11-240 to  
MH11-244

REF	Description
MH11-2	JANSEN Retractor, 4-1/4" (10.8 cm), spread 2" (5.1 cm), 3 X 3 blunt prongs 10 mm deep, standard model
MH11-102	CUSHING Vein Retractor, 9" (22.9 cm), blade 11 X 13 mm
MH11-240	MEYERDING Retractor, 9-1/2" (24.1 cm), large blade, 3-1/2" (8.9 cm) X 2" (5.1 cm)
MH11-242	MEYERDING Retractor, 9-1/2" (24.1 cm), medium blade, 3" (7.6 cm) X 1" (2.5 cm)
MH11-244	MEYERDING Retractor, 9-1/2" (24.1 cm), small blade, 2" (5.1 cm) X 5/8" (1.6 cm)



**MH11-364 to  
MH11-366**



**MH11-630-BL to  
MH11-630-SH**



**MH11-684**



**MH11-700  
MH11-702**

REF	Description
MH11-364	GELPI Retractor, 3-1/2" (8.9 cm), with grip lock, Mini
MH11-365	GELPI Retractor, Pediatric, 5-1/2" (14 cm), with grip lock
MH11-366	GELPI Retractor, 7-1/2" (19.1 cm), with grip lock
MH11-630-BL	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth, with hinged blades
MH11-630-SH	BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth, with hinged blades
MH11-684	BECKMAN-ADSON Retractor, 12-1/2" (31.8 cm), hinged, sharp, 4 X 4 teeth, 3/4" (1.9 cm) X 1" (2.5 cm)
MH11-700	ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth, 3/4" (1.9 cm) X 3/4" (1.9 cm), straight arms
MH11-702	ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth, 3/4" (1.9 cm) X 3/4" (1.9 cm), angled arms




MH19-570 to  
MH19-577



MH19-580 to  
MH19-584

REF	Description
MH19-570	FRAZIER Suction Tube, 6 French, (2 mm), with finger cut-off, angled
MH19-572	FRAZIER Suction Tube, 7 French, (2.3 mm), with finger cut-off, angled
MH19-573	FRAZIER Suction Tube, 8 French, (2.6 mm), with finger cut-off, angled
MH19-574	FRAZIER Suction Tube, 9 French, (3 mm), with finger cut-off, angled
MH19-575	FRAZIER Suction Tube, 10 French, (3.3 mm), with finger cut-off, angled
MH19-576	FRAZIER Suction Tube, 11 French, (3.6 mm), with finger cut-off, angled
MH19-577	FRAZIER Suction Tube, 12 French, (4 mm), with finger cut-off, angled
MH19-580	BARON Suction Tube, 3 French (1 mm), with finger cut-off, working length 7.5 cm
MH19-582	BARON Suction Tube, 5 French (1.7 mm), with finger cut-off, working length 7.5 cm
MH19-584	BARON Suction Tube, 7 French (2.3 mm), with finger cut-off, working length 7.5 cm




MH19-700 to  
MH19-720

**SPRATT (BRUN) MASTOID CURETTES, OVAL CUPS**

REF	Description
MH19-700	6-1/2" (16.5 cm), size 5/0
MH19-702	6-1/2" (16.5 cm), size 4/0
MH19-704	6-1/2" (16.5 cm), size 3/0
MH19-706	6-1/2" (16.5 cm), size 2/0
MH19-708	6-1/2" (16.5 cm), size 0
MH19-710	6-1/2" (16.5 cm), size 1

REF	Description
MH19-712	6-1/2" (16.5 cm), size 2
MH19-714	6-1/2" (16.5 cm), size 3
MH19-716	6-1/2" (16.5 cm), size 4
MH19-718	6-1/2" (16.5 cm), size 5
MH19-720	6-1/2" (16.5 cm), size 6



MH26-1627 to  
MH26-1640

### HIBBS SPRATT SPINAL FUSION CURETTES, OVAL CUPS

REF	Description
MH26-1627	9" (22.9 cm), size 3/0
MH26-1628	9" (22.9 cm), size 2/0
MH26-1629	9" (22.9 cm), size 0
MH26-1630	9" (22.9 cm), size 1
MH26-1632	9" (22.9 cm), size 2

REF	Description
MH26-1634	9" (22.9 cm), size 3
MH26-1636	9" (22.9 cm), size 4
MH26-1638	9" (22.9 cm), size 5
MH26-1640	9" (22.9 cm), size 6





MH19-844



MH19-848  
MH19-849



MH19-850  
MH19-852

REF	Description
MH19-844	BEYER Rongeur, 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide
MH19-848	RUSKIN Rongeur, 6" (15.2 cm), double action, straight 3 mm bite jaws
MH19-849	RUSKIN Rongeur, 6" (15.2 cm), double action, curved 3 mm bite jaws
MH19-850	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, straight 5 mm bite jaws
MH19-852	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, curved 5 mm bite jaws



MH25-454  
MH25-456



MH25-458



MH25-510  
MH25-511

REF	Description
MH25-454	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), straight jaws 9 X 15 mm
MH25-456	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), curved jaws 9 X 15 mm
MH25-458	STILLE-LUER Duckbill Rongeur, 9" (22.9 cm), angled on side, jaws 6 X 12 mm
MH25-510	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 3 mm bite
MH25-511	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 8 mm bite



MH26-119



MH26-124



MH26-336  
MH26-339A

REF	Description
MH26-119	GIGLI Saw Handle, standard pattern, T-shaped, turnable
MH26-124	GIGLI-STRULLY Saw Handle, standard pattern, loop grip, rigid, with snap-lock
MH26-336	SPURLING-KERRISON Cervical Rongeur, 8" (20.3 cm) shaft, up-bite, 3 mm
MH26-339A	SPURLING-KERRISON Cervical Rongeur, 8" (20.3 cm) shaft, 40 degree forward angle, 3 mm



MH26-394  
MH26-397



MH26-400  
MH26-402



MH26-430  
MH26-435

REF	Description
MH26-394	CUSHING Pituitary Rongeur, 5" (12.7 cm) shank, straight, 2 X 10 mm cup jaws
MH26-397	CUSHING Pituitary Rongeur, 7" (17.8 cm) shank, straight, 2 X 10 mm cup jaws
MH26-400	LOVE-GRUENWALD Pituitary Rongeur, 5" (12.7 cm) shank, straight, 3 X 10 mm cup jaws
MH26-402	LOVE-GRUENWALD Pituitary Rongeur, 7" (17.8 cm) shank, straight, 3 X 10 mm cup jaws
MH26-430	SPURLING Pituitary Rongeur, 5" (12.7 cm) shank, straight, 4 X 10 mm cup jaws
MH26-435	SPURLING Pituitary Rongeur, 7" (17.8 cm) shank, straight, 4 X 10 mm cup jaws



MH26-1084



MH26-1090




MH26-1410



MH26-1412

REF	Description
MH26-1084	FRAZIER Dura Hook, 5" (12.7 cm), sharp point
MH26-1090	ADSON Dura and Dissecting Hook, 8" (20.3 cm), sharp blade 5 mm at right angle
MH26-1410	LANGENBECK Periosteal Elevator, 7-1/2" (19.1 cm), narrow model, 7 mm wide blunt blade
MH26-1412	LANGENBECK Periosteal Elevator, 7-1/2" (19.1 cm), wide model, 17 mm wide sharp blade



MH18-1968


MH26-1450

MH26-1451  
MH26-1452

MH26-1453

MH26-1454

REF	Description
MH18-1968	FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
MH26-1450	PENFIELD Dissector, style No. 1, 7-1/4" (18.4 cm), double end, 6 mm sharp cup and 6 mm blunt dissector
MH26-1451	PENFIELD Dissector, style No. 2, 7-3/4" (19.7 cm), double end, wax packer and 6 mm blunt dissector, slightly curved
MH26-1452	PENFIELD Dissector, style No. 3, 7-1/2" (19.1 cm), double end, wax packer and 6 mm blunt dissector, fully curved
MH26-1453	PENFIELD Dissector, style No. 4, 8" (20.3 cm), single end, light pattern, 3 mm blunt dissector, slightly curved
MH26-1454	PENFIELD Dissector, style No. 5, 11-1/2" (29.2 cm), double end, blunt dissector blades 7 mm and 8 mm wide, slightly curved



MH11-530 to  
MH11-544



MH11-580  
MH11-582

**VOLKMAN RETRACTORS, 8-1/2" (21.6 cm) ring handle**

Sharp	Description	Blunt
MH11-530	2 prongs	MH11-538
MH11-532	3 prongs	MH11-540
MH11-534	4 prongs	MH11-542
MH11-536	6 prongs	MH11-544

REF	Description
MH11-580	BENNETT Tibia Retractor, 9-1/2" (24.1 cm), small blade 1-3/4" (4.4 cm) wide
MH11-582	BENNETT Tibia Retractor, 9-1/2" (24.1 cm), large blade 2-1/2" (6.4 cm) wide




MH27-213  
MH27-214



MH27-215  
MH27-219

REF	Description
MH27-213	MINI-HOHMANN Retractor, 6-1/4" (15.9 cm), blade 6 mm wide
MH27-214	MINI-HOHMANN Retractor, 6-1/4" (15.9 cm), blade 8 mm wide
MH27-215	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 10 mm wide
MH27-216	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 18 mm wide
MH27-217	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 25 mm wide
MH27-218	HOHMANN Retractor, 9-1/4" (23.5 cm), blade 43 mm wide
MH27-219	HOHMANN Retractor, 9-3/4" (24.8 cm), blade 70 mm wide






MH19-700 to  
MH19-720

**SPRATT (BRUN) MASTOID CURETTES, OVAL CUPS**

REF	Description
MH19-700	6-1/2" (16.5 cm), size 5/0
MH19-702	6-1/2" (16.5 cm), size 4/0
MH19-704	6-1/2" (16.5 cm), size 3/0
MH19-706	6-1/2" (16.5 cm), size 2/0
MH19-708	6-1/2" (16.5 cm), size 0
MH19-710	6-1/2" (16.5 cm), size 1

REF	Description
MH19-712	6-1/2" (16.5 cm), size 2
MH19-714	6-1/2" (16.5 cm), size 3
MH19-716	6-1/2" (16.5 cm), size 4
MH19-718	6-1/2" (16.5 cm), size 5
MH19-720	6-1/2" (16.5 cm), size 6



MH26-1627 to  
MH26-1640

**HIBBS SPRATT SPINAL FUSION CURETTES, OVAL CUPS**

REF	Description
MH26-1627	9" (22.9 cm), size 3/0
MH26-1628	9" (22.9 cm), size 2/0
MH26-1629	9" (22.9 cm), size 0
MH26-1630	9" (22.9 cm), size 1
MH26-1632	9" (22.9 cm), size 2

REF	Description
MH26-1634	9" (22.9 cm), size 3
MH26-1636	9" (22.9 cm), size 4
MH26-1638	9" (22.9 cm), size 5
MH26-1640	9" (22.9 cm), size 6



MH19-844



MH19-848

MH19-849



MH19-850

MH19-852

REF	Description
MH19-844	BEYER Rongeur, 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide
MH19-848	RUSKIN Rongeur, 6" (15.2 cm), double action, straight 3 mm bite jaws
MH19-849	RUSKIN Rongeur, 6" (15.2 cm), double action, curved 3 mm bite jaws
MH19-850	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, straight 5 mm bite jaws
MH19-852	RUSKIN Rongeur, 7-1/4" (18.4 cm), double action, curved 5 mm bite jaws



MH21-626  
MH21-627



MH25-384



MH25-390



MH25-392

REF	Description
MH21-626	KLEINERT-KUTZ (RUSKIN) Bone Splitting Forceps, 6" (15.2 cm), straight
MH21-627	KLEINERT-KUTZ (RUSKIN) Bone Splitting Forceps, 6" (15.2 cm), curved
MH25-384	RUSKIN Bone Cutting Forceps, 7-1/2" (19.1 cm), straight standard blades
MH25-390	STILLE-LISTON Bone Cutting Forceps, 10-3/4" (27.3 cm), straight
MH25-392	STILLE-LISTON Bone Cutting Forceps, 10-3/4" (27.3 cm), curved


**MH25-396**
**MH25-454  
MH25-456**
**MH25-458**
**MH25-510 to  
MH25-511**

REF	Description
MH25-396	STILLE-HORSLEY Bone Cutting Forceps, 10-1/2" (26.7 cm), angled blades
MH25-454	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), straight jaws 9 X 15 mm
MH25-456	STILLE-LUER Rongeur, 8-3/4" (22.2 cm), curved jaws 9 X 15 mm
MH25-458	STILLE-LUER Duckbill Rongeur, 9" (22.9 cm), angled on side, jaws 6 X 12 mm
MH25-510	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 3 mm bite
MH25-511	LEKSELL Rongeur, 9-1/2" (24.1 cm), standard curve, 8 mm bite



MH26-119



MH26-124



MH27-4  
MH27-7



MH27-16

REF	Description
MH26-119	GIGLI Saw Handle, standard pattern, T-shaped, turnable
MH26-124	GIGLI-STRULLY Saw Handle, standard pattern, loop grip, rigid, with snap-lock
MH27-4	KERN Bone Holding Forceps, 9-1/2" (24.1 cm), standard size, with ratchet
MH27-7	KERN Bone Holding Forceps, 6" (15.2 cm), pediatric size, with ratchet
MH27-16	LEWIN Bone Holding Forceps, 7" (17.8 cm), serrated jaws



MH27-30 to  
MH27-34



MH27-160



MH27-3070



MH27-3100

REF	Description
MH27-30	LOWMAN Bone Clamp, 5" (12.7 cm), calibrated in 1/8" (.32 cm), 1 X 2 prong jaws 3/4" (1.9 cm) wide
MH27-32	LOWMAN Bone Clamp, 7" (17.8 cm), calibrated in 1/8" (.32 cm), 1 X 2 prong jaws 1" (2.5 cm) wide
MH27-34	LOWMAN Bone Clamp, 8" (20.3 cm), calibrated in 1/8" (.32 cm), 1 X 2 prong jaws 1-1/4" (3.2 cm) wide
MH27-160	DIAMOND Pin & Wire Cutter, 6-1/2" (16.5 cm), cuts up to 2 mm (5/64") diameter hard wire, stainless steel
MH27-3070	BRUNS Plaster Shears, 9-1/4" (23.5 cm), serrated blade
MH27-3100	Three Prong Cast Spreader, 9" (22.9 cm), with spring action for one-hand operation



MH27-438 to  
MH27-448



MH27-438C to  
MH27-448C



MH27-488 to  
MH27-498



MH27-499 to  
MH27-508

**HIBBS OSTEOTOMES 9" (22.9 cm)**

Straight	Size	Curved
MH27-438	1/4" (.64 cm)	MH27-438C
MH27-440	3/8" (.95 cm)	MH27-440C
MH27-442	1/2" (1.3 cm)	MH27-442C
MH27-444	5/8" (1.6 cm)	MH27-444C
MH27-446	3/4" (1.9 cm)	MH27-446C
MH27-447	7/8" (2.2 cm)	MH27-447C
MH27-448	1" (2.5 cm)	MH27-448C

**LAMBOTTE OSTEOTOMES 9" (22.9 cm)**

Straight	Width	Curved
MH27-488	6 mm (1/4")	MH27-499
MH27-490	13 mm (1/2")	MH27-500
MH27-492	19 mm (3/4")	MH27-502
MH27-494	25 mm (1")	MH27-504
MH27-496	32 mm (1-1/4")	MH27-506
MH27-498	38 mm (1-1/2")	MH27-508





MH19-793



MH19-796




MH26-1410



MH26-1412

REF	Description
MH19-793	Mallet, 7-1/2" (19.1 cm), lightweight 7.5 oz. (213 g), stainless steel head 1" (2.5 cm) diameter with 2 replaceable nylon caps
MH19-796	MEAD Mallet, 7-1/2" (19.1 cm), head 8 oz. (227 g), with 2 replaceable nylon caps
MH26-1410	LANGENBECK Periosteal Elevator, 7-1/2" (19.1 cm), narrow model, 7 mm wide blunt blade
MH26-1412	LANGENBECK Periosteal Elevator, 7-1/2" (19.1 cm), wide model, 17 mm wide sharp blade



MH27-820 to  
MH27-826



MH27-700



MH27-702

REF	Description
MH27-820	KEY Periosteal Elevator, 7" (17.8 cm), 1/4" (.64 cm) wide
MH27-822	KEY Periosteal Elevator, 7-1/2" (19.1 cm), 1/2" (1.3 cm) wide
MH27-824	KEY Periosteal Elevator, 8" (20.3 cm), 3/4" (19. cm) wide
MH27-826	KEY Periosteal Elevator, 8-1/4" (21 cm), 1" (2.5 cm) wide
MH27-700	PUTTI Bone Rasp, 10-1/2" (26.7 cm), flat ends, one side curved on flat, other side curved sideways
MH27-702	PUTTI Bone Rasp, 12" (30.5 cm), round and half-round tapered curved ends



MH27-990



MH27-1000



MH27-1016



MH7-106



MH7-220  
MH7-221



MH16-10



MH16-30

REF	Description
MH27-990	MARTIN Cartilage Clamp, 7" (17.8 cm)
MH27-1000	MARTIN Cartilage Scissors, 8" (20.3 cm), curved serrated non-slip cutting edges
MH27-1016	Tendon Pulling Forceps, 5" (12.7 cm), curved shaft, alligator type
MH7-106	PRATT T-Shape Forceps, 6" (15.2 cm)
MH7-220	PENNINGTON Tissue Grasping Forceps, 6" (15.2 cm), triangular jaws
MH7-221	PENNINGTON Tissue Grasping Forceps, 8" (20.3 cm), triangular jaws
MH16-10	ALLIS Tissue Forceps, 7-1/4" (18.4 cm), 5 X 6 teeth, straight
MH16-30	ALLIS-ADAIR Tissue Forceps, 6-1/2" (16.5 cm), 9 X 10 teeth



MH40-210-SS



MH40-212-SS



MH40-224



MH40-225-SS



MH40-226A



MH40-226-SS



MH40-227-SS



MH40-230



MH40-245-SS



MH40-250-SS

REF	Description
MH40-210-SS	Nail Nipper, 5-1/2" (14.0 cm), concave, stainless steel
MH40-212-SS	Nail Nipper, 5-1/2" (14.0 cm), straight, stainless steel
MH40-224	Nail Nipper, 4" (10.2 cm), straight, stainless steel
MH40-225-SS	Nail Nipper, 4-1/2" (11.4 cm), straight, stainless steel
MH40-226A	Nail Nipper, 5" (12.7 cm), delicate, straight, stainless steel
MH40-226-SS	Nail Nipper, 5" (12.7 cm), straight, stainless steel
MH40-227-SS	Nail Nipper, 6" (15.2 cm), heavy, straight, stainless steel
MH40-230	Nail Splitter, 5" (12.7 cm), english anvil, stainless steel
MH40-245-SS	Tissue & Cuticle Nipper, 4" (10.2 cm), convex, stainless steel
MH40-250-SS	Tissue & Cuticle Nipper, 4-1/2" (11.4 cm), stainless steel



MH16-38



MH28-196



MH28-1



MH28-2

REF	Description
MH16-38	THOMS-ALLIS Tissue Forceps, 8" (20.3 cm), 6 X 7 teeth
MH28-196	BUIE Pile Clamp, 8-1/2" (21.6 cm), curved angiotribe jaws
MH28-1	PRATT Rectal Speculum, 8-1/2" (21.6 cm) long, with set screw, blades 3-1/2" (8.9 cm) long X 1" (2.5 cm) wide at distal end
MH28-2	SIMS Rectal Speculum, 6" (15.2 cm) long, with set screw, fenestrated blades, 3-1/2" (8.9 cm) long X 5/8" (1.6 cm) wide at distal end



MH28-20 to  
MH28-24



MH28-68




MH28-69

REF	Description
MH28-20	HIRSCHMAN Anoscope, small size, 2-1/4" (5.7 cm) long X 9/16" (1.4 cm) I.D.
MH28-22	HIRSCHMAN Anoscope, medium size, 2-5/8" (6.7 cm) long X 11/16" (1.7 cm) I.D.
MH28-24	HIRSCHMAN Anoscope, large size, 3-1/2" (8.9 cm) long X 7/8" (2.2 cm) I.D.
MH28-68	HILL-FERGUSON Retractor, 8-1/2" (21.6 cm), small blade 2.2 X 6.4 cm
MH28-69	HILL-FERGUSON Retractor, 8-3/4" (22.2 cm), large blade 3 X 7.6 cm




MH29-6-8 to MH29-6-36

REF	Description
MH29-6-8	VAN BUREN Urethral Sound, 11" (27.9 cm), 8 Fr. (2.6 mm)
MH29-6-10	VAN BUREN Urethral Sound, 11" (27.9 cm), 10 Fr. (3.3 mm)
MH29-6-12	VAN BUREN Urethral Sound, 11" (27.9 cm), 12 Fr. (4 mm)
MH29-6-14	VAN BUREN Urethral Sound, 11" (27.9 cm), 14 Fr. (4.6 mm)
MH29-6-16	VAN BUREN Urethral Sound, 11" (27.9 cm), 16 Fr. (5.3 mm)
MH29-6-18	VAN BUREN Urethral Sound, 11" (27.9 cm), 18 Fr. (6 mm)
MH29-6-20	VAN BUREN Urethral Sound, 11" (27.9 cm), 20 Fr. (6.6 mm)
MH29-6-22	VAN BUREN Urethral Sound, 11" (27.9 cm), 22 Fr. (7.3 mm)
MH29-6-24	VAN BUREN Urethral Sound, 11" (27.9 cm), 24 Fr. (8 mm)
MH29-6-26	VAN BUREN Urethral Sound, 11" (27.9 cm), 26 Fr. (8.6 mm)
MH29-6-28	VAN BUREN Urethral Sound, 11" (27.9 cm), 28 Fr. (9.2 mm)
MH29-6-30	VAN BUREN Urethral Sound, 11" (27.9 cm), 30 Fr. (10 mm)
MH29-6-32	VAN BUREN Urethral Sound, 11" (27.9 cm), 32 Fr. (10.6 mm)
MH29-6-34	VAN BUREN Urethral Sound, 11" (27.9 cm), 34 Fr. (11.2 mm)
MH29-6-36	VAN BUREN Urethral Sound, 11" (27.9 cm), 36 Fr. (11.9 mm)



REF	Description
MH29-284	RANDALL Kidney Stone Forceps, 9-1/4" (23.5 cm), quarter curved
MH29-286	RANDALL Kidney Stone Forceps, 8-1/4" (21 cm), half curved
MH29-288	RANDALL Kidney Stone Forceps, 7-3/4" (19.7 cm), three quarter curved
MH29-290	RANDALL Kidney Stone Forceps, 7-1/2" (19.1 cm), full curved
MH29-302	STILLE Pattern Vessel Clamp, 9" (22.9 cm)
MH29-308	MAYO (GUYON) Vessel Clamp, 9-1/4" (23.5 cm), double curved jaws, serrated





MH7-124 to  
MH7-144

MH7-600  
MH7-602

MH7-622

REF	Description
MH7-124	ROCHESTER-PEAN Forceps, 8" (20.3 cm), straight
MH7-126	ROCHESTER-PEAN Forceps, 9" (22.9 cm), straight
MH7-142	ROCHESTER-PEAN Forceps, 8" (20.3 cm), curved
MH7-144	ROCHESTER-PEAN Forceps, 9" (22.9 cm), curved
MH7-600	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), straight, serrated jaws
MH7-602	FOERSTER Sponge Forceps, 9-1/2" (24.1 cm), straight, smooth jaws
MH7-622	BOZEMAN Sponge Forceps, 10-1/2" (26.7 cm), double curved



MH11-180

MH11-182

MH11-184

MH11-185

MH30-395 to  
MH30-415

REF	Description
MH11-180	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1" (2.5 cm) X 1-3/4" (4.4 cm)
MH11-182	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1-1/8" (2.8 cm) X 2-7/8" (7.3 cm)
MH11-184	HEANEY Retractor, 10-1/2" (26.7 cm), blade 1-1/8" (2.8 cm) X 3-1/4" (8.3 cm)
MH11-185	HEANEY-SIMON Retractor, 11-1/2" (29.2 cm), blade 1" (2.5 cm) X 4-1/2" (11.4 cm)
MH30-395	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 1-3/4" (4.4 cm) X 2-1/2" (6.4 cm)
MH30-406	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 2-1/8" (5.4 cm) X 3" (7.6 cm)
MH30-415	DOYEN Vaginal Retractor, 9" (22.9 cm), blade 2-1/4" (5.7 cm) X 3-1/2" (8.9 cm)



MH11-364 to  
MH11-366

MH12-80

MH16-38

MH16-51

REF	Description
MH11-364	GELPI Retractor, 3-1/2" (8.9 cm), with grip lock, Mini
MH11-365	GELPI Retractor, Pediatric, 5-1/2" (14 cm), with grip lock
MH11-366	GELPI Retractor, 7-1/2" (19.1 cm), with grip lock
MH12-80	LAHEY Traction Forceps, 6-1/4" (15.9 cm), 3 X 3 teeth
MH16-38	THOMS-ALLIS Tissue Forceps, 8" (20.3 cm), 6 X 7 teeth
MH16-51	ADAIR Tissue Holding Tenaculum Forceps, 6-1/2" (16.5 cm)



MH30-10 to  
MH30-20



MH30-50 to  
MH30-62

REF	Description
MH30-10	GRAVES Vaginal Speculum, original model, small size, 3/4" (1.9 cm) X 3" (7.6 cm)
MH30-15	GRAVES Vaginal Speculum, original model, medium size, 1-3/8" (3.5 cm) X 4" (10.2 cm)
MH30-20	GRAVES Vaginal Speculum, original model, large size, 1-1/2" (3.8 cm) X 4-1/2" (11.4 cm)
MH30-30	GRAVES Open Sided Vaginal Speculum, medium size, 1-3/8" (3.5 cm) X 4" (10.2 cm), wide angle blades
MH30-50	PEDERSON Vaginal Speculum, narrow blades, small size, 1/2" (1.3 cm) X 3" (7.6 cm)
MH30-55	PEDERSON Vaginal Speculum, narrow blades, medium size, 7/8" (2.2 cm) X 4" (10.2 cm)
MH30-60	PEDERSON Vaginal Speculum, narrow blades, large size, 1" (2.5 cm) X 4-3/4" (12.1 cm)
MH30-62	PEDERSON Vaginal Speculum, extra narrow blades, 5/8" (1.6 cm) X 4" (10.2 cm)



MH30-500



MH30-505-910 to MH30-505-1920

REF	Description
MH30-500	HANK Uterine Dilators, 10-1/2" (26.7 cm), double end, set of 6 Dilators in Khaki Roll, Sizes 9-10, (3-3.3 mm), 11-12, (3.6-4 mm), 13-14 (4.3-4.6 mm), 15-16 (5-5.3 mm), 17-18 (5.6-6 mm) and 19-20 (6.3-6.6 mm) French
MH30-505-910	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 9-10 Fr. (3-3.3 mm)
MH30-505-1112	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 11-12 Fr. (3.6-4 mm)
MH30-505-1314	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 13-14 Fr. (4.3-4.6 mm)
MH30-505-1516	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 15-16 Fr. (5-5.3 mm)
MH30-505-1718	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 17-18 Fr. (5.6-6 mm)
MH30-505-1920	HANK Uterine Dilator, 10-1/2" (26.7 cm), double end, 19-20 Fr. (6.3-6.6 mm)



MH30-550

MH30-555-34 to  
MH30-555-1718

REF	Description
MH30-550	HEGAR Uterine Dilators, 7-3/4" (19.7 cm), double end, set of 8 Dilators in Khaki Roll, sizes 3-4 mm to 17-18 mm
MH30-555-34	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 3-4 mm
MH30-555-56	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 5-6 mm
MH30-555-78	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 7-8 mm
MH30-555-910	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 9-10 mm
MH30-555-1112	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 11-12 mm
MH30-555-1314	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 13-14 mm
MH30-555-1516	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 15-16 mm
MH30-555-1718	HEGAR Uterine Dilator, 7-3/4" (19.7 cm), double end, 17-18 mm



MH30-561




MH30-560-1315 to  
MH30-560-4143

**PRATT UTERINE DILATORS, DOUBLE END, 11-1/2" (29.2 cm)**

REF	Description
MH30-561	Set of 8 Dilators in Khaki Roll, sizes 13-15 Fr. (4.3 to 4.6 mm) to 41-43 Fr. (13.5-14.2 mm)
MH30-560-1315	13-15 Fr. (4.3-4.6 mm)
MH30-560-1719	17-19 Fr. (5.6-6.3 mm)
MH30-560-2123	21-23 Fr. (7-7.6 mm)

REF	Description
MH30-560-2527	25-27 Fr. (8.3-9 mm)
MH30-560-2931	29-31 Fr. (9.6-10.2 mm)
MH30-560-3335	33-35 Fr. (10.9-11.6 mm)
MH30-560-3739	37-39 Fr. (12.2-12.9 mm)
MH30-560-4143	41-43 Fr. (13.5-14.2 mm)



REF	Description
MH30-965	SCHROEDER-BRAUN Uterine Tenaculum Forceps, 9-3/4" (24.8 cm), straight
MH30-966	SCHROEDER Uterine Tenaculum Forceps, 10" (25.4 cm), round jaw
MH30-967	HULKA Tenaculum Forceps, 11" (27.9 cm), with 3.5 mm probe
MH30-985	BARRETT Uterine Tenaculum Forceps, 7-1/4" (18.4 cm), straight
MH30-965ATR	BRAUN Tenaculum Forceps, 10" (25.4 cm), straight/squared jaw with non-overlapping atraumatic tips
MH30-966ATR	SCHROEDER Tenaculum Forceps, 10" (25.4 cm), rounded jaw with non-overlapping atraumatic tips





MH30-997




MH30-1000



MH30-1025

REF	Description
MH30-997	TEALE Vulsellum Forceps, 10" (25.4 cm), curved sideways, 3 X 4 teeth, 5 mm wide
MH30-1000	SCHROEDER Uterine Vulsellum Forceps, 10" (25.4 cm), straight
MH30-1025	JACOBS Uterine Vulsellum Forceps, 8-1/2" (21.6 cm), straight, serrated jaws, 2 X 2 teeth



MH30-1205-0 to  
MH30-1205-6

REF	Description
MH30-1205-0	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 0
MH30-1205-1	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 1
MH30-1205-2	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 2
MH30-1205-3	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 3
MH30-1205-4	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 4
MH30-1205-5	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 5
MH30-1205-6	SIMS Uterine Curette, 11" (27.9 cm), sharp blades on malleable shank, size 6



MH30-1382

MH30-1383

MH30-1385


MH30-1390

MH30-1395

REF	Description
MH30-1382	KEVORKIAN-YOUNGE Endocervical Biopsy Curette, 12" (30.5 cm), loop size 3 X 12 mm, with basket
MH30-1383	KEVORKIAN-YOUNGE Endocervical Biopsy Curette, 12" (30.5 cm), loop size 3 X 12 mm, without basket
MH30-1385	HEANEY Uterine Biopsy Curette, 9" (22.9 cm), serrated loop
MH30-1390	NOVAK Endometrial Suction Biopsy Curette, 9-3/4" (24.8 cm), with luer hub, 4 mm diameter, standard 12 Fr. size
MH30-1395	RANDALL Endometrial Suction Biopsy Curette, 9-3/4" (24.8 cm), with luer hub, 4 mm diameter



MH30-1405




MH30-1442



MH30-1443

REF	Description
MH30-1405	THOMS-GAYLOR Uterine Biopsy Forceps, 8-1/2" (21.6 cm), angled shanks, curved jaws, 5 mm bite, interfitted cups
MH30-1442	TISCHLER Cervical Biopsy Punch Forceps, 8" (20.3 cm) shaft, 7 X 3 X 1.5 mm bite
MH30-1443	BURKE (Modified Baby TISCHLER) Cervical Biopsy Punch Forceps, 7-3/4" (19.7 cm) shaft, 4.5 X 3 X 1.5 mm bite



MH30-1482



MH30-1485

REF	Description
MH30-1482	KEVORKIAN-YOUNGE Uterine Biopsy Forceps, 9-1/2" (24.1 cm) shaft, 3.5 X 8 mm bite, teeth on lower jaw
MH30-1485	EPPENDORFER Uterine Biopsy Forceps, 8-1/2" (21.6 cm) shaft



MH30-1700  
MH30-1710



MH30-1716  
MH30-1718




MH30-2190  
MH30-2195

REF	Description
MH30-1700	HEANEY Hysterectomy Forceps, 8-1/4" (21 cm), heavy pattern, single tooth, curved
MH30-1710	HEANEY Hysterectomy Forceps, 8-1/4" (21 cm), heavy pattern, double tooth, curved
MH30-1716	HEANEY-BALLENTINE Hysterectomy Forceps, 8-1/2" (21 cm), longitudinal serrations, single tooth, curved
MH30-1718	BALLENTINE Hysterectomy Forceps, 8-1/2" (21.6 cm), longitudinal serrations, single tooth, straight
MH30-2190	BRAUN Episiotomy Scissors, 5-1/2" (14 cm), angled to side, guarded lower blade
MH30-2195	BRAUN Episiotomy Scissors, 8-1/2" (21.6 cm), angled to side, guarded lower blade



MH5-128



MH5-130



MH6-142 to  
MH6-145

REF	Description
MH5-128	MAYO Dissecting Scissors, 9" (22.9 cm), straight, standard beveled blades
MH5-130	MAYO Dissecting Scissors, 9" (22.9 cm), curved, standard beveled blades
MH6-142	RUSSIAN Tissue Forceps, 6" (15.2 cm)
MH6-144	RUSSIAN Tissue Forceps, 8" (20.3 cm)
MH6-145	RUSSIAN Tissue Forceps, 10" (25.4 cm)



MH1-226



MH3-10

REF	Description
MH1-226	WARTENBERG Neurological Pin Wheel, 7" (17.8 cm)
MH3-10	Three Prong Sterilizer Forceps, 8" (20.3 cm)





**A**

Adair Tissue Forceps	42, 43, 226
Adenoid Curettes	120
Adenotome Blades	120
Adenotomes	120
Adson Dressing Forceps	33, 34, 36, 44, 86, 108, 113, 133, 160, 190, 192
Adson Dura & Dissecting Hook	85, 204
Adson Hemostatic Forceps	50, 194
Adson Micro Dressing Forceps	34, 86, 133, 160, 190
Adson Micro Tissue Forceps	34, 86, 133, 160, 190
Adson Retractors	76, 196
Adson Tissue Forceps	33-34, 186, 133, 160, 190
Alligator Forceps	88, 109, 110, 114
Allis Tissue Forceps	41-42, 185, 218
Allis-Adair Tissue Forceps	42, 218
Alm Retractors	67, 104, 163
Andrews-Pynchon Suction Tube	65, 168
Anoscope	221
Ash Septum Forceps	151
Aufricht Glabella Rasps	150
Aufricht Retractor	78, 148

**B**

Babcock Forceps	43, 186
Backhaus Towel Clamps	55
Ballenger V-Shaped Chisels	118, 146
Ballentine Hysterectomy Forceps	54, 237
Bandage Scissors	24-25
Bandage Shears	25
Barnhill Adenoid Curettes	120
Baron Suction Tubes	66, 111, 197
Barraquer Eye Specula	90
Barrett Uterine Tenaculum Forceps	231
Baumgartner Needle Holders	59, 62, 140, 162
Bayonet Dressing Forceps	36, 40, 108, 113, 191-192
Beckman Nasal Speculum	117
Beckman-Adson Retractor	76, 196
Beckman-Weitlaner Retractors	76, 105, 196
Bennett Tibia Retractors	74, 206
Beyer Rongeur	112, 172, 200, 210
Billeau Flexible Ear Loops	106
Biopsy Curettes, Endocervical	234
Biopsy Curettes, Endometrial	234
Biopsy Curettes, Uterine	234
Biopsy Forceps, Uterine	235, 236
Biopsy Punch Forceps, Cervical	235
Biopsy Punch Forceps, Uterine	235, 236

Bishop-Harmon Dressing Forceps	38, 94, 135, 159
Bishop-Harmon Tissue Forceps	38, 88, 94, 135, 159
Bistouri, Converse	146
Blades, Adenotome	120
Blades, Mouth Gags	169-170
Bone Clamps	214
Bone Cutting Forceps	186-187, 211-212
Bone Holding Forceps	213
Bone Rasps	217
Bone Splitting Forceps	151, 211
Bonn Micro Suturing Forceps	38, 87, 94
Bonney Tissue Forceps	34
Bowman Lacrimal Probes	92
Boys-Allis Tissue Forceps	42
Bozeman Sponge Forceps	56, 224
Braun Episiotomy Scissors	20, 237
Braun Atraumatic Tenaculum Forceps	231
Brown Needle Holder	58, 139
Brown-Adson Tissue Forceps	33, 133, 160, 190
Brun Mastoid Curettes	107, 198, 208
Bruns Plaster Shears	25, 214
Buck Ear Curettes	106
Buie Pile Clamp	220
Bulldog Clamps	178
Burke Cervical Biopsy Punch Forceps	235
Button End Knife	146

**C**

Capsulotomy Scissors	17, 100
Cartilage Clamp	218
Cartilage Crusher	149
Cartilage Scissors	20, 218
Cast Shears	25, 214
Cast Spreader	214
Cast Spreaders	214
Castroviejo Corneal Scissors	17, 100
Castroviejo Needle Holders	60, 64, 103, 140, 162
Castroviejo Suturing Forceps	39, 95
Cervical Biopsy Punch Forceps	235
Cervical Rongeur	202
Chisels	118, 146
Clamps, Bone	214
Clamps, Bulldog	178
Clamps, Cartilage	218
Clamps, Coarctation	183
Clamps, Peripheral Vascular	178
Clamps, Pile	220
Clamps, Towel	55
Clamps, Vascular	178-180

Clamps, Vascular, Pediatric	179
Clamps, Vessel	223
Coarctation Clamp	183
Collier Needle Holder	57, 138
Columella Forceps	135
Comedone Extractors	167
Converse Bistouri	146
Converse Needle Holder	58, 139
Converse Skin Hooks	82, 143, 166
Cooley Pediatric Vascular Clamps	179
Cooley Vascular Clamp	179-180
Cooley-Satinsky Clamp	180
Corneal Scissors	17, 100
Cottle Cartilage Crusher	149
Cottle Columella Forceps	135
Cottle Dorsal Scissors	19, 131
Cottle Knife Guide & Retractor	78, 144
Cottle Nasal Knife	146
Cottle Osteotomes	148
Cottle Rasp	150
Cottle Retractors	77-78, 142, 144
Cottle Septum Elevator	147
Cottle Skin Elevator	147
Cottle Skin Hooks	83, 142-143, 166
Crafoord Coarctation Clamp	183
Crile Forceps	47, 137
Crile Retractor	68
Crile-Wood Needle Holders	59, 63, 140, 184
Crusher, Cartilage	149
Currettes, Adenoid	120
Currettes, Bone	107, 198-199, 208, 209
Currettes, Ear	106-107
Currettes, Endocervical	234
Currettes, Endometrial	234
Currettes, Mastoid	107, 198, 208
Currettes, Neurosurgical	198-199
Currettes, Spinal Fusion	199, 209
Currettes, Suction	234
Currettes, Uterine	233-234
Currettes, Uterine Biopsy	234
Cushing Dressing Forceps	36, 113, 191
Cushing Pituitary Rongeur	203
Cushing Vein Retractor	69, 195

**D**

Dandy Forceps	48, 193
Davis Mouth Gag	169
Day Ear Hooks	82, 106
Dean Dissecting Scissors	15, 173

Deaver Retractors	71
DeBakey Bulldog Clamps	178
DeBakey Peripheral Vascular Clamps	178
Delicate Scissors	11, 122
Derf Needle Holders	57, 61, 138, 161
Desmarres Lid Retractors	77, 90, 141, 165
Diamond Pin & Wire Cutter	214
Dilators, Tracheal	177
Dilators, Uterine	228, 229, 230
Dissecting Hooks	85, 204
Dissectors	78, 175, 205
Dissectors, Neurosurgical	205
Dissectors, Tonsil	78, 175
Dorsal Scissors	18-19, 130-131
Doyen Intestinal Forceps	51
Doyen Vaginal Retractors	80, 225
Doyen-Jansen Mouth Gags	171
Dressing Forceps	32-40, 44, 86, 92, 94, 108, 113, 132-135, 158-160, 181, 189-192
Dressing Forceps, Bayonet	36, 40, 108, 113, 191-192
Duckbill Rongeur	187, 201, 212
Dura Hooks	85, 142, 165, 204

**E**

Ear Currettes	106-107
Ear Forceps	36, 40, 88, 108-110
Ear Hooks	82, 106
Ear Loops	106
Elevators	90, 118, 147, 171, 204-205, 216-217
Elevators, Bone	216, 217
Elevators, Neurosurgical	204-205
Elevators, Periosteal	90, 118, 147, 171, 204-205, 216, 217
Elevators, Septum	118, 147
Elevators, Skin	147
Endocervical Biopsy Curette	234
Endometrial Biopsy Suction Curette	234
Episiotomy Scissors	20, 237
Eppendorfer Uterine Biopsy Forceps	236
Esmarch Bandage Shears	25
Ethmoid Forceps	119
Evans Dressing Forceps	44
Eve Tonsil Snare	175
Eye Dressing Forceps	37-38, 92, 134, 158
Eye Specula	90
Eye Suture Scissors	23, 102, 129, 157

**F**

Farabeuf Retractor	69
--------------------	----

- Feilchenfeld Splinter Forceps ..... 45  
 Ferris-Smith Tissue Forceps ..... 40, 191  
 Finger Ring Cutter ..... 167  
 Finger Ring Cutter, Blade ..... 167  
 Fisher Tonsil Knife-Dissector ..... 175  
 Fixation Hook ..... 81, 91, 143, 165  
 Flexible Ear Loops ..... 106  
 Foerster Sponge Forceps ..... 56, 173, 224  
 Fomon Dorsal Scissors ..... 18, 130  
 Fomon Lower Lateral Scissors ..... 18, 130  
 Fomon Rasps ..... 150  
 Fomon Retractor ..... 77, 142  
 Fomon Saber-Back Scissors ..... 18, 130  
 Fomon Scissors ..... 18, 130  
 Forceps, Adair ..... 42, 43, 226  
 Forceps, Adson ..... 33-34, 36, 44, 50, 86, 108,  
     113, 133, 160, 190, 192, 194  
 Forceps, Alligator ..... 88, 109, 110, 114  
 Forceps, Allis ..... 41-42, 185, 218  
 Forceps, Ash ..... 151  
 Forceps, Babcock ..... 43, 186  
 Forceps, Baby Babcock ..... 43  
 Forceps, Baby Crile ..... 47, 137  
 Forceps, Ballentine ..... 54, 237  
 Forceps, Barrett ..... 231  
 Forceps, Bayonet ..... 36, 40, 108, 113, 191-192  
 Forceps, Biopsy Punch ..... 235, 236  
 Forceps, Bishop-Harmon ..... 38, 88, 94, 135, 159  
 Forceps, Bone Cutting ..... 186-187, 211-212  
 Forceps, Bone Holding ..... 213  
 Forceps, Bone Splitting ..... 151, 211  
 Forceps, Bonn ..... 38, 87, 94  
 Forceps, Bonney ..... 34  
 Forceps, Boys-Allis ..... 42  
 Forceps, Bozeman ..... 56, 224  
 Forceps, Braun Atraumatic Tenaculum ..... 231  
 Forceps, Brown-Adson ..... 33, 133, 160, 190  
 Forceps, Burke ..... 235  
 Forceps, Castroviejo ..... 39, 95  
 Forceps, Cervical Biopsy ..... 235  
 Forceps, Columella ..... 135  
 Forceps, Cottle ..... 135  
 Forceps, Crile ..... 47, 137  
 Forceps, Cushing ..... 36, 113, 191  
 Forceps, Dandy ..... 48, 193  
 Forceps, Doyen ..... 51  
 Forceps, Dressing ..... 32-40, 44, 86, 92, 94, 108, 113,  
     132-135, 158-160, 181, 189-192  
 Forceps, Ear ..... 36, 40, 88, 108-110  
 Forceps, Ear Dressing ..... 36, 40, 108, 113, 158  
 Forceps, Eppendorfer ..... 236  
 Forceps, Ethmoid ..... 119  
 Forceps, Evans ..... 44  
 Forceps, Eye Dressing ..... 37-38, 92, 134, 158  
 Forceps, Feilchenfeld ..... 45  
 Forceps, Ferris-Smith ..... 40, 191  
 Forceps, Foerster ..... 56, 173, 224  
 Forceps, Gall Duct ..... 51  
 Forceps, Gemini-Mixer ..... 51, 182  
 Forceps, Gerald ..... 35, 191  
 Forceps, Gruenewald ..... 36, 108, 113, 192  
 Forceps, Halsted ..... 46, 89, 136, 193  
 Forceps, Hartman ..... 46, 89, 109, 136  
 Forceps, Heaney ..... 54, 237  
 Forceps, Heaney-Ballentine ..... 54, 237  
 Forceps, Hemostatic ..... 46-54, 86, 89, 136-137,  
     174, 182, 193-194, 218, 224  
 Forceps, Hulka ..... 231  
 Forceps, Hysterectomy ..... 54, 237  
 Forceps, Intestinal ..... 43, 51, 186  
 Forceps, Intra-Thoracic Artery ..... 52, 183  
 Forceps, Iris ..... 37, 39, 87, 93, 95, 135, 159  
 Forceps, Jacobs ..... 54, 232  
 Forceps, Jacobson ..... 52, 86, 89, 137, 194  
 Forceps, Jansen ..... 36, 108, 113, 192  
 Forceps, Jansen-Middleton ..... 119  
 Forceps, Kantrowitz ..... 53, 183  
 Forceps, Kelly ..... 46  
 Forceps, Kern ..... 213  
 Forceps, Kevorkian-Young ..... 236  
 Forceps, Kidney Stone ..... 223  
 Forceps, Kleinert-Kutz ..... 151, 211  
 Forceps, Knight ..... 115  
 Forceps, Kocher ..... 47  
 Forceps, Lahey ..... 51, 176, 226  
 Forceps, Lewin ..... 213  
 Forceps, Lucae ..... 40, 108  
 Forceps, McPherson ..... 39, 87, 95  
 Forceps, Micro ..... 34, 38, 39, 86-89, 95, 110,  
     133, 135, 137, 159-160, 190, 194  
 Forceps, Micro Adson ..... 34, 86, 133, 160, 190  
 Forceps, Micro Alligator ..... 88, 110  
 Forceps, Micro Iris ..... 39, 87, 95, 135, 159  
 Forceps, Micro Suture ..... 38-39, 87, 94-95  
 Forceps, Micro Tissue ..... 34, 38-39, 86-88,  
     94-95, 133, 135, 159-160, 190  
 Forceps, Mixer ..... 50, 51, 53, 182-183  
 Forceps, Mosquito ..... 46, 52, 86, 89, 136-137, 193-194  
 Forceps, Nasal Polyp ..... 114-115  
 Forceps, Noyes ..... 109, 114  
 Forceps, Pelvic Dressing ..... 44  
 Forceps, Pennington ..... 50, 218

Forceps, Potts-Smith . . . . .35, 44, 181  
 Forceps, Pratt . . . . .48, 218  
 Forceps, Randall . . . . .223  
 Forceps, Rankin . . . . .48  
 Forceps, Rochester-Carmalt . . . . .49  
 Forceps, Rochester-Ochsner . . . . .49  
 Forceps, Rochester-Pean . . . . .48, 224  
 Forceps, Ruskin . . . . .151, 186, 211  
 Forceps, Russian . . . . .34, 181, 238  
 Forceps, Sarot . . . . .52, 183  
 Forceps, Sawtell . . . . .52, 174  
 Forceps, Schmidt . . . . .52, 174  
 Forceps, Schroeder . . . . .231-232  
 Forceps, Schroeder Atraumatic Tenaculum . . . . .231  
 Forceps, Schroeder-Braun . . . . .231  
 Forceps, Semken . . . . .32, 132, 158, 189  
 Forceps, Septum . . . . .115, 119, 151  
 Forceps, Singley . . . . .35, 176  
 Forceps, Splinter . . . . .45  
 Forceps, Sponge . . . . .56, 173, 224  
 Forceps, Sterilizer . . . . .239  
 Forceps, Stille-Horsley . . . . .187, 212  
 Forceps, Stille-Liston . . . . .187, 211  
 Forceps, Suturing . . . . .33, 38-39, 87, 94-95, 133, 160  
 Forceps, Takahashi . . . . .119  
 Forceps, Teale . . . . .232  
 Forceps, Tenaculum . . . . .43, 226, 231  
 Forceps, Tendon Pulling . . . . .218  
 Forceps, Thoms-Allis . . . . .42, 220, 226  
 Forceps, Thoms-Gaylor . . . . .235  
 Forceps, Thoracic . . . . .52-53, 183  
 Forceps, Tischler . . . . .235  
 Forceps, Tissue . . . . .32-44, 86-88, 93-95,  
 . . . . .132-133, 135, 158-160, 176,  
 . . . . .181, 185, 189-191, 218, 220, 226, 238  
 Forceps, Tonsil . . . . .52, 174  
 Forceps, Tonsil Hemostat . . . . .52, 174  
 Forceps, Traction . . . . .176, 226  
 Forceps, Tungsten Carbide . . . . .44  
 Forceps, Uterine . . . . .54, 231, 232  
 Forceps, Uterine Biopsy . . . . .235, 236  
 Forceps, Vulsellum . . . . .54, 232  
 Forceps, Watson-Williams . . . . .115  
 Forceps, White . . . . .174  
 Forceps, Wilde . . . . .40, 108  
 Frazier Dura Hook . . . . .85, 142, 165, 204  
 Frazier Suction Tubes . . . . .66, 111, 116, 197  
 Freer Elevator . . . . .90, 118, 147, 171, 205

Freer Septum Knife . . . . .118

**G**

Gall Duct Forceps . . . . .51  
 Gelpi Retractors . . . . .75, 196, 226  
 Gemini-Mixter Hemostatic Forceps . . . . .51, 182  
 Gerald Dressing Forceps . . . . .35, 191  
 Gerald Tissue Forceps . . . . .35, 191  
 Gigli Saw Handles . . . . .202, 213  
 Gigli-Strully Saw Handle . . . . .202, 213  
 Gillies Converse Skin Hooks . . . . .82, 143, 166  
 Glabella Rasp . . . . .150  
 Goelet Retractor . . . . .69  
 Gradle Scissors . . . . .23, 102, 129, 157  
 Graves Vaginal Specula . . . . .227  
 Green Retractor . . . . .68  
 Gruenwald Adson Forceps . . . . .36, 108, 113, 192  
 Gruenwald Forceps . . . . .36, 108, 113, 192  
 Gruenwald Jansen Forceps . . . . .36, 108, 113, 192  
 Guthrie Fixation Hook . . . . .81, 91, 143, 165  
 Guyon Vessel Clamp . . . . .223

**H**

Hajek-Ballenger Septum Elevator . . . . .118  
 Halsey Needle Holders . . . . .57, 61, 138, 161  
 Halsted Mosquito Forceps . . . . .46, 89, 136, 193  
 Hank Uterine Dilators . . . . .228  
 Harrington Retractors . . . . .73, 185  
 Hartman Ear Dressing Forceps . . . . .109  
 Hartman Mosquito forceps . . . . .46, 89, 136  
 Heaney Hysterectomy Forceps . . . . .54, 237  
 Heaney Needle Holders . . . . .60, 63  
 Heaney Retractors . . . . .71, 225  
 Heaney Uterine Biopsy Curette . . . . .234  
 Heaney-Ballentine Hysterectomy Forceps . . . . .54, 237  
 Heaney-Simon Retractor . . . . .71, 225  
 Hegar Uterine Dilators . . . . .229  
 Heiss Retractors . . . . .67, 104, 163  
 Hemostatic Forceps . . . . .46-54, 86, 89, 136-137,  
 . . . . .174, 182, 193-194, 218, 224  
 Hibbs Osteotomes . . . . .215  
 Hibbs Retractor . . . . .74, 185  
 Hibbs Spratt Spinal Fusion Curettes . . . . .199, 209  
 Hill-Ferguson Retractors . . . . .80, 221  
 Hirschman Anoscopes . . . . .221  
 Hohmann Retractors . . . . .79, 207  
 Hooks, Adson . . . . .85, 204  
 Hooks, Cottle . . . . .83, 142-143, 166

Hooks, Day	82, 106
Hooks, Dissecting	85, 204
Hooks, Dura	85, 142, 165, 204
Hooks, Ear	82, 106
Hooks, Fixation	81, 91, 143, 165
Hooks, Frazier	85, 142, 165, 204
Hooks, Gillies Converse	82, 143, 166
Hooks, Guthrie	81, 91, 143, 165
Hooks, Jackson	84, 177
Hooks, Jameson	81, 91
Hooks, Joseph	84, 144, 167
Hooks, News	84, 177
Hooks, Skin	81-85, 91, 142-144, 165-167
Hooks, Stevens	81, 91
Hooks, Strabismus	81-82, 91-92
Hooks, Tenaculum	84, 177
Hooks, Tenotomy	81, 91
Hooks, Tyrell	81, 91, 143, 165
Hooks, Von Graefe	82, 92
Hulka Tenaculum Forceps	231
Hurd Dissector-Retractor	78, 175
Hysterectomy Forceps	54, 237

**I**

Infant Nasal Speculum	117
Intestinal Forceps	43, 51, 186
Intra-Thoracic Artery Forceps	52, 183
Iris Scissors	12-14, 28, 31, 96-98, 123-125, 152-154
Iris Suturing Forceps	39, 87, 95
Iris Tissue Forceps	37, 39, 87, 93, 95, 135, 159

**J**

Jackson Tenaculum Hook Retractor	84, 177
Jackson Tracheal Retractor	78, 142, 176
Jacobs Uterine Vulsellum Forceps	54, 232
Jacobson Micro Mosquito Forceps	52, 86, 89, 137, 194
Jameson Strabismus Hook	81, 91
Jansen Bayonet Dressing Forceps	36, 108, 113, 192
Jansen Retractor	67, 104, 195
Jansen-Middleton Septum Forceps	119
Jennings Mouth Gags	171
Joseph Button End Knife	146
Joseph Hooks	84, 144, 167
Joseph Nasal Saw	149
Joseph Periosteal Elevator	147
Joseph Rasp	149

Joseph Scissors	19, 131, 156
-----------------	--------------

**K**

Kantrowitz Thoracic Forceps	53, 183
Kelly Forceps	46
Kelly Retractors	73
Kern Bone Holding Forceps	213
Kevoorkian-Young Endocervical Biopsy Curettes	234
Kevoorkian-Young Uterine Biopsy Punch Forceps	236
Key Periosteal Elevators	217
Kidney Stone Forceps	223
Killian Septum Specula	117
Kleinert-Kutz Bone Splitting Forceps	151, 211
Knapp Iris Scissors	14, 98, 125, 154
Knapp Retractors	77, 90, 141, 164
Knapp Rollet Retractors	77, 90, 141, 164
Knife Guide	78, 144
Knife Handle	6
Knife, Button End	146
Knife, Nasal	118, 146
Knife, Septum	118
Knife, Tonsil	175
Knight Nasal Scissors	18, 114
Knight Septum Forceps	115
Knowles Bandage Scissors	24
Kocher Forceps	47
Krause Nasal Snare	117

**L**

LaBorde Tracheal Dilator	177
Lacrimal Probes	92
LaForce Adenotome Blades	120
LaForce Adenotomes	120
Lahey Gall Duct Forceps	51
Lahey Scissors	10, 121, 152
Lahey Traction Forceps	176, 226
Lambotte Osteotomes	215
Langenbeck Periosteal Elevators	204, 216
Langenbeck Retractors	70
Leksell Rongeurs	188, 201, 212
Lempert Rongeurs	112
Lewin Bone Holding Forceps	213
Liberator Splinter	45
Lid Retractors	77, 90, 141, 165
Lister Bandage Scissors	24
Littauer Stitch Scissors	21, 127, 157
Littler Suture Carrying Scissors	23, 129
Lorna Towel Clamps	55

Love-Gruenwald Pituitary Rongeurs . . . . . 203  
 Lower Lateral Scissors . . . . . 18, 130  
 Lowman Bone Clamps . . . . . 214  
 Lucae Ear Forceps . . . . . 40, 108

**M**

Mallet . . . . . 107, 145, 216  
 Maltz Nasal Rasp . . . . . 150  
 Martin Cartilage Clamp . . . . . 218  
 Martin Cartilage Scissors . . . . . 20, 218  
 Masson Needle Holder . . . . . 60, 64  
 Mastoid Currettes . . . . . 107, 198, 208  
 Mayo (Guyon) Vessel Clamp . . . . . 223  
 Mayo Dissecting Scissors . . . . . 9, 26, 29-30, 181, 238  
 Mayo-Hegar (Masson) Needle Holders . . . . . 60, 64  
 Mayo-Hegar Needle Holders . . . . . 59-60, 63-64, 140  
 McIvor Mouth Gag . . . . . 169  
 McPherson Micro Iris Suturing Forceps . . . . . 39, 87, 95  
 Mead Mallet . . . . . 107, 145, 216  
 Metzbaum Scissors . . . . . 10, 27, 30, 121, 152, 173  
 Meyerding Retractors . . . . . 72, 195  
 Micro Adson Dressing Forceps . . . . . 34, 86, 133, 160, 190  
 Micro Adson Tissue Forceps . . . . . 34, 86, 133, 160, 190  
 Micro Dressing Forceps . . . . . 34, 38-39, 86, 94, 133, 135, 159-160, 190  
 Micro Ear Forceps . . . . . 88, 110  
 Micro Mosquito Forceps . . . . . 52, 86, 89, 137, 194  
 Micro Suturing Forceps . . . . . 38-39, 87, 94-95  
 Micro Tissue Forceps . . . . . 34, 38-39, 86-88, 94-95, 133, 135, 159-160, 190  
 Mini-Hohmann Retractors . . . . . 79, 207  
 Mixer Forceps . . . . . 50, 51, 53, 182-183  
 Mixer Thoracic Forceps . . . . . 53, 183  
 Morselizer, Septal . . . . . 119  
 Mosquito Forceps . . . . . 46, 52, 86, 89, 136-137, 193-194  
 Mouth Gags . . . . . 169, 171

**N**

Nail Nippers . . . . . 219  
 Nail Splitters . . . . . 219  
 Nasal Forceps . . . . . 113-115  
 Nasal Knives . . . . . 146  
 Nasal Polyp Forceps . . . . . 115  
 Nasal Rasp . . . . . 149-150  
 Nasal Saw . . . . . 149  
 Nasal Scissors . . . . . 18, 114  
 Nasal Septum Forceps . . . . . 115  
 Nasal Snare . . . . . 117

Nasal Specula . . . . . 116-117  
 Needle Holders . . . . . 57-64, 103, 138-140, 161-162, 184  
 Needle Holders, Micro Surgery . . . . . 60, 64, 103, 140, 162  
 Needle Holders, Tungsten Carbide . . . . . 61-64, 103, 140, 184  
 Neivert Osteotome . . . . . 148  
 Neivert-Anderson Osteotome . . . . . 148  
 New's Suture Scissors . . . . . 21  
 News Hook . . . . . 84, 177  
 Novak Endometrial Suction Biopsy Curettes . . . . . 234  
 Noyes Alligator Forceps . . . . . 109, 114

**O**

O'Brien Stitch Scissors . . . . . 21, 101, 127, 157  
 Olsen-Hegar Needle Holders . . . . . 58, 62, 139, 162  
 Operating Scissors . . . . . 7-8  
 Osteotomes . . . . . 148, 215

**P**

Parker Retractors . . . . . 69  
 Pederson Vaginal Specula . . . . . 227  
 Pediatric Vascular Clamps . . . . . 179  
 Peers Towel Clamps . . . . . 55  
 Pelvic Dressing Forceps . . . . . 44  
 Penfield Dissectors . . . . . 205  
 Pennington Hemostatic Forceps . . . . . 50, 218  
 Pennington Septum Elevator . . . . . 118  
 Pennington Tissue Forceps . . . . . 50, 218  
 Periosteal Elevators . . . . . 90, 118, 147, 171, 204-205, 216, 217  
 Peripheral Vascular Clamp . . . . . 178  
 Pile Clamp . . . . . 220  
 Pillar Retractor . . . . . 78, 175  
 Pin & Wire Cutter . . . . . 214  
 Pin Wheel . . . . . 239  
 Pituitary Rongeurs . . . . . 203  
 Plaster Shears . . . . . 25, 214  
 Plastic Surgery Scissors . . . . . 11, 121-131  
 Poole Suction Tubes . . . . . 65  
 Potts-Smith Dressing Forceps . . . . . 35, 44, 181  
 Potts-Smith Tissue Forceps . . . . . 35, 181  
 Pratt Rectal Speculum . . . . . 220  
 Pratt T-Shape Forceps . . . . . 48, 218  
 Pratt Uterine Dilators . . . . . 230  
 Probes, Lacrimal . . . . . 92  
 Putti Bone Rasps . . . . . 217

**R**

Ragnell Retractor . . . . . 68, 141, 176

Randall Endometrial Suction Biopsy Curette	234
Randall Kidney Stone Forceps	223
Rankin Forceps	48
Rasps	149, 150, 217
Rectal Speculum	220
Retractor Blades	170
Retractor-Speculum, Aufricht	78, 148
Retractors, Adson	76, 196
Retractors, Alm	67, 104, 163
Retractors, Aufricht	78, 148
Retractors, Beckman-Adson	76, 196
Retractors, Beckman-Weitlaner	76, 105, 196
Retractors, Bennett	74, 206
Retractors, Cottle	77-78, 142, 144
Retractors, Crile	68
Retractors, Cushing	69, 195
Retractors, Deaver	71
Retractors, Desmarest	77, 90, 141, 165
Retractors, Doyen	80, 225
Retractors, Farabeuf	69
Retractors, Fomon	77, 142
Retractors, Gelpi	75, 196, 226
Retractors, Goelet	69
Retractors, Green	68
Retractors, Harrington	73, 185
Retractors, Heaney	71, 225
Retractors, Heaney-Simon	71, 225
Retractors, Heiss	67, 104, 163
Retractors, Hibbs	74, 185
Retractors, Hill-Ferguson	80, 221
Retractors, Hohmann	79, 207
Retractors, Hurd	78, 175
Retractors, Jackson	78, 142, 176-177
Retractors, Jansen	67, 104, 195
Retractors, Kelly	73
Retractors, Knapp	77, 90, 141, 164
Retractors, Langenbeck	70
Retractors, Lid	77, 90, 141, 165
Retractors, Meyerding	72, 195
Retractors, Mini-Hohmann	79, 207
Retractors, Parker	69
Retractors, Pillar	78, 175
Retractors, Ragnell	68, 141, 176
Retractors, Rectal	221
Retractors, Ribbon	70
Retractors, Richardson	72, 184
Retractors, Richardson-Eastman	73
Retractors, Ring Tongue	170
Retractors, Senn	68, 141, 164, 176-177
Retractors, Tibia	74, 206

Retractors, Tongue	170
Retractors, Tracheal	78, 142, 176-177
Retractors, U.S. Army	70
Retractors, Vaginal	80, 225
Retractors, Volkman	74, 206
Retractors, Weitlaner	75-76, 105, 196
Ribbon Retractors	70
Richardson Retractors	72, 184
Richardson-Eastman Retractors	73
Ring Cutter	167
Ring Cutter, Blade	167
Ring Tongue Retractor	170
Rochester-Carmalt Forceps	49
Rochester-Ochsner Forceps	49
Rochester-Pean Forceps	48, 224
Rongeurs, Beyer	112, 172, 200, 210
Rongeurs, Cervical	202
Rongeurs, Cushing	203
Rongeurs, Duckbill	187, 201, 212
Rongeurs, Leksell	188, 201, 212
Rongeurs, Lempert	112
Rongeurs, Love-Gruenwald	203
Rongeurs, Pituitary	203
Rongeurs, Ruskin	112, 186, 200, 210
Rongeurs, Sauerbruch	188
Rongeurs, Spurling	203
Rongeurs, Spurling-Kerrison	202
Rongeurs, Still-Luer	187, 201, 212
Rubin Septal Morselizer	119
Ruskin Bone Cutting Forceps	186, 211
Ruskin Bone Splitting Forceps	151, 211
Ruskin Rongeurs	112, 186, 200, 210
Russian Tissue Forceps	34, 181, 238
Ryder Needle Holders	64, 184

## S

Sarot Intra-Thoracic Artery Forceps	52, 183
Sauerbruch Rongeur	188
Saw	149
Saw Handles	202, 213
Sawtell Tonsil Hemostatic Forceps	52, 174
Scalpel Handles	6
Schamberg Comedones	167
Schmidt Tonsil Hemostatic Forceps	52, 174
Schroeder Uterine Tenaculum Forceps	231
Schroeder Uterine Vulsellum Forceps	232
Schroeder-Braun Uterine Tenaculum Forceps	231
Scissors, Bandage	24-25
Scissors, Braun	20, 237



- Scissors, Capsulotomy . . . . .17, 100  
 Scissors, Cartilage . . . . .20, 218  
 Scissors, Castroviejo . . . . .17, 100  
 Scissors, Corneal . . . . .17, 100  
 Scissors, Cottle . . . . .19, 131  
 Scissors, Dean . . . . .15, 173  
 Scissors, Delicate . . . . .11, 122  
 Scissors, Dorsal . . . . .18-19, 130-131  
 Scissors, Episiotomy . . . . .20, 237  
 Scissors, Esmarch . . . . .25  
 Scissors, Eye Suture . . . . .23, 102, 129, 157  
 Scissors, Fomon . . . . .18, 130  
 Scissors, Gradle . . . . .23, 102, 129, 157  
 Scissors, Heavy Duty . . . . .25  
 Scissors, Iris . . . . .12-14, 28, 31, 96-98, 123-125, 152-154  
 Scissors, Joseph . . . . .19, 131, 156  
 Scissors, Knapp . . . . .14, 98, 125, 154  
 Scissors, Knight . . . . .18, 114  
 Scissors, Knowles . . . . .24  
 Scissors, Lahey . . . . .10, 121, 152  
 Scissors, Lister . . . . .24  
 Scissors, Littauer . . . . .21, 127, 157  
 Scissors, Littler . . . . .23, 129  
 Scissors, Lower Lateral . . . . .18, 130  
 Scissors, Martin . . . . .20, 218  
 Scissors, Mayo . . . . .9, 26, 29-30, 181, 238  
 Scissors, Metzenbaum . . . . .10, 27, 30, 121, 152, 173  
 Scissors, Nasal . . . . .18, 114  
 Scissors, New's . . . . .21  
 Scissors, O'Brien . . . . .21, 101, 127, 157  
 Scissors, Operating . . . . .7-8  
 Scissors, Plastic Surgery . . . . .11, 121-131  
 Scissors, Saber-Back . . . . .18, 130  
 Scissors, Spencer . . . . .21, 101, 127, 157  
 Scissors, Stevens . . . . .16, 99, 126, 155  
 Scissors, Stitch . . . . .21-23, 101-102, 127-129, 156-157, 162  
 Scissors, Strabismus . . . . .15, 28, 96, 125, 152  
 Scissors, SuperCut . . . . .29-31  
 Scissors, Suture . . . . .21-23, 58, 62, 101-102, 127-129, 156-157, 162  
 Scissors, Tenotomy . . . . .16, 99, 126, 155  
 Scissors, Tungsten Carbide . . . . .26-28, 128  
 Scissors, Vannas . . . . .17, 100  
 Scissors, Westcott . . . . .22, 102, 128, 156  
 Scissors, Wire Cutting . . . . .22, 28, 128  
 Semken Dressing Forceps . . . . .32, 132, 158, 189  
 Semken Tissue Forceps . . . . .32, 132, 158, 189  
 Senn Retractors . . . . .68, 141, 164, 176  
 Septal Morselizer . . . . .119  
 Septal Morselizer, Slip-On Guard . . . . .119  
 Septum Elevators . . . . .118, 147  
 Septum Forceps . . . . .115, 119, 151  
 Septum Knife . . . . .118  
 Septum Speculum . . . . .117  
 Shears . . . . .25, 214  
 Sheehan Osteotomes . . . . .148  
 Sims Rectal Speculum . . . . .220  
 Sims Uterine Curettes . . . . .233  
 Singley Tissue Forceps . . . . .35, 176  
 Skin Elevator . . . . .147  
 Skin Hooks . . . . .81-85, 91, 143-144, 165-167  
 Snare, Nasal . . . . .117  
 Snare, Tonsil . . . . .175  
 Sounds . . . . .222  
 Specula, Barraquer . . . . .90  
 Specula, Eye . . . . .90  
 Specula, Graves . . . . .227  
 Specula, Infant . . . . .117  
 Specula, Killian . . . . .117  
 Specula, Nasal . . . . .116, 117  
 Specula, Pederson . . . . .227  
 Specula, Pratt . . . . .220  
 Specula, Rectal . . . . .220  
 Specula, Septum . . . . .117  
 Specula, Sims . . . . .220  
 Specula, Tieck-Halle . . . . .117  
 Specula, Vaginal . . . . .227  
 Specula, Vienna . . . . .116  
 Spencer Stitch Scissors . . . . .21, 101, 127, 157  
 Spinal Fusion Curette . . . . .199, 209  
 Splinter Forceps . . . . .45  
 Splinter Liberator . . . . .45  
 Sponge Forceps . . . . .56, 173, 224  
 Spratt Curettes . . . . .107, 198-199, 208-209  
 Spurling Pituitary Rongeurs . . . . .203  
 Spurling-Kerrison Cervical Rongeurs . . . . .202  
 Sterilizer Forceps . . . . .239  
 Stevens Tenotomy Hook . . . . .81, 91  
 Stevens Tenotomy Scissors . . . . .16, 99, 126, 155  
 Stille-Luer Rongeurs . . . . .187, 201, 212  
 Stille Vessel Clamp . . . . .223  
 Stille-Horsley Bone Cutting Forceps . . . . .187, 212  
 Stille-Liston Bone Cutting Forceps . . . . .187, 211  
 Stitch Scissors . . . . .21-23, 101-102, 127-129, 156-157  
 Strabismus Hooks . . . . .81-82, 91-92  
 Strabismus Scissors . . . . .15, 28, 96, 125, 152  
 Suction Tube Tip . . . . .65, 168  
 Suction Tubes . . . . .65-66, 111, 116, 168, 197  
 SuperCut Scissors . . . . .29-31

Suture Scissors . . . 21-23, 58, 62, 101-102, 127-129, 156-157  
 Suturing Forceps . . . . . 33, 38-39, 87, 94-95, 133, 160

**T**

Takahashi Ethmoid Forceps . . . . . 119  
 Teale Vulsellum Forceps . . . . . 232  
 Tenaculum Forceps, Atraumatic . . . . . 231  
 Tenaculum Forceps, Tissue . . . . . 43, 226, 231  
 Tenaculum Forceps, Uterine . . . . . 231  
 Tenaculum Hook . . . . . 84, 177  
 Tendon Pulling Forceps . . . . . 218  
 Tenotomy Hook . . . . . 81, 91  
 Tenotomy Scissors . . . . . 16, 99, 126, 155  
 Thoms-Allis Forceps . . . . . 42, 220, 226  
 Thoms-Gaylor Uterine Biopsy Forceps . . . . . 235  
 Thoracic Forceps . . . . . 52, 53, 183  
 Three Prong Cast Spreader . . . . . 214  
 Three Prong Sterilizer Forceps . . . . . 239  
 Thyroid Retractor . . . . . 68, 76  
 Tibia Retractor . . . . . 74, 206  
 Tieck-Halle Infant Speculum . . . . . 117  
 Tischler Cervical Biopsy Punch Forceps . . . . . 235  
 Tissue Forceps . . . . . 32-44, 86-88, 93-95,  
 . . . . . 132-133, 135, 158-160, 176,  
 . . . . . 181, 185, 189-191, 218, 220, 226, 238  
 Tongue Depressor . . . . . 172  
 Tongue Retractor . . . . . 170  
 Tongue Retractor Blades . . . . . 170  
 Tonsil Dissectors . . . . . 78, 175  
 Tonsil Hemostatic Forceps . . . . . 52, 174  
 Tonsil Knife . . . . . 175  
 Tonsil Seizing Forceps . . . . . 52, 174  
 Tonsil Snare . . . . . 175  
 Towel Clamps . . . . . 55  
 Towel Clamps, Non-Perforating . . . . . 55  
 Tracheal Dilator . . . . . 177  
 Tracheal Retractor . . . . . 78, 142, 176  
 Traction Forceps . . . . . 176, 226  
 Trousseau Tracheal Dilator . . . . . 177  
 Tungsten Carbide Forceps . . . . . 44  
 Tungsten Carbide Needle Holders . . . . . 61-64, 103, 140, 184  
 Tungsten Carbide Scissors . . . . . 26-28, 128  
 Tyrell Hooks . . . . . 81, 91, 143, 165

**U**

U.S. Army Retractor . . . . . 70  
 Urethral Sound . . . . . 222  
 Uterine Biopsy Curettes . . . . . 234

Uterine Biopsy Forceps . . . . . 236  
 Uterine Biopsy Punch Forceps . . . . . 235-236  
 Uterine Curettes . . . . . 233, 234  
 Uterine Dilators . . . . . 228, 229, 230  
 Uterine Tenaculum Forceps . . . . . 231-232  
 Uterine Vulsellum Forceps . . . . . 54, 232

**V**

Vaginal Retractors . . . . . 80, 225  
 Vaginal Specula . . . . . 227  
 Van Buren Urethral Sounds . . . . . 222  
 Vannas Capsulotomy Scissors . . . . . 17, 100  
 Vascular Clamps . . . . . 178-180  
 Vein Retractor . . . . . 69, 195  
 Vessel Clamp . . . . . 223  
 Vienna Nasal Specula . . . . . 116  
 Volkman Retractors . . . . . 74, 206  
 Von Graefe Strabismus Hooks . . . . . 82, 92  
 Vulsellum Forceps . . . . . 54, 232

**W**

Wartenberg Pin Wheel . . . . . 239  
 Watson-Williams Nasal Polyp Forceps . . . . . 115  
 Webster Needle Holders . . . . . 57, 61, 138, 161  
 Weitlaner Retractors . . . . . 75-76, 105, 196  
 Westcott Stitch Scissors . . . . . 22, 102, 128, 156  
 White Tonsil Seizing Forceps . . . . . 174  
 Wieder Tongue Depressors . . . . . 172  
 Wilde Ear Forceps . . . . . 40, 108  
 Wire Cutting Scissors . . . . . 22, 28, 128

**Y**

Yankauer Suction Tube . . . . . 65, 168



REF	Page
MH1-226.....	239
MH2-104SS.....	65, 168
MH2-106SS.....	65, 168
MH2-109.....	65, 168
MH2-110SS.....	65, 168
MH2-120.....	169
MH2-122.....	169
MH2-132.....	169
MH2-134.....	169
MH2-136.....	169
MH2-140.....	170
MH2-141.....	170
MH2-142.....	170
MH2-143.....	170
MH2-144.....	170
MH2-145.....	170
MH2-146.....	170
MH2-147.....	170
MH2-149.....	171
MH2-150.....	171
MH2-160.....	171
MH2-161.....	171
MH2-162.....	171
MH2-164.....	171
MH2-180.....	169
MH3-10.....	239
MH4-7.....	6
MH4-8.....	6
MH4-9.....	6
MH4-10.....	6
MH5-2.....	7
MH5-4.....	7
MH5-6.....	7
MH5-8.....	7
MH5-12.....	7
MH5-14.....	7
MH5-16.....	7
MH5-18.....	7
MH5-22.....	7
MH5-24.....	7
MH5-26.....	7
MH5-36.....	8
MH5-46.....	8
MH5-56.....	8
MH5-120.....	9
MH5-120TC.....	26
MH5-122.....	9

REF	Page
MH5-122TC.....	26
MH5-124.....	9
MH5-124TC.....	26
MH5-126.....	9
MH5-126TC.....	26
MH5-128.....	9, 181, 238
MH5-128TC.....	26
MH5-130.....	9, 181, 238
MH5-130TC.....	26
MH5-179.....	10, 121, 152
MH5-179TC.....	27
MH5-180.....	10, 121, 152
MH5-180TC.....	27
MH5-181.....	10, 121
MH5-182.....	10, 121, 173
MH5-182TC.....	27
MH5-182A-TC.....	27
MH5-264.....	15, 173
MH5-270.....	11, 122
MH5-272.....	11, 122
MH5-274.....	11, 122
MH5-276.....	11, 122
MH5-278.....	11, 122
MH5-280.....	11, 122
MH5-288.....	10, 121
MH5-300.....	12, 96, 123, 152
MH5-302.....	12, 96, 123, 152
MH5-304.....	12, 96, 123, 152
MH5-304TC.....	28
MH5-306.....	12, 96, 123, 152
MH5-306TC.....	28
MH5-312.....	15, 96, 125, 152
MH5-312TC.....	28
MH5-314.....	15, 96, 125, 152
MH5-314TC.....	28
MH5-512.....	24
MH5-514.....	24
MH5-516.....	24
MH5-550.....	24
MH5-560.....	24
MH5-568.....	25
MH5-SC-120.....	29
MH5-SC-122.....	29
MH5-SC-124.....	29
MH5-SC-126.....	29
MH5-SC-128.....	30
MH5-SC-130.....	30

REF	Page
MH5-SC-180.....	30
MH5-SC-182.....	30
MH5-SC-304.....	31
MH5-SC-306.....	31
MH6-4.....	32
MH6-6.....	32
MH6-8.....	32
MH6-10.....	32
MH6-14.....	32
MH6-40.....	32
MH6-42.....	32
MH6-44.....	32
MH6-46.....	32
MH6-48.....	32
MH6-104.....	32, 132, 158, 189
MH6-106.....	32, 132, 158, 189
MH6-118.....	33, 133, 160, 190
MH6-118TC.....	44
MH6-120.....	33, 133, 160, 190
MH6-121.....	33, 133, 160, 190
MH6-122.....	33, 133, 160, 190
MH6-123.....	33, 133, 160, 190
MH6-124.....	33, 133, 160, 190
MH6-142.....	34, 181, 238
MH6-144.....	34, 181, 238
MH6-145.....	34, 181, 238
MH6-148.....	34
MH6-150.....	34
MH6-154.....	35, 181
MH6-154TC.....	44
MH6-160.....	35, 181
MH6-165TC.....	44
MH6-180.....	35, 191
MH6-182.....	35, 191
MH6-190.....	36, 113, 191
MH6-192.....	36, 113, 191
MH6-196.....	36, 108, 113, 192
MH6-197.....	36, 108, 113, 192
MH6-198.....	36, 108, 113, 192
MH6-214.....	35, 176
MH6-300.....	45
MH6-304.....	45
MH6-316.....	45
MH6-318.....	45
MH6-400.....	45
MH7-2.....	46, 89, 136
MH7-4.....	46, 89, 136

<b>REF</b>	<b>Page</b>	<b>REF</b>	<b>Page</b>	<b>REF</b>	<b>Page</b>
MH7-8.....	46, 89, 136, 193	MH7-258.....	51, 182	MH8-101TC.....	64, 184
MH7-10.....	46, 89, 136, 193	MH7-259.....	51, 182	MH8-103TC.....	64, 184
MH7-24.....	46, 89, 136	MH7-504.....	55	MH8-104TC.....	64, 184
MH7-26.....	46, 89, 136	MH7-506.....	55	MH8-105TC.....	64, 184
MH7-36.....	46	MH7-514.....	55	MH8-130.....	60
MH7-38.....	46	MH7-515.....	55	MH8-130TC.....	64
MH7-42.....	47	MH7-516.....	55	MH9-98.....	21
MH7-44.....	47	MH7-600.....	56, 224	MH9-100.....	21, 101, 127, 157
MH7-46.....	47	MH7-602.....	56, 224	MH9-102.....	21, 127, 157
MH7-48.....	47	MH7-604.....	56, 173	MH9-104.....	21, 127, 157
MH7-50.....	47, 137	MH7-610.....	56, 173	MH9-110.....	21, 101, 127, 157
MH7-52.....	47, 137	MH7-614.....	56, 173	MH9-124.....	22, 128
MH7-66.....	47	MH7-622.....	56, 224	MH9-124TC.....	28, 128
MH7-68.....	47	MH8-2.....	57, 138	MH10-310.....	65
MH7-92.....	48	MH8-4.....	57, 138, 161	MH10-312.....	65
MH7-94.....	48	MH8-4TC.....	61	MH11-2.....	67, 104, 195
MH7-102.....	48, 193	MH8-6.....	57, 138, 161	MH11-12.....	67, 104, 163
MH7-106.....	48, 218	MH8-6TC.....	61	MH11-13.....	67, 104, 163
MH7-118.....	48	MH8-8.....	57, 138, 161	MH11-14.....	67, 104, 163
MH7-120.....	48	MH8-8A-TC.....	61	MH11-16.....	67, 104, 163
MH7-122.....	48	MH8-14.....	58, 139, 162	MH11-22.....	68
MH7-124.....	48, 224	MH8-15.....	58, 139	MH11-73.....	68, 141, 176
MH7-126.....	48, 224	MH8-15TC.....	62	MH11-74.....	68, 141, 164, 176
MH7-136.....	48	MH8-16.....	58, 139	MH11-76.....	68, 141, 164, 176
MH7-138.....	48	MH8-16TC.....	62	MH11-78.....	68
MH7-140.....	48	MH8-17.....	58	MH11-102.....	69, 195
MH7-142.....	48, 224	MH8-17TC.....	62	MH11-110.....	69
MH7-144.....	48, 224	MH8-20.....	58, 139	MH11-114.....	69
MH7-150.....	49	MH8-22.....	58, 139	MH11-120.....	69
MH7-152.....	49	MH8-40.....	59, 140, 162	MH11-122.....	70
MH7-154.....	49	MH8-40TC.....	62	MH11-132.....	70
MH7-160.....	49	MH8-42.....	59, 140	MH11-136.....	70
MH7-162.....	49	MH8-42TC.....	63	MH11-140.....	70
MH7-164.....	49	MH8-44.....	59, 140	MH11-144.....	70
MH7-170.....	49	MH8-44TC.....	63	MH11-146.....	70
MH7-172.....	49	MH8-46.....	59, 140	MH11-162.....	70
MH7-180.....	49	MH8-46TC.....	63	MH11-164.....	70
MH7-182.....	49	MH8-48.....	59	MH11-166.....	70
MH7-202.....	50	MH8-48TC.....	63	MH11-180.....	71, 225
MH7-206.....	50, 182	MH8-50.....	59, 140	MH11-182.....	71, 225
MH7-220.....	50, 218	MH8-50TC.....	63	MH11-184.....	71, 225
MH7-221.....	50, 218	MH8-51.....	59, 140	MH11-185.....	71, 225
MH7-240.....	50, 194	MH8-51TC.....	63, 184	MH11-202.....	71
MH7-242.....	50, 194	MH8-52.....	59	MH11-206.....	71
MH7-255.....	51, 182	MH8-52TC.....	63, 184	MH11-208.....	71
MH7-256.....	51, 182	MH8-100.....	60	MH11-210.....	71
MH7-257.....	51, 182	MH8-100TC.....	63	MH11-219.....	71

REF	Page	REF	Page	REF	Page
MH11-220.....	71	MH16-10.....	41, 218	MH18-867....	38,88, 94, 135, 159
MH11-240.....	72, 195	MH16-12.....	41, 185	MH18-868.....	38, 94, 135, 159
MH11-242.....	72, 195	MH16-30.....	42, 218	MH18-948.....	39, 87, 95
MH11-244.....	72, 195	MH16-31.....	42	MH18-949.....	39, 87, 95
MH11-246.....	72, 184	MH16-38.....	42, 220, 226	MH18-954.....	39, 95
MH11-247.....	72, 184	MH16-42.....	43	MH18-955.....	39, 95
MH11-248.....	72, 184	MH16-44.....	43	MH18-1396....	12, 96, 123, 153
MH11-249.....	72, 184	MH16-46.....	43, 186	MH18-1398....	12, 96, 123, 153
MH11-260.....	73	MH16-48.....	43, 186	MH18-1415....	13, 97, 124, 153
MH11-260A.....	73	MH16-51.....	43, 226	MH18-1416....	13, 97, 124, 153
MH11-260B.....	73	MH16-160.....	51	MH18-1417....	13, 97, 124, 154
MH11-266.....	73	MH16-162.....	51	MH18-1418....	13, 97, 124, 154
MH11-268.....	73	MH17-2500... 34, 86, 133, 160, 190		MH18-1420....	14, 98, 125, 154
MH11-269.....	73	MH17-2510... 34, 86, 133, 160, 190		MH18-1422....	14, 98, 125, 154
MH11-314.....	73, 185	MH17-2600... 52, 86, 89, 137, 194		MH18-1426....	14, 98, 125, 154
MH11-316.....	73, 185	MH17-2602... 52, 86, 89, 137, 194		MH18-1428....	14, 98, 125, 154
MH11-364.....	75, 196, 226	MH18-35.....	90	MH18-1460....	16, 99, 126, 155
MH11-365.....	75, 196, 226	MH18-36.....	90	MH18-1462....	16, 99, 126, 155
MH11-366.....	75, 196, 226	MH18-100....	77, 90, 141, 164	MH18-1464....	16, 99, 126, 155
MH11-530.....	74, 206	MH18-102....	77, 90, 141, 164	MH18-1466....	16, 99, 126, 155
MH11-532.....	74, 206	MH18-110....	77, 90, 141, 165	MH18-1472....	16, 99, 126, 155
MH11-534.....	74, 206	MH18-112....	77, 90, 165	MH18-1476....	16, 99, 126, 155
MH11-536.....	74, 206	MH18-114....	77, 90, 165	MH18-1486....	22, 102, 128, 156
MH11-538.....	74, 206	MH18-116....	77, 90, 165	MH18-1572....	17, 100
MH11-540.....	74, 206	MH18-412....	81, 91, 143, 165	MH18-1576....	17, 100
MH11-542.....	74, 206	MH18-414....	81, 91, 143, 165	MH18-1578....	17, 100
MH11-544.....	74, 206	MH18-416....	81, 91, 143, 165	MH18-1620....	17, 100
MH11-580.....	74, 206	MH18-440....	81, 91	MH18-1622....	17, 100
MH11-582.....	74, 206	MH18-446....	81, 91	MH18-1652....	23, 102, 129, 157
MH11-590.....	74, 185	MH18-450....	82, 92	MH18-1653....	23, 102, 129, 157
MH11-600.....	75, 105	MH18-452....	82, 92	MH18-1828... 60, 103, 140, 162	
MH11-602.....	75, 105	MH18-454....	82, 92	MH18-1828TC	64, 103, 140
MH11-608.....	75, 105	MH18-740....	92	MH18-1832... 60, 103, 140, 162	
MH11-610.....	75, 105	MH18-742....	92	MH18-1968... 90, 118, 147, 171, 205	
MH11-614.....	75, 105	MH18-744....	92	MH19-220.....	82, 106
MH11-616.....	75, 105	MH18-746....	92	MH19-222....	82, 106
MH11-630-BL.....	76, 105, 196	MH18-748....	92	MH19-280.....	106
MH11-630-SH.....	76, 105, 196	MH18-750....	92	MH19-282....	106
MH11-684.....	76, 196	MH18-780....	37, 92, 134, 158	MH19-284....	106
MH11-700.....	76, 196	MH18-782....	37, 92, 134, 158	MH19-286....	106
MH11-702.....	76, 196	MH18-784....	37, 92, 134, 158	MH19-288....	106
		MH18-786....	37, 93, 135, 159	MH19-290....	106
MH12-80.....	176, 226	MH18-788....	37, 93, 135, 159	MH19-292....	106
MH14-48.....	51	MH18-790....	37, 93, 135, 159	MH19-294....	106
MH16-4.....	41	MH18-835....	38, 87, 94	MH19-296....	106
MH16-6.....	41	MH18-836....	38, 87, 94	MH19-298....	106
MH16-8.....	41	MH18-866....	38, 94, 135, 159	MH19-314....	106

REF	Page	REF	Page	REF	Page
MH19-316	106	MH20-32	117	MH21-153	84, 144, 167
MH19-318	106	MH20-140	114	MH21-154	84, 144, 167
MH19-350	40, 108	MH20-142	114	MH21-156	84, 144, 167
MH19-370	40, 108	MH20-160	18, 114	MH21-158	84, 144, 167
MH19-410	109	MH20-180	117	MH21-160	84, 144, 167
MH19-411	109	MH20-200	118	MH21-165	78, 144
MH19-434	109	MH20-246	118	MH21-188	78, 148
MH19-570	66, 111, 197	MH20-258	118	MH21-205-2	148
MH19-572	66, 111, 197	MH20-278	118, 146	MH21-205-4	148
MH19-573	66, 111, 197	MH20-280	118, 146	MH21-205-6	148
MH19-574	66, 111, 116, 197	MH20-282	118, 146	MH21-205-8	148
MH19-575	66, 111, 197	MH20-284	118, 146	MH21-205-10	148
MH19-576	66, 111, 197	MH20-504	115	MH21-205-12	148
MH19-577	66, 111, 197	MH20-506	115	MH21-206-4	148
MH19-580	66, 111, 197	MH20-522	115	MH21-206-7	148
MH19-582	66, 111, 197	MH20-552	119	MH21-206-9	148
MH19-584	66, 111, 197	MH20-553	119	MH21-206-12	148
MH19-700	107, 198, 208	MH20-554	119	MH21-232	148
MH19-702	107, 198, 208	MH20-556	119	MH21-233	148
MH19-704	107, 198, 208	MH20-570	119	MH21-234	148
MH19-706	107, 198, 208	MH20-800	120	MH21-280	149
MH19-708	107, 198, 208	MH20-802	120	MH21-300	149
MH19-710	107, 198, 208	MH20-804	120	MH21-330	149
MH19-712	107, 198, 208	MH20-806	120	MH21-338	150
MH19-714	107, 198, 208	MH20-808	120	MH21-342	150
MH19-716	107, 198, 208	MH20-852	120	MH21-346	150
MH19-718	107, 198, 208	MH20-854	120	MH21-348	150
MH19-720	107, 198, 208	MH20-856	120	MH21-368	150
MH19-793	107, 145, 216	MH20-858-1	120	MH21-369	150
MH19-796	107, 145, 216	MH20-858-2	120	MH21-370	150
MH19-844	112, 172, 200, 210	MH20-858-3	120	MH21-440	135
MH19-848	112, 200, 210	MH21-20	146	MH21-536	23, 129
MH19-849	112, 200, 210	MH21-24	146	MH21-600	18, 130
MH19-850	112, 186, 200, 210	MH21-32	146	MH21-602	18, 130
MH19-852	112, 186, 200, 210	MH21-60	147	MH21-604	18, 130
MH19-1230	112	MH21-64	147	MH21-610	19, 131
MH19-1232	112	MH21-72	147	MH21-626	151, 211
MH19-2081	88, 110	MH21-88	82, 143, 166	MH21-627	151, 211
MH19-2110	88, 110	MH21-89	82, 143, 166	MH21-630	19, 131, 156
MH20-2	116	MH21-90	83, 143, 166	MH21-632	19, 131, 156
MH20-3	116	MH21-92	83, 143, 166	MH21-660	151
MH20-4	116	MH21-94	83, 143, 166	MH22-2	172
MH20-12	117	MH21-144	83, 142, 166	MH22-4	172
MH20-24	117	MH21-145	77, 83, 142, 166	MH22-550	174
MH20-28	117	MH21-146	77, 83, 142, 166	MH22-552	174
MH20-29	117	MH21-148	77, 142	MH22-662	52, 174
MH20-30	117	MH21-152	83, 91, 144, 167	MH22-668	52, 174

REF	Page	REF	Page	REF	Page
MH22-808.....	175	MH26-430.....	203	MH27-447C.....	215
MH22-854.....	78, 175	MH26-435.....	203	MH27-448.....	215
MH22-1000.....	175	MH26-956.....	40, 191	MH27-448C.....	215
MH23-1049.....	78, 142, 176	MH26-958.....	40, 191	MH27-488.....	215
MH23-1050.....	84, 177	MH26-1084.....	85, 142, 165, 204	MH27-490.....	215
MH23-1051.....	84, 177	MH26-1090.....	85, 204	MH27-492.....	215
MH23-1420.....	177	MH26-1410.....	204, 216	MH27-494.....	215
MH23-1430.....	177	MH26-1412.....	204, 216	MH27-496.....	215
MH24-800.....	178	MH26-1450.....	205	MH27-499.....	215
MH24-802.....	178	MH26-1451.....	205	MH27-500.....	215
MH24-804.....	178	MH26-1452.....	205	MH27-502.....	215
MH24-806.....	178	MH26-1453.....	205	MH27-504.....	215
MH24-808.....	178	MH26-1454.....	205	MH27-506.....	215
MH24-1180.....	178	MH26-1627.....	199, 209	MH27-508.....	215
MH24-1190.....	178	MH26-1628.....	199, 209	MH27-700.....	217
MH24-1200.....	178	MH26-1629.....	199, 209	MH27-702.....	217
MH24-3200.....	179	MH26-1630.....	199, 209	MH27-820.....	217
MH24-3210.....	179	MH26-1632.....	199, 209	MH27-822.....	217
MH24-3220.....	179	MH26-1634.....	199, 209	MH27-824.....	217
MH24-3230.....	179	MH26-1636.....	199, 209	MH27-826.....	217
MH24-3240.....	179	MH26-1638.....	199, 209	MH27-990.....	218
MH24-3250.....	179	MH26-1640.....	199, 209	MH27-1000.....	20, 218
MH24-3260.....	180	MH27-4.....	213	MH27-1016.....	218
MH24-3270.....	180	MH27-7.....	213	MH27-3070.....	25, 214
MH25-384.....	186, 211	MH27-16.....	213	MH27-3100.....	214
MH25-390.....	187, 211	MH27-30.....	214	MH28-1.....	220
MH25-392.....	187, 211	MH27-32.....	214	MH28-2.....	220
MH25-396.....	187, 212	MH27-34.....	214	MH28-20.....	221
MH25-454.....	187, 201, 212	MH27-160.....	214	MH28-22.....	221
MH25-456.....	187, 201, 212	MH27-213.....	79, 207	MH28-24.....	221
MH25-458.....	187, 201, 212	MH27-214.....	79, 207	MH28-68.....	80, 221
MH25-480.....	188	MH27-215.....	79, 207	MH28-69.....	80, 221
MH25-510.....	188, 201, 212	MH27-216.....	79, 207	MH28-196.....	220
MH25-511.....	188, 201, 212	MH27-217.....	79, 207	MH29-6-8.....	222
MH25-800.....	52, 183	MH27-218.....	79, 207	MH29-6-10.....	222
MH25-814.....	53, 183	MH27-219.....	79, 207	MH29-6-12.....	222
MH25-836.....	53, 183	MH27-438.....	215	MH29-6-14.....	222
MH25-837.....	53, 183	MH27-438C.....	215	MH29-6-16.....	222
MH25-1160.....	183	MH27-440.....	215	MH29-6-18.....	222
MH26-119.....	202, 213	MH27-440C.....	215	MH29-6-20.....	222
MH26-124.....	202, 213	MH27-442.....	215	MH29-6-22.....	222
MH26-336.....	202	MH27-442C.....	215	MH29-6-24.....	222
MH26-339A.....	202	MH27-444.....	215	MH29-6-26.....	222
MH26-394.....	203	MH27-444C.....	215	MH29-6-28.....	222
MH26-397.....	203	MH27-446.....	215	MH29-6-30.....	222
MH26-400.....	203	MH27-446C.....	215	MH29-6-32.....	222
MH26-402.....	203	MH27-447.....	215	MH29-6-34.....	222



<b>REF</b>	<b>Page</b>	<b>REF</b>	<b>Page</b>	<b>REF</b>	<b>Page</b>
MH29-6-36	222	MH30-555-1112	229	MH30-1383	234
MH29-284	223	MH30-555-1314	229	MH30-1385	234
MH29-286	223	MH30-555-1516	229	MH30-1390	234
MH29-288	223	MH30-555-1718	229	MH30-1395	234
MH29-290	223	MH30-560-1315	230	MH30-1405	235
MH29-302	223	MH30-560-1719	230	MH30-1442	235
MH29-308	223	MH30-560-2123	230	MH30-1443	235
MH30-10	227	MH30-560-2527	230	MH30-1482	236
MH30-15	227	MH30-560-2931	230	MH30-1485	236
MH30-20	227	MH30-560-3335	230	MH30-1700	54, 237
MH30-30	227	MH30-560-3739	230	MH30-1710	54, 237
MH30-50	227	MH30-560-4143	230	MH30-1716	54, 237
MH30-55	227	MH30-561	230	MH30-1718	54, 237
MH30-60	227	MH30-965	231	MH30-2190	20, 237
MH30-62	227	MH30-965ATR	231	MH30-2195	20, 237
MH30-395	80, 225	MH30-966	231	MH33-140	167
MH30-406	80, 225	MH30-966ATR	231	MH33-142	167
MH30-415	80, 225	MH30-967	231	MH33-200	167
MH30-500	228	MH30-985	231	MH33-201	167
MH30-505-910	228	MH30-997	232	MH33-202	167
MH30-505-1112	228	MH30-1000	232	MH40-210-SS	219
MH30-505-1314	228	MH30-1025	54, 232	MH40-212-SS	219
MH30-505-1516	228	MH30-1205-0	233	MH40-224	219
MH30-505-1718	228	MH30-1205-1	233	MH40-225-SS	219
MH30-505-1920	228	MH30-1205-2	233	MH40-226A	219
MH30-550	229	MH30-1205-3	233	MH40-226-SS	219
MH30-555-34	229	MH30-1205-4	233	MH40-227-SS	219
MH30-555-56	229	MH30-1205-5	233	MH40-230	219
MH30-555-78	229	MH30-1205-6	233	MH40-245-SS	219
MH30-555-910	229	MH30-1382	234	MH40-250-SS	219



[www.miltex.com](http://www.miltex.com)  
toll-free phone 866-854-8300  
toll-free fax 866-854-8400  
phone 717-840-9335  
fax 717-840-9347

**Please call Customer Service today for additional information, an in-office demonstration, to place an order, or to receive other Miltex catalogs.**

Miltex, Miltex Redefining Excellence and MeisterHand are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation.

The Integra wave logo is a trademark of Integra LifeSciences Corporation.

Copyright© 2009 Integra LifeSciences Corporation

CAT-MEISTHAND      MPM0401

1009/6/C/4



\*+H834CATMEISTHAND0/\*