

MILTEX® HI-LIGHT

ORTHOPEDIC INSTRUMENTS

EXPANDED OFFERING!



Miltex announces the expansion of our Orthopedic Instrumentation line. These additions, which are detailed on the reverse side, include Bone Reduction Forceps, Cobb Curettes, Bone Tamps, Retractors, Synthetic Cast Cutters and more!


See reverse side for more Orthopedic Instruments

Illustrations and content provide general description only and may be subject to change. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra Life Sciences Corporation. Copyright © 2010 Integra LifeSciences Corporation

CURETTES

REF	Description
19-721	COBB Curette 11" (27.9 cm), straight, oval cup, Size 3/0
19-722	COBB Curette 11" (27.9 cm), straight, oval cup, Size 2/0
19-723	COBB Curette 11" (27.9 cm), straight, oval cup, Size 0
19-724	COBB Curette 11" (27.9 cm), straight, oval cup, Size 1
19-725	COBB Curette 11" (27.9 cm), straight, oval cup, Size 2
19-726	COBB Curette 11" (27.9 cm), straight, oval cup, Size 3
19-727	COBB Curette 11" (27.9 cm), straight, oval cup, Size 4
19-728	COBB Curette 11" (27.9 cm), straight, oval cup, Size 5
19-729	COBB Curette 11" (27.9 cm), straight, oval cup, Size 6
19-737	COBB Curette 11" (27.9 cm), angled, oval cup, Size 3/0
19-738	COBB Curette 11" (27.9 cm), angled, oval cup, Size 2/0
19-730	COBB Curette 11" (27.9 cm), angled, oval cup, Size 0
19-731	COBB Curette 11" (27.9 cm), angled, oval cup, Size 1
19-732	COBB Curette 11" (27.9 cm), angled, oval cup, Size 2
19-733	COBB Curette 11" (27.9 cm), angled, oval cup, Size 3
19-734	COBB Curette 11" (27.9 cm), angled, oval cup, Size 4
19-735	COBB Curette 11" (27.9 cm), angled, oval cup, Size 5
19-736	COBB Curette 11" (27.9 cm), angled, oval cup, Size 6
26-1728	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 00
26-1729	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 0
26-1730	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 1
26-1731	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 2
26-1732	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 3
26-1733	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 4
26-1734	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 5
26-1735	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 6
26-1736	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 2/0
26-1737	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 3/0
26-1738	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 4/0
26-1739	SPRATT Curette 9" (22.9 cm), angled, oval cup, Size 5/0
26-1700	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 000, 2.8 mm
26-1701	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 00, 3.3 mm
26-1702	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 0, 3.7 mm
26-1703	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 1, 4.3 mm
26-1704	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 2, 4.8 mm
26-1705	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 3, 5.6 mm
26-1706	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 4, 6.1 mm
26-1707	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 5, 6.7 mm
26-1708	BRUN (HIBBS-SPRATT) Curette 9" (22.9 cm), reverse angle, oval cup, size 6, 8.7 mm




19-721 to 19-729

19-731 to 19-738


26-1728 to 26-1739

26-1700 to 26-1708



BONE TAMPS

27-1040	Bone Tamp 6-1/4" (15.9 cm), 2 mm wide, serrated tip
27-1042	Bone Tamp 6-1/4" (15.9 cm), 3 mm wide, serrated tip
27-1044	Bone Tamp 6-1/4" (15.9 cm), 4 mm wide, serrated tip
27-1046	Bone Tamp 6-1/4" (15.9 cm), 5 mm wide, serrated tip
27-1048	Bone Tamp 6-1/4" (15.9 cm), 6 mm wide, serrated tip
27-1050	Bone Tamp 6-1/4" (15.9 cm), 8 mm wide, serrated tip
27-1052	Bone Tamp 6-1/4" (15.9 cm), 10 mm wide, serrated tip
27-2000	Bone Graft Impactor 6-1/2" (16.5 cm), double-ended



27-1040 to 27-1052



BONE HOLDING

27-80	VERBRUGGE Bone Forceps 10" (25.4 cm), 10 mm jaws
27-70	VERBRUGGE Bone Forceps 7" (17.8 cm), 5 mm jaws
27-26	Bone Reduction Forceps 5-1/4" (13.3 cm), ratchet, pointed jaw with step
27-28	Bone Reduction Forceps 8" (20.3 cm), ratchet, pointed jaw
27-93	Bone Reduction Forceps 5-1/2" (14 cm), curved, with ratchet
27-95	Bone Reduction Forceps 8" (20.3 cm), curved, with ratchet*



27-80

MISCELLANEOUS

27-2040	COBB Elevator 9-1/2" (24.1 cm), 10 mm wide
27-2030	COBB Elevator 9-1/2" (24.1 cm), 13 mm wide
27-2020	COBB Elevator 9-1/2" (24.1 cm), 19 mm wide
27-2010	COBB Elevator 9-1/2" (24.1 cm), 25 mm wide
27-2001	COBB Elevator 9-1/2" (24.1 cm), 32 mm wide
27-3000	COBB Gouge 11" (27.9 cm), straight
27-995	AUFRANC Cobra Retractor 10-1/4" (26 cm), 32 mm wide pointed tip, serrated
27-996	AUFRANC Cobra Retractor 10-1/4" (26 cm), 32 mm wide blunt tip, serrated
26-1350	D'ERRICO Root Retractor 9" (22.9 cm), curved
21-200	JOHNSON Skin Hook 6-1/2" (16.5 cm), sharp
25-512	LEKSELL Rongeur Forceps 9-1/2" (24.1 cm), 5 mm bite, double action*
27-652	Wire Pulling Forceps 6-1/2" (16.5 cm)*
27-162	Pin Cutter, Heavy Duty 22" (25.8 cm)*
9-129	Double Action Pliers & Wire Cutter 5" (12.7 cm), stainless steel*
27-3018	Synthetic Cast Cutter 6-3/4" (17 cm), sharp*
27-3019	Synthetic Cast Cutter 7-1/4" (18.4 cm), blunt*
27-83	HATT Spoon 9" (22.9 cm), 10 mm oval cup
27-84	HATT Spoon 9" (22.9 cm), 15 mm oval cup
27-85	HATT Spoon 9" (22.9 cm), 20 mm oval cup



27-2001 to 27-2040

27-995 to 27-996

21-200

27-83 to 27-85

* see image on reverse side